

## Operator Preventive Maintenance Checks and Services

6630-01-376-9823  
Analyzer, Clinical Chemistry, DT60

[B-Before Operation, D-During Operation, A-After Operation, Q-Quarterly, and S-Semiannually]

ITEM NO	INTERVAL	ITEM TO BE INSPECTED AND PROCEDURE	IS NOT MISSION CAPABLE IF:
1	B	<p><b>Analyzer, Clinical Chemistry</b></p> <p>a. Conduct an inventory to ensure that the items listed on the Equipment Parts and Accessories List are on hand.</p> <p>b. Inspect the unit for damage, discoloration, or excessively worn components.</p> <p>c. Verify that the Verification/Certification sticker (DD Form 2163) has a current date.</p>	<p>Missing components or accessories prevent the operation of the unit.</p> <p>Unserviceable components prevent the use of the unit.</p> <p>The unit has not been verified within the last six (6) months.</p>
2	B, D	<p><b>Installation</b></p> <p>a. Set up equipment as directed by the manufacturer's literature.</p> <p>b. Check power source and plug the analyzer in as directed by the manufacturer's literature.</p>	
3	B, D	<p><b>Operating Instructions</b></p> <p>a. Perform startup procedures as directed by the manufacturer's operator's manual.</p> <p>b. Perform testing procedure as directed by the manufacturer's operator's manual.</p> <p style="padding-left: 20px;">(1) Conduct "Pipetting Techniques" as directed by the operator's manual.</p> <p style="padding-left: 20px;">(2) Perform the "Steps for Analysis on the Vitros DT II System" as directed.</p> <p>c. Perform the "Calibration Data Module and Chemistry Language Module" as directed by the operator's manual.</p> <p>d. Perform the normal shutdown procedure as directed by the operator's manual.</p> <p style="padding-left: 20px;">(1) Check incubator for tests in progress.</p> <p style="padding-left: 20px;">(2) Turn analyzer off.</p> <p><b>NOTE: Conducting an EMERGENCY SHUTDOWN will result in having to repeat the analysis for all tests that remained in the incubator at the time of the shutdown.</b></p>	<p>The unit cannot be assembled or will not start up properly.</p>





# Operator Preventive Maintenance Checks and Services

6630-01-376-9823  
Analyzer, Clinical Chemistry, DT60

[B-Before Operation, D-During Operation, A-After Operation, Q-Quarterly, and S-Semiannually]

ITEM NO	INTERVAL	ITEM TO BE INSPECTED AND PROCEDURE	IS NOT MISSION CAPABLE IF:
		<p><b>CAUTION: NEVER LEAVE THE CHARGER CONNECTED TO THE DT PIPETTE FOR MORE THAN 72 HOURS.</b></p> <p>e. Perform other cleaning as directed by the manufacturer's operator's manual: "Vitros DT60 II Chemistry System FORS Head."</p>	