

SB 8-75-S8

DEPARTMENT OF THE ARMY SUPPLY BULLETIN

Army Medical Department Supply Information

Headquarters, Department of the Army, Washington, DC 20310-2300

20 August 2013

Effective until rescinded or superseded

DISTRIBUTION STATEMENT A: Approval for public release; distribution is unlimited

Printing of this information does not constitute endorsement of the product, manufacturer, or the use of any related item thereof by the US Army.

REPORTING ERRORS AND RECOMMENDING IMPROVEMENTS

You can improve this publication. If you find any mistakes, or you know a way to improve it, please let us know. Mail your letter or DA Form 2028 (Recommended Changes to Publications and Blank Forms) directly to the US Army Medical Materiel Agency, ATTN: MCMR-MMP-M, Fort Detrick, MD 21702-5001.

| Table of Contents | | Page |
|---------------------------|-----------------------------------------------------------------------|------|
| CHAPTER 1 | - Preventive Maintenance Checks and Services..... | 1-1 |
| CHAPTER 2 | - Supportability Strategy Overview | 2-1 |
| APPENDIX A | - Operator PMCS | A-1 |
| APPENDIX B | - Repairer PMCS | B-1 |
| APPENDIX C | - Maintenance Allocation Chart | C-1 |
| APPENDIX D | - Equipment Parts and Accessories List | D-1 |
| APPENDIX E | - Tools and TMDE Code Listing for Maintenance Allocation Charts | E-1 |
| 2012 Index for SB 8-75-S8 | - | IN-1 |

NOTICE

This issue of the DA SB 8-75-S8 is devoted entirely to the Medical Equipment Supportability Information developed by the Medical Devices Program Management Office, USAMMA, and supersedes any previous editions.

CHAPTER 1. PREVENTIVE MAINTENANCE CHECKS AND SERVICES

1-1. PREVENTIVE MAINTENANCE CHECKS AND SERVICES (PMCS)

a. The call for stringent application of PMCS procedures is issued throughout the Army Medical Department (AMEDD). Your vigilance is needed in following the required procedures in maintaining medical equipment. Medical equipment must be inspected and serviced systematically and periodically to ensure that it is ready for operation at all times. Inspections can reveal defects that can be corrected before they result in serious damage or failure.

b. Medical equipment will be maintained to 10/20 standards as published in equipment technical manuals (TMs) or to the operational standards published in the manufacturer's operator or service literature. In the absence of a TM, the standards established in Appendix A and Appendix B of this supply bulletin, in conjunction with the requirements of the manufacturer's literature, will be used.

c. Complete PMCS will be performed before placing Army equipment in administrative storage. When equipment is removed from storage, perform PMCS to ensure its operational readiness. Results of the equipment inspection will be documented on DA Form 2404, *Equipment Inspection and Maintenance Worksheet*. All discrepancies will be annotated with corrective action required and steps taken to correct the deficiency.

d. This SB 8-75-S8 is used to update 10/20 standards for medical equipment and will include additional standards for reportable items and newly fielded medical equipment. Each new publication of SB 8-75-S8 supersedes the previous year's edition. Equipment specific 10/20 PMCS and Maintenance Allocation Chart (MAC) updates will also be available through the USAMMA website at www.usamma.amedd.army.mil. The 10/20 PMCS and MAC update is an ongoing project. New releases will be posted as they are developed.

e. SB 8-75-S8 will be used when discrepancies are encountered with the 10/20 standards published in the equipment's TM. Standards published in the manufacturer's literature take precedence over SB 8-75-S8.

1-2. OPERATOR LEVEL PMCS

a. Appendix A contains a list of tasks to be performed by unit level operator/user personnel. These PMCS tables are also referred to as "10 PMCS" requirements. Preventive maintenance by operator/user personnel is not limited to performing the checks and services listed in the tables. There are tasks the operator/user must do any time the equipment is used or stored, such as checking for general cleanliness, observing for improper operational indicators, and maintaining the proper quantities of operating supplies.

b. The following is a list of PMCS table column headings with a description of the information found in each column:

(1) Item Number. This column shows the sequence to perform the PMCS, and is used to identify the equipment area on the DA Form 2404.

(2) Interval. This column shows when each PMCS item is to be serviced. B, D, and A should be performed with daily use of the equipment.

| | |
|----------------------|----------------------|
| B – Before Operation | D – During Operation |
| A – After Operation | Q – Quarterly |
| S - Semiannually | |

(3) Item to be Inspected and Procedure. This column identifies the general area or specific part to be checked or serviced.

(4) Equipment is not Mission Capable If. This column lists conditions that make the equipment unavailable or unusable.

c. When the equipment must be kept in continuous operation, check and service only those items that will not disrupt operation. Perform the complete daily checks and services when the equipment can be shut down.

d. Operator/user personnel will report problems with medical equipment discovered during their "10 PMCS" that they are unable to correct. Refer to TB 38-750-2, *Maintenance Management Procedures for Medical Equipment*, and report the deficiency using the proper forms. Consult with your unit's Medical Equipment Repairer if you need assistance.

e. Table 1-1 lists the "10 PMCS" standards located in Appendix A.

TABLE 1-1. 10 PMCS STANDARDS LOCATED IN APPENDIX A

| NSN | EQUIPMENT NAME AND MODEL | PAGE |
|------------------|------------------------------------------------------|------|
| N/A | Generic Operator Standards | A-1 |
| 6520-01-272-4531 | Dental Operating Unit, ADEC Model 3406 Porta-Cart | A-2 |
| 6520-01-398-4613 | Compressor Dehydrator, Dental, Model PAC 6.7 | A-4 |
| 6525-01-312-6411 | X-Ray Apparatus, Radiographic/Fluoroscopic, CS-8952 | A-6 |
| 6525-01-325-3740 | Portable X-Ray System, Model 1200 | A-9 |
| 6525-01-370-7552 | Portable Dental X-Ray System, Model ALPHA MPDX | A-12 |
| 6525-01-384-9296 | X-Ray Apparatus, Model LCROKS | A-14 |
| 6530-00-926-2151 | Sterilizer, Surgical Dressing 16x36 in., Model M-138 | A-15 |

1-3. REPAIRER LEVEL PMCS

a. Appendix B contains a list of tasks to be performed by the unit level repairer. These PMCS tables are also referred to as "20 PMCS" requirements.

b. The following is a list of 20 PMCS table column headings with a description of the information found in each column:

(1) Item Number. This column shows the sequence to perform the PMCS, and is used to identify the equipment area on DA Form 2404.

(2) Interval. This column shows when each PMCS item is to be serviced:
M – Monthly, Q – Quarterly, S – Semiannually, and A – Annually.

(3) Item to be Inspected and Procedure. This column identifies the general area or specific part to be checked or serviced.

(4) Equipment is not Mission Capable If. This column lists conditions that make the equipment unavailable or unusable.

c. When the equipment must be kept in continuous operation, check and service only those items that will not disrupt operation. Perform the complete daily checks and services when the equipment can be shut down.

d. The following list in Table 1-2 is the "20 PMCS" charts located in Appendix B. This table identifies the NSN, the name and model of the equipment, and the page number.

TABLE 1-2. 20 PMCS CHARTS LOCATED IN APPENDIX B

| NSN | EQUIPMENT NAME AND MODEL | PAGE |
|------------------|-----------------------------------------------------------|------|
| N/A | Generic Repairer Standards | B-1 |
| 6520-01-272-4531 | Dental Operating Unit, ADEC Model 3406 Porta-Cart | B-2 |
| 6520-01-398-4613 | Compressor Dehydrator, Dental, Model PAC 6.7 | B-4 |
| 6525-01-312-6411 | X-Ray Apparatus, Radiographic/Fluoroscopic, Model CS-8952 | B-5 |
| 6525-01-325-3740 | Portable X-Ray System, Model 1200 | B-7 |
| 6525-01-370-7552 | Portable Dental X-Ray System, Model ALPHA MPDX | B-9 |
| 6525-01-384-9296 | X-Ray Apparatus, Model LCROKS | B-11 |
| 6530-00-926-2151 | Sterilizer, Surgical Dressing 16x36 in., Model M-138 | B-14 |

1-4. MAINTENANCE ALLOCATION CHART (MAC)

a. Appendix C provides a general explanation of all maintenance and repair functions authorized at various maintenance levels.

b. The following is a list of MAC table column headings with a description of the information found in each column:

(1) Group Number. This column is a numerical group assigned to each assembly. The applicable assembly groups are listed in the MAC in disassembly sequence beginning with the first assembly removed in a top down disassembly sequence.

(2) Assembly Group. This column contains a brief description of the components of each assembly group.

(3) Maintenance Function. This column lists the various maintenance functions authorized to be performed. These maintenance functions are defined as follows:

(a) Inspect. To determine serviceability of an item by comparing its physical, mechanical, and electrical characteristics with established standards.

(b) Test. To verify serviceability and to detect electrical or mechanical failure using test equipment.

(c) Service. To clean, to preserve, to charge, and to add lubricants, cooling agents, and air. If it is desired that elements, such as painting and lubricating, be defined separately, they may be so listed.

(d) Adjust. To rectify to the extent necessary to bring device into proper operation range.

(e) Align. To adjust specified variable elements of an item to bring it to optimum performance.

(f) Calibrate. To determine the corrections to be made in the readings of instruments or test equipment used in precise measurement. Consists of the comparison of two instruments, one of which is a certified standard of known accuracy, to detect and adjust any discrepancy in the accuracy of the instrument being compared with the certified standard.

(g) Install. To set for use in an operational environment such as tents or International Standards Organization shelters.

(h) Replace. To replace unserviceable items with comparable serviceable items.

(i) Repair. Those maintenance operations necessary to restore an item to serviceable condition through correction of material damage to a specific failure. Repair may be accomplished at each level of maintenance.

(j) Overhaul. Normally the highest degree of maintenance performed by the Army in order to minimize time work in process consistent with quality and economy of operation. It consists of that maintenance necessary to restore an item to completely serviceable condition as prescribed by a maintenance standard in technical publications for each item of equipment. Overhaul normally does not return an item to like-new condition.

(k) Rebuild. The highest degree of material maintenance, it consists of restoring equipment as nearly as possible to new condition in accordance with original manufacturing standards. Rebuild is performed only when required by operational considerations or other paramount factors and then only at the depot maintenance level.

(4) Maintenance Level. This column indicates the maintenance level authorized to perform the maintenance functions.

| | |
|----------------------------------|-------------------------------------|
| "FM" is crew or unit maintenance | "SM" is MEDLOG or depot maintenance |
|----------------------------------|-------------------------------------|

(5) Tools and Equipment. This column corresponds to tools and, test, measurement, and diagnostic equipment (TMDE) listed in the chart in Appendix E.

(6) Remarks. This column is provided for information pertinent to the maintenance functions

c. The following is a list of the MACs located in Appendix C.

TABLE 1-3. LIST OF THE MACS LOCATED IN APPENDIX C

| NSN | EQUIPMENT NAME AND MODEL | PAGE |
|------------------|-----------------------------------------------------------|------|
| 6520-01-272-4531 | Dental Operating Unit, ADEC Model 3406 Porta-Cart | C-1 |
| 6520-01-398-4613 | Compressor Dehydrator, Dental, Model PAC 6.7 | C-3 |
| 6525-01-312-6411 | X-Ray Apparatus, Radiographic/Fluoroscopic, Model CS-8952 | C-5 |
| 6525-01-325-3740 | Portable X-Ray System, Model 1200 | C-8 |
| 6525-01-370-7552 | Portable Dental X-Ray System, Model ALPHA MPDX | C-13 |
| 6525-01-384-9296 | X-Ray Apparatus, Model LCROKS | C-15 |
| 6530-00-926-2151 | Sterilizer, Surgical Dressing 16x36 in., Model M-138 | C-18 |

d. Appendix D contains the equipment parts and accessories list for each item of equipment. The following is a list of the items located in Appendix D.

TABLE 1-4. LIST OF PARTS AND ACCESSORIES LOCATED IN APPENDIX D

| NSN | EQUIPMENT NAME AND MODEL | PAGE |
|------------------|-----------------------------------------------------------|------|
| 6520-01-272-4531 | Dental Operating Unit, ADEC Model 3406 Porta-Cart | D-1 |
| 6520-01-398-4613 | Compressor Dehydrator, Dental, Model PAC 6.7 | D-3 |
| 6525-01-312-6411 | X-Ray Apparatus, Radiographic/Fluoroscopic, Model CS-8952 | D-4 |
| 6525-01-325-3740 | Portable X-Ray System, Model 1200 | D-5 |
| 6525-01-370-7552 | Portable Dental X-Ray System, Model ALPHA MPDX | D-6 |
| 6525-01-384-9296 | X-Ray Apparatus, Model LCROKS | D-7 |
| 6530-00-926-2151 | Sterilizer, Surgical Dressing 16x36 in., Model M-138 | D-9 |

e. Appendix E contains the tools and TMDE code listing for MACs

CHAPTER 2. SUPPORTABILITY STRATEGY OVERVIEW

2-1. REQUIREMENTS

In accordance with *AR 700-127, Integrated Logistics Support*, the Materiel Developer (USAMMA for COTS items) is responsible for developing, coordinating, and maintaining the supportability strategy throughout the life of the system. It is required for all equipment in our unit assemblages.

2-2. RESPONSIBILITY

Key areas of the USAMMA's Acquisition and Supportability Strategy include the acquisition strategy, clinical application, requisition and supply information, and the 12 integrated logistics support elements: maintenance planning and management; training and training support; manpower and personnel; supply support; technical data; facilities; product support management; sustaining engineering; packaging, handling, storage, and transportation; design interface; support equipment; and computer resources.

a. The Acquisition Strategy notes what the item is, the unit assemblages and types of units it is assigned to, the cost, and life expectancy.

b. The Clinical Application section provides a high-level overview regarding the clinical use of the system.

c. Requisition/Supply notes the cataloging terms regarding how the item is procured, stored, reviewed and resupplied.

d. The 12 integrated logistics support elements are described below.

(1) Maintenance Planning includes high-level user and depot level procedures and requirements. There is a separate maintenance plan which carries maintenance and inspection procedures, tools and TMDE requirements, 10/20 standards, maintenance allocation charts (MAC), MEDSTEP requirements, repair and spare parts listing, man-hour requirements, warranty procedures, and DA 2406 reportable information.

(2) Training and Training Support notes how both maintainer and operator training are performed.

(3) Manpower and Personnel provides an overview of the military occupation specialties (MOS) required to operate and maintain the system.

(4) Supply Support provides effective and efficient end-to-end customer service to meet operational requirements for all classes of supply.

(5) Support Equipment provides a list of the support items and how to order these items.

(6) Technical Data provides information on the operator and service manuals.

(7) Computer Resource Support notes whether a computer or special software is required separately from the system for updates or operation.

(8) The Facilities section explains if any additional facility requirements are required for storage or operation.

(9) Packaging, Handling, Storage, and Transportation annotates if the system is procured and transported in a hardened case, if there are any special package or storage requirements, item shelf life, and special characteristics.

(10) Design Interface comments on any special interfaces required for system operation.

(11) Product Support Management's purpose is to plan, manage, and fund the system product support across all Integrated Product Support (IPS) elements.

(12) Sustaining Engineering spans those technical tasks (engineering and logistics).

2-3. CUSTOMERS

Customers for these documents include, but are not limited to the logistician, the maintainer, and the operator.

2-4. EQUIPMENT SUPPORT PLANS AVAILABILITY

Equipment support plans are available on the USAMMA website under "Equipment Supportability Plans" which are listed by NSN and name. Check the ones required and enter your email address to immediately receive an email with the support plan attached. Only personnel with a ".mil" address are authorized to receive the support plans. You can access this page by going to http://www.usamma.amedd.army.mil/assets/apps/email_doc/email_doc_index.cfm. The USAMMA PMs are developing equipment maintenance plans which will be posted to the USAMMA website starting in FY14. As the maintenance plans are posted the corresponding support plans will be removed from the website.

Appendix A. Operator PMCS

Generic Standards

[B-Before Operation, D-During Operation, A-After Operation, Q-Quarterly, and S-Semiannually]

| ITEM NO | INTERVAL | ITEM TO BE INSPECTED AND PROCEDURE | IS NOT MISSION CAPABLE IF: |
|---------|----------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------|
| 1 | B, A | Ensure that all supplies, reagents, and ancillary components necessary to operate the equipment or system are on hand. | Supplies, reagents, or ancillary components are missing. |
| 2 | B | Inspect for corrosion, rust, physically damaged parts, deteriorated materials, and damage to protective coatings. | Rust on outer surface parts determined by the Infection Control Nurse to be a health hazard. |
| 3 | B, A | Ensure the operator manual or documentation is on hand. Identify the location of such material if it is not packed with the equipment. | Operator manual is not readily available. |
| 4 | B, A | Verify that the equipment or system has no broken parts or accessories, i.e., switches, knobs, casters, plastic coverings, hoses, casings, etc. | Equipment is not functional due to broken parts. |
| 5 | B, D, A | Ensure that fluid levels, lubricants, physical limits or settings for operation are correct. | Levels are below those established in the TM or manufacturer's literature. |
| 6 | B | Verify date of last electrical safety test, PMCS, or CVC services from DA Form 2163 or other record (typically, annual inspection for patient care, laboratory and incidental use; semi-annual inspection for critical care and anesthetizing locations). If beyond designated period, arrange for CVC services. | Performance of CVC cannot be verified or CVCs are past due. |
| 7 | B, D, A | Verify operation of the equipment or system in accordance with published TMs and manufacturer's literature. | Equipment fails to operate according to TM or manufacturer's specifications. |
| 8 | B, D, A | Inspect for unusual operation, noises, leakage, or other unexpected results. | Noticeable fluid leaks or unexpected noises are detected. |
| 9 | A | Shut down equipment and clean and dry parts or components that were subjected to liquid contact. | Unit or components are not clean or dry. |
| 10 | A | Locate and store components, accessories, and operator documentation with the equipment or in appropriate location. | Items are not stored with the equipment or are not readily available. |
| 11 | B, A | Check the electrical power cord for cuts, fraying, or deterioration. | Electrical plug is missing a pin/blade or the cord insulation is cut through the outer coating. |
| 12 | B, D | Ensure that alarms and visual indicators are functioning properly. | Alarms and indicators are not functioning properly. |
| 13 | B, A | Inspect storage container for damage to case, hinges, latches, and seals. | Storage container will not latch or could leak. |

6520-01-272-4531
 Dental Operating Unit, ADEC Model 3406 Porta-Cart

[B-Before Operation, D-During Operation, A-After Operation, Q-Quarterly, and S-Semiannually]

| ITEM NO | INTERVAL | ITEM TO BE INSPECTED AND PROCEDURE | IS NOT MISSION CAPABLE IF: |
|---------|----------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1 | B, A | <p>Dental Unit</p> <p>a. Conduct an inventory to ensure that the items listed on the Equipment Parts and Accessories List are on hand.</p> <p>b. Inspect components for damage, discoloration, or excessively worn components.</p> | <p>Missing components or accessories prevent operation of the dental unit.</p> <p>Unserviceable components prevent the use of the dental unit.</p> |
| 2 | B | <p>Dental Unit Operational Test</p> <p>a. Ensure the unit is assembled by performing the equipment setup procedure as directed by the manufacturer's literature.</p> <p>b. Verify the function of the controls according to the manufacturer's literature.</p> <p>c. Prepare the dental unit for use according to the manufacturer's literature.</p> <p>d. Verify the function of the syringe according to the manufacturer's literature.</p> <p>e. Verify the function of the air vacuum system according to the manufacturer's literature.</p> <p>f. Verify the function of the water tank according to the manufacturer's literature.</p> | <p>Missing components prevent the assembly of the unit.</p> <p>Broken controls prevent effective patient care.</p> <p>The dental unit does not maintain air pressure between 60 psi to 80 psi or water pressure between 30 psi to 40 psi.</p> <p>The syringe does not pass water and/or air.</p> <p>The air vacuum system does not create sufficient vacuum.</p> <p>The water tank cannot be pressurized.</p> |
| 3 | B | <p>Dental Handpieces</p> <p>a. Adjust the maximum dynamic air pressure according to the handpiece manufacturer's literature.</p> <p>b. Adjust the water coolant flow according to the manufacturer's literature.</p> | <p>The maximum dynamic air pressure cannot be reached for the particular manufacturer's handpiece.</p> <p>The handpiece coolant water cannot be adjusted.</p> |
| 4 | A | <p>Dental Unit Care</p> <p>After each patient, clean and disinfect all surfaces to include the air vacuum system.</p> | |
| 5 | A | <p>Dental Unit Shut Down</p> <p>Perform the "System Shut-Down" according to the manufacturer's literature.</p> | |

(continued) Appendix A. Operator PMCS

6520-01-272-4531
 Dental Operating Unit, ADEC Model 3406 Porta-Cart

[B-Before Operation, D-During Operation, A-After Operation, Q-Quarterly, and S-Semiannually]

| ITEM NO | INTERVAL | ITEM TO BE INSPECTED AND PROCEDURE | IS NOT MISSION CAPABLE IF: |
|---------|----------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------|
| 6 | A | <p>Dental Unit Storage and Transportation</p> <p>a. Prepare the unit for storage or transportation according to the manufacturer's literature.</p> <p>b. Repack the unit according to the manufacturer's literature.</p> | |
| 7 | B, A | <p>Storage Case</p> <p>Inspect the storage case for cracks, dents, or broken latches.</p> | |

6520-01-398-4613
Compressor Dehydrator, Dental, Model PAC 6.7

[B-Before Operation, D-During Operation, A-After Operation, Q-Quarterly, and S-Semiannually]

| ITEM NO | INTERVAL | ITEM TO BE INSPECTED AND PROCEDURE | IS NOT MISSION CAPABLE IF: |
|---------|----------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1 | B, D | <p>Compressor Dehydrator</p> <p>a. Conduct an inventory to ensure that the items listed on the Equipment Parts and Accessories List are on hand.</p> <p>b. Inspect the unit for any damaged or deteriorated hoses, tubes, cables, and other components.</p> <p>c. Inspect the unit for an excessive accumulation of dust or dirt. Ensure intake filter elements are clean and serviceable.</p> <p>d. Inspect muffler on the water separator for serviceability.</p> | <p>Missing interconnecting air hoses, with appropriate connectors, which connect compressor to dental operating and treatment unit.</p> <p>Damaged or deteriorated components prevent operation of the unit.</p> <p>Intake filter elements are unserviceable.</p> |
| 2 | B, D | <p>Operational Checkout</p> <p>a. Observe the pressure gauge.</p> <p>b. Observe that the unloader valve switches and compressor vents to atmosphere.</p> <p>c. Observe that when pressure decreases to 60 psi the unloader valve switches back and compressor pumps for approximately 8 seconds to reach 80 psi.</p> <p>d. Observe that the cycle repeats.</p> <p>e. Verify that the dryness indicator is blue.</p> <p>f. Rotate the four transit cover supports. Place transit case cover on supports.</p> | <p>The unit fails to operate.</p> <p>Pressure does not increase to 80psi in approximately 40 seconds.</p> <p>The unloader valve does not switch or the pressure does not decrease to 60 psi.</p> <p>The unloader valve does not switch or the pressure does not reach 80psi.</p> <p>The cycle does not repeat.</p> <p>The dryness indicator is other than blue.</p> |
| 3 | B, D | <p>Air Storage Tank</p> <p>a. Inspect air tank for leaks, damage, or excessive rust.</p> <p>b. Inspect hoses and ensure that the hose(s) can be properly connected.</p> <p>c. Ensure pressure relief/drain valve opens and closes properly.</p> | <p>Air tank leaks or damage or rust accumulation precludes operation.</p> <p>The hose(s) cannot be connected to the storage tank.</p> <p>The valve cannot be opened or it leaks when closed.</p> |
| 4 | B, D, A | <p>Pressure Gauge</p> <p>Check for dents, a cracked or broken dial cover, or gauge indications beyond the normal range.</p> | <p>The pressure gauge does not function.</p> |

(continued) Appendix A. Operator PMCS

6520-01-398-4613
Compressor Dehydrator, Dental, Model PAC 6.7

[B-Before Operation, D-During Operation, A-After Operation, Q-Quarterly, and S-Semiannually]

| ITEM NO | INTERVAL | ITEM TO BE INSPECTED AND PROCEDURE | IS NOT MISSION CAPABLE IF: |
|---------|----------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------|
| 5 | B, D, A | <p>Dryness Indicator</p> <p>a. Inspect for dents, a cracked or missing indicator cover, or the lack of any color indication.</p> <p>b. Ensure that the indicator is blue.</p> | <p>The damaged indicator prevents operation of the unit.</p> <p>The dryness indicator is other than blue.</p> |
| 6 | S | <p>Case</p> <p>a. Inspect the case for signs of excessive wear.</p> <p>b. Check the air relief valve.</p> | |

6525-01-312-6411

X-Ray Apparatus, Radiographic/Fluoroscopic, Model CS-8952

[B-Before Operation, D-During Operation, A-After Operation, Q-Quarterly, and S-Semiannually]

| ITEM NO | INTERVAL | ITEM TO BE INSPECTED AND PROCEDURE | IS NOT MISSION CAPABLE IF: |
|---------|----------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1 | B, A | <p>X-Ray Apparatus</p> <p>a. Conduct an inventory to ensure that the items listed on the Equipment Parts and Accessories List are on hand.</p> <p>b. Inspect unit for physical damage, rust, or excessively worn components.</p> <p>c. Verify that the Verification/Certification sticker (DD Form 2163) has a current date.</p> | <p>Missing components prevent the use of the unit.</p> <p>Unserviceable components prevent the use of the unit.</p> <p>The x-ray apparatus has not been verified within the last 12 months.</p> |
| 2 | B, A | <p>X-Ray Operational Test</p> <p>a. Perform daily pre-operational system check as directed by manufacture's literature.</p> <p>NOTE: Ensure that personal protective apron, lead blockers, and suitable radiation protection measures are taken.</p> <p>(1) Turn power on and adjust line set as needed.</p> <p>(2) Perform table check.</p> <p>(a) Press and hold the longitudinal switch on spot film device (SFD) until the tabletop reaches its limit of travel.</p> <p>(b) Press and hold the table longitudinal foot switch until the tabletop reaches its limit of travel.</p> <p>(c) Press and hold the table center switch until the tabletop stops.</p> <p>(d) Press and hold the Trendelenburg tilt switch until the table reaches its maximum tilt and stops.</p> <p>(e) Press and hold the vertical tilt switch. The table should stop at the horizontal position. Release the switch, and press and hold the switch again. The table should rotate to its maximum tilt of 88 degrees, proving the tabletop is on "center."</p> <p>(f) Press and hold the Trendelenburg tilt switch until the table stops at horizontal. Release the switch.</p> | <p>The line adjustment cannot be accomplished.</p> <p>There are any malfunctions or unusual noises.</p> <p>The tabletop does not move approximately 30" from its center position before it stops.</p> <p>The tabletop does not move approximately 30" from its center position before it stops.</p> <p>The tabletop does not move to its center position from either of the above mentioned longitudinal positions, before stopping.</p> <p>The table does not reach its maximum 12 degrees before stopping.</p> <p>The table does not reach center or if it does not rotate to 88 degrees.</p> <p>The table does not reach horizontal position.</p> |

(continued) Appendix A. Operator PMCS

6525-01-312-6411

X-Ray Apparatus, Radiographic/Fluoroscopic, Model CS-8952

[B-Before Operation, D-During Operation, A-After Operation, Q-Quarterly, and S-Semiannually]

| ITEM NO | INTERVAL | ITEM TO BE INSPECTED AND PROCEDURE | IS NOT MISSION CAPABLE IF: |
|---------|----------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| | | <p>(g) On the spot film device, disengage the carriage locks and the compression locks and move the spot film device in all directions.</p> <p>(3) Perform the tube stand check</p> <p>(a) Verify that the tube stand is energized by operating the locks and moving it through its various motions.</p> <p>(b) Check the collimator to ensure that all lamps will light.</p> <p>(4) Perform the spot film device (SFD) check.</p> <p>(a) Observe the spot film device for the presence of power. All push buttons should be lit.</p> <p>(b) Insert an empty 9" x 9" cassette into the SFD tunnel. Cycle the cassette carriage by pressing the PARK/LOAD switch. The carriage should alternate between its park and load positions.</p> <p>(c) Verify that various pictorial representations can be set on the display (i.e., 2 on 1, 3 on 1, and 9 on 1).</p> <p>(5) Perform the warm-up procedure.</p> <p>NOTE: Always perform the warm-up procedure no more than one hour before the first case of the day or if the system has been idle for one hour or longer.</p> <p>(a) Warm up the over-table tube.</p> <p>[1] Disable autotiming and close the collimator blades. Select 70 kVp, 100 mA, 1.0 second.</p> <p>[2] Warm up the over-table x-ray tube by making four (4) exposures at 15-second intervals.</p> <p>(b) Make a fluoroscopic exposure by performing the following steps:</p> <p>[1] Press the 200L SPOT push button switch on the generator front panel.</p> <p>[2] On the fluoroscopic controls section of the generator panel, select mA station B and rotate the "minutes" dial to the 5 (minute) position.</p> | <p>The device requires more than 20 pounds of force to move it.</p> <p>The locks do not work or if the tube stand cannot be moved into various positions.</p> <p>All lamps do not energize.</p> <p>The buttons are not lit.</p> <p>The carriage does not alternate between park and load positions.</p> <p>The display does not indicate the correct selection or the cassette is not motor powered into the correct position.</p> <p>The selections cannot be made.</p> <p>The unit will not make the exposures.</p> <p>Fluoroscopic exposures cannot be made.</p> |

6525-01-312-6411

X-Ray Apparatus, Radiographic/Fluoroscopic, Model CS-8952

[B-Before Operation, D-During Operation, A-After Operation, Q-Quarterly, and S-Semiannually]

| ITEM NO | INTERVAL | ITEM TO BE INSPECTED AND PROCEDURE | IS NOT MISSION CAPABLE IF: |
|---------|----------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| | | <p>[3] Rotate the fluoroscopic kVp control until 70 kV is indicated on the fluoroscopic kVp meter. Ensure that the spot film device cassette carriage is in the park position.</p> <p>[4] Place a suitable fluoroscopic kVp phantom on the tabletop in the in-beam position.</p> <p>[5] Depress either the footswitch or x-ray push button on the spot film device.</p> <p>[6] Observe the imaging system mirror. A sharp, clear x-ray image of the grid chamber mechanism should be displayed.</p> <p>NOTE: Under-table (UT) shutters must always be visible and mechanically coned down as necessary.</p> <p>[7] Place a 9" x 9" cassette into the SFD (this should activate the system to make radiographic exposure). Locate the footswitch behind the operator barrier. Select an under table (UT) exposure of 70 kVp, 0.1 second. Depress footswitch, make fluoro exposure. From SFD location, make radiographic exposure.</p> <p>[8] Repeat above procedure with "Autotiming" set "ON." Select "Table" and "Normal density." Set radiographic exposure to about 50% more time than expected.</p> <p>NOTE: Phototiming failure does not deadline the system, but does reduce overall capability.</p> <p>b. Clean x-ray unit as directed by the manufacturer's literature.</p> | <p>The unit does not produce a clear image.</p> <p>The system will not transition from "fluoro" imaging to radiographic mode, with actual radiographic exposure.</p> <p>The system will not transition from "fluoro" mode to radiographic mode with exposure.</p> |

(continued) Appendix A. Operator PMCS

6525-01-325-3740

Portable X-Ray System, Model 1200

[B-Before Operation, D-During Operation, A-After Operation, Q-Quarterly, and S-Semiannually]

| ITEM NO | INTERVAL | ITEM TO BE INSPECTED AND PROCEDURE | IS NOT MISSION CAPABLE IF: |
|---------|----------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1 | B | <p>X-Ray System</p> <p>a. Conduct an inventory to ensure that the items listed on the Equipment Parts and Accessories List are on hand.</p> <p>b. Inspect unit for damage, discoloration, or excessively worn components.</p> <p>c. Verify that the Medical Equipment Verification/Certification sticker (DD Form 2163) has a current date.</p> | <p>Missing components prevent the use of the x-ray apparatus.</p> <p>Unserviceable components prevent the use of the x-ray apparatus.</p> <p>The unit has not been verified within the last 12 months.</p> |
| 2 | B | <p>X-Ray Operational Test</p> <p>a. Perform the "Assembly/Setup Procedure" as directed by the manufacturer's literature.</p> <p>(1) Remove reusable storage container from the wooden shipping crate, release leg clips.</p> <p>(2) Open the reusable container.</p> <p>(3) Remove the stand frame assembly, position on floor, engage rear wheel brakes, fold out legs and insert locking pins to frame/leg holes to lock legs.</p> <p>(4) Remove the pipe assembly, lower section, and position locking handles down (to the horizontal unlocked position).</p> <p>(5) Position the pipe assembly, lower section, with the gear rack toward the rear of the stand. Align the four "T" head bolts on the bottom of the pipe assembly, lower section with the four key-slots on the stand frame assembly. Lower into place, being sure the "T" bolts fit into the key-slots.</p> <p>(6) Slide the pipe assembly, lower section, forward (approximately 1 inch) and lift, the two locking handles up (to the vertical locking position). Ensure that both locking clips fit into locking clip slots.</p> <p>(7) Remove pipe assembly, upper section, from the reusable container.</p> <p>(8) Position pipe assembly, upper section, locking handle to the up (unlocked) position.</p> | <p>The unit cannot be set up.</p> |

6525-01-325-3740

Portable X-Ray System, Model 1200

[B-Before Operation, D-During Operation, A-After Operation, Q-Quarterly, and S-Semiannually]

| ITEM NO | INTERVAL | ITEM TO BE INSPECTED AND PROCEDURE | IS NOT MISSION CAPABLE IF: |
|---------|----------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------|
| | | <p>(9) Position the pipe assembly, upper section, on top of the pipe assembly, lower section, with the rack gear facing the rear of the stand.</p> <p>(10) Pull locking handle down (from the unlocked position) to the perpendicular position and rotate 90 degrees to the lock position, secure locking handle with spring clip.</p> <p>(11) Crank gearbox assembly up to a comfortable working height.</p> <p>(12) Remove cross arm assembly from side of stand frame assembly.</p> <p>(13) Press cross arm horizontal travel release brake and slide cross arm into gearbox assembly.</p> <p>(14) Lift x-ray generator assembly out of reusable container, remove safety pin, position x-ray generator yoke assembly into end of cross arm assembly, secure safety pin.</p> <p>(15) Lift control assembly out of reusable container. Position the control arm assembly on the stand bracket pull out on the end clips, and snap in to place.</p> <p>(16) Attach line cord to control assembly "LINE IN" connector, attach exposure switch cable to control assembly, "HAND SWITCH" connector and connect one end of the interconnect cable to the control assembly "LINE OUT" connector and the remaining end to the x-ray generator assembly connector.</p> <p>b. Verify the "Assembly Check Out" procedure as directed by the manufacturer's literature.</p> <p>(1) Verify that the stand foldout leg locks pins are installed.</p> <p>(2) Verify that the pipe assembly lower section locking handles are in the up position and that the handle locking clips are engaged.</p> <p>(3) Verify that pipe assembly upper section locking handle is in the locked position and the spring clip is engaged.</p> <p>(4) Verify that the x-ray generator safety pin is installed and locked.</p> | <p>The assembly cannot be accomplished.</p> |

(continued) Appendix A. Operator PMCS

6525-01-325-3740

Portable X-Ray System, Model 1200

[B-Before Operation, D-During Operation, A-After Operation, Q-Quarterly, and S-Semiannually]

| ITEM NO | INTERVAL | ITEM TO BE INSPECTED AND PROCEDURE | IS NOT MISSION CAPABLE IF: |
|---------|----------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------|
| 3 | B, A | <p>(5) Verify that the line cord, the exposure switch and the interconnect cable are properly installed.</p> <p>c. Perform pre-operational checkout procedures as directed by the manufacturer's literature.</p> <p>(1) Check power cord, and all interconnecting cables.</p> <p>(2) Verify that the 50/60 Hz switch is set correctly.</p> <p>(3) Turn on power switch; verify the correct line set.</p> <p>(4) Close the collimator shutters.</p> <p>(5) Select the 60kVp/40mA station.</p> <p>(6) Set timer for 0.1 seconds.</p> <p>(7) Step back from the unit with the exposure switch.</p> <p>(8) Press the exposure switch to the first position; verify the ready lamp goes off and on after about a one second delay.</p> <p>(9) Press for second trigger position; verify the x-ray on lamp and audio tone operate.</p> <p>Periodic Maintenance</p> <p>a. Perform Operator Maintenance as directed by manufacture's literature.</p> <p>b. Inspect and clean the unit as directed by the manufacturer's literature.</p> | <p>The checkout cannot be accomplished.</p> <p>The scheduled maintenance cannot be completed.</p> <p>The unit is unsafe or hazardous.</p> |

6525-01-370-7552

Portable Dental X-Ray System, Model ALPHA MPDX

[B-Before Operation, D-During Operation, A-After Operation, Q-Quarterly, and S-Semiannually]

| ITEM NO | INTERVAL | ITEM TO BE INSPECTED AND PROCEDURE | IS NOT MISSION CAPABLE IF: |
|---------|----------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1 | B | <p>X-Ray System</p> <p>a. Conduct an inventory to ensure that the items listed on the Equipment Parts and Accessories List are on hand.</p> <p>b. Inspect components for damage, discoloration, or excessively worn components.</p> <p>c. Verify the date on the Medical Equipment Verification/Certification sticker (DD Form 2163) is current.</p> | <p>Missing components or accessories prevent the operation of the x-ray system.</p> <p>Unserviceable components prevent the use of x-ray.</p> <p>The unit has not been verified within the last 12 months.</p> |
| 2 | B, D, A | <p>X-Ray System Operational Check Out</p> <p>a. Perform the "Assembly/Setup Procedure" as directed by the manufacturer's literature.</p> <p>b. Inspection after assembly as directed by the manufacturer's literature.</p> <p>(1) Ensure that all quick release pins are fully inserted and locked in place.</p> <p>(2) Ensure that all locking knobs are hand-tight (full clockwise position).</p> <p>(3) Verify the security of the cone installed on the x-ray control assembly.</p> <p>(4) Check security of all electrical connectors.</p> <p>(5) Verify that all labels are securely affixed and legible.</p> <p>(6) Thoroughly inspect the assembled x-ray system for tight fittings, possible missing parts (including the operation and maintenance manuals), frayed electrical cords, cracks, chips, excessive wear, or other signs of deterioration.</p> <p>(7) Using a lint-free cloth, remove any noticeable dirt or excess dust from the assembled unit.</p> <p>(8) Check x-ray head subassembly 1A2A1 in all working positions for possible drift.</p> <p>(9) Check scissor arm assembly 1A3 in all working positions for ease of motion.</p> | <p>The unit cannot be assembled.</p> <p>Knobs cannot be tightened sufficiently to prevent drift or to keep unit from falling.</p> <p>Loose connectors prevent the operation of the x-ray system.</p> <p>Loose fittings, missing parts, or frayed cords prevent operation of the x-ray system.</p> <p>X-ray head drift prevents the operation of the x-ray system.</p> <p>Scissor arm assembly is unable to hold position prevents the operation of the x-ray system.</p> |

(continued) Appendix A. Operator PMCS

6525-01-370-7552

Portable Dental X-Ray System, Model ALPHA MPDX

[B-Before Operation, D-During Operation, A-After Operation, Q-Quarterly, and S-Semiannually]

| ITEM NO | INTERVAL | ITEM TO BE INSPECTED AND PROCEDURE | IS NOT MISSION CAPABLE IF: |
|---------|----------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| | | <p>c. Perform the pre-operational checkout procedure as directed by the manufacturer's literature.</p> <p>(1) Check the line power plug connection to the line power receptacle.</p> <p>(2) Check security of the electrical connection between x-ray control assembly 1A1 and x-ray unit.</p> <p>(3) Cover the cone port with lead shielding.</p> <p>(4) Position the tube head away from the x-ray unit mounting post (scissor arm assembly fully extended).</p> <p>(5) Set the rotary "TIMER" switch to 0.1 second.</p> <p>(6) While holding exposure switch, STEP BACK FROM THE UNIT APPROXIMATELY SIX (6) FEET.</p> <p>(7) Depress and hold down the exposure switch. The x-ray Indicator light will illuminate and the buzzer will emit as audible tone. The exposure switch will automatically switch off when the time set on the "TIMER" switch expires.</p> <p>d. Perform routine cleaning as directed by the manufacturer's literature.</p> | <p>Improper fit prevents the operation of the x-ray system.</p> <p>Improper fit prevents the operation of the x-ray system.</p> <p>The unit does not shut off.</p> |

6525-01-384-9296

X-Ray Apparatus, Model LCROKS

[B-Before Operation, D-During Operation, A-After Operation, Q-Quarterly, and S-Semiannually]

| ITEM NO | INTERVAL | ITEM TO BE INSPECTED AND PROCEDURE | IS NOT MISSION CAPABLE IF: |
|---------|----------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1 | B, D, A | <p>X-Ray Apparatus</p> <p>a. Conduct an inventory to ensure that the items listed on the Equipment Parts and Accessories List are on hand.</p> <p>b. Inspect unit for damage, discoloration, or excessively worn components.</p> <p>c. Verify that the Verification/Certification sticker (DD Form 2163) has a current date.</p> | <p>Missing components prevent the use of the x-ray.</p> <p>Unserviceable components prevent the use of x-ray.</p> <p>The unit has not been verified within the last 12 months.</p> |
| 2 | B | <p>X-Ray Operational Test</p> <p>Conduct "Operator Maintenance" as directed by manufacturer's literature.</p> <p>a. Check control panel for nicks, scratches, and/or dents.</p> <p>b. Ensure proper seating of APR labels.</p> <p>c. Inspect unit for all warning labels, serial tags, UL, and CSA tags.</p> | <p>The labels are missing, unreadable, or outdated.</p> |

(continued) Appendix A. Operator PMCS

6530-00-926-2151

Sterilizer, Surgical Dressing 16x36 in., Model M-138

[B-Before Operation, D-During Operation, A-After Operation, Q-Quarterly, and S-Semiannually]

| ITEM NO | INTERVAL | ITEM TO BE INSPECTED AND PROCEDURE | IS NOT MISSION CAPABLE IF: |
|---------|----------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1 | B, A | <p>Sterilizer</p> <p>a. Conduct an inventory to ensure that the items listed on the Equipment Parts and Accessories List are on hand.</p> <p>b. Inspect the unit for obvious signs of damage such as cracks, dents, leaks, or broken components.</p> | <p>The shelves are missing.</p> <p>Gasket is torn, sight-glass is broken, or legs cannot be locked into their supporting position.</p> |
| 2 | B | <p>Sterilizer Operational Check</p> <p>a. Ensure that the unit is set up and assembled properly as directed by the manufacturer's literature.</p> <p>b. Remove the chamber drain-plug and inspect for lint and sediment from the strainer.</p> <p>c. Inspect and clean the interior surfaces of the chamber, with mild detergent and water, before heating. Clean the shelves in the same manner.</p> <p>CAUTION: DO NOT USE STEEL WOOL OR ABRASIVE CLEANERS.</p> <p>d. Inspect door and door gasket</p> <p>e. Inspect sight glass for mineral deposits.</p> <p>f. Inspect fill washer.</p> | <p>The sterilizer cannot be assembled properly.</p> <p>Built-up sediment cannot be removed and prevents the chamber from draining.</p> <p>Chamber does not hold pressure.</p> <p>Door does not seal.</p> <p>Sight glass is broken or mineral deposits obscure water level in sight glass.</p> <p>Fill washer is missing.</p> |
| 3 | B, D | <p>Electrical Operations</p> <p>Ensure that the frame of the sterilizer is adequately grounded before operating on electrical power as directed by the manufacturer's literature. Seek assistance from unit medical maintenance section if necessary.</p> | <p>Unit is not grounded.</p> |
| 4 | B | <p>Sterilizer Jacket</p> <p>a. Turn operating valve to sterilize. This opens an escape route for trapped air.</p> <p>b. Open drain valve to allow for a lower air escape route.</p> <p>c. Fill jacket with water to about ½ mark.</p> <p>d. Close drain valve when water flows freely without burping.</p> | <p>Jacket cannot be filled with water.</p> |

6530-00-926-2151

Sterilizer, Surgical Dressing 16x36 in., Model M-138

[B-Before Operation, D-During Operation, A-After Operation, Q-Quarterly, and S-Semiannually]

| ITEM NO | INTERVAL | ITEM TO BE INSPECTED AND PROCEDURE | IS NOT MISSION CAPABLE IF: |
|---------|----------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| | | <p>e. Ensure that the water level viewed in the sight glass is at least at the ¼ mark as directed by the manufacturer's literature.</p> <p>WARNING: LIFT THE RELIEF HANDLE OF THE SAFETY VALVE OR TURN OPERATING VALVE TO THE DRY POSITION TO RELEASE ANY PRESSURE IN THE JACKET BEFORE REMOVING THE PLUG FROM THE FILLING FUNNEL. FILL THE STERILIZER JACKET WITH THE PUREST WATER AVAILABLE AND INSPECT FOR WATER LEAKS. ENSURE THE WATER IN THE SIGHT GLASS IS AT LEAST AT THE ¼ MARK.</p> <p>f. Verify operation of the pressure control switch knob. Turn the pressure control knob to the maximum clockwise position.</p> <p>g. Verify operation of the operating valve. Ensure that the operating valve is in the OFF position.</p> <p>h. Turn the heat switch on and verify that the red pilot light is glowing.</p> <p>i. Turn pressure valve fully counterclockwise to open the low-pressure relief valve.</p> <p>j. When pressure reaches 18 – 20 psi, the low-pressure valve should release pressure.</p> <p>k. Turn pressure relief valve fully clockwise to take the low-pressure valve out of the system.</p> <p>l. Verify the increase in pressure and test the safety valve by depressing the safety lever.</p> <p>m. When pressure reaches 27 – 32 psi, the high-pressure valve should release pressure.</p> <p>n. Verify that the pressure gauge indicates the desired pressure of 18 psi for 250 degrees F or 29 psi for 270 degrees F.</p> <p>o. Turn the pressure control switch knob slowly counterclockwise until the pilot light goes out. Verify that the pressure control cycles and maintains the selected pressure.</p> | <p>Jacket cannot be filled with water.</p> <p>Jacket leaks or cannot be filled with water.</p> <p>Pressure control switch does not operate.</p> <p>Operating valve does not function.</p> <p>Heating elements do not energize.</p> <p>Safety valve does not activate.</p> <p>Safety valve does not activate.</p> <p>Safety valve does not activate.</p> <p>Desired steam pressure cannot be reached or pressure gauge is faulty.</p> |

(continued) Appendix A. Operator PMCS

6530-00-926-2151

Sterilizer, Surgical Dressing 16x36 in., Model M-138

[B-Before Operation, D-During Operation, A-After Operation, Q-Quarterly, and S-Semiannually]

| ITEM NO | INTERVAL | ITEM TO BE INSPECTED AND PROCEDURE | IS NOT MISSION CAPABLE IF: |
|---------|----------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| | | <p>NOTE: A pre-heat period of 10 to 15 minutes is recommended to allow the pressure to stabilize. There are no markings or calibration on the pressure control switch since temperature is a function of absolute pressure rather than gauge pressure. Depending on altitude and atmospheric conditions, reaching 270° F may require between 27 and 32 psi gauge pressure. The pressure switch must be adjusted to the pressure, which will give the desired temperature.</p> <p>p. Load the sterilizer and verify proper operation.</p> <p>CAUTION: IN THE EVENT THAT WATER IN THE JACKET RUNS LOW, THE WATER CUT-OFF WILL INTERRUPT THE POWER SUPPLY TO THE HEATERS. IF THIS OCCURS, LIFT THE RELIEF HANDLE ON THE SAFETY VALVE TO RELEASE ANY PRESSURE IN THE JACKET BEFORE REMOVING PLUG FROM FILLING FUNNEL. WAIT UNTIL INTERNAL PARTS COOL BELOW THE BOILING POINT AND REFILL THE JACKET WITH WATER AND PRESS THE RESET BUTTON (LOCATED UNDER THE HEATER BOX). PROCEED WITH THE REGULAR OPERATING CYCLE FROM THE BEGINNING.</p> <p>q. Close the chamber door.</p> <p>r. Turn the operating valve to sterilize.</p> <p>s. Let the chamber pressurize.</p> <p>t. Check for leaks. The steam trap may stick open (rap with a solid object to release it).</p> <p>u. Set the timer.</p> <p>v. Check that the pilot light cycles on and off.</p> <p>w. Check that the chamber maintains pressure.</p> <p>x. When the timer goes off, turn the operating valve to "DRY."</p> <p>y. Check that the pressure goes to about -5 psi for about 15 minutes before the pressure releases and the door can be opened.</p> | <p>Pressure control does not operate.</p> <p>Door does not seal.</p> <p>Chamber leaks or trap fails to close.</p> <p>Desired pressure cannot be maintained.</p> <p>Sterilizer chamber does not release pressure.</p> <p>Sterilizer chamber does not pull a vacuum.</p> |

6530-00-926-2151

Sterilizer, Surgical Dressing 16x36 in., Model M-138

[B-Before Operation, D-During Operation, A-After Operation, Q-Quarterly, and S-Semiannually]

| ITEM NO | INTERVAL | ITEM TO BE INSPECTED AND PROCEDURE | IS NOT MISSION CAPABLE IF: |
|---------|----------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------|
| 5 | B, D | Gasoline Heat Note: No longer authorized for Department of the Defense use. | |
| 6 | B, D | Direct Steam Operation a. Conduct direct steam operation as directed by TM 8-6530-004-24&P. b. Load the sterilizer and verify proper operation. WARNING: TO PREVENT POSSIBLE INJURY TO PERSONNEL RESULTING FROM BURSTING BOTTLES AND HOT FLUID, USE ONLY BOROSILICATE (PYREX) FLASKS WITH VENTED CLOSURES FOR STERILIZING LIQUIDS. | Sterilizer does not operate. |

Appendix B. Repairer PMCS

Generic Standards

[M-Monthly, Q-Quarterly, S-Semiannually, and A-Annually]

| ITEM NO | INTERVAL | ITEM TO BE INSPECTED AND PROCEDURE | IS NOT MISSION CAPABLE IF: |
|---------|----------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------|
| 1 | A | Ensure that all ancillary components necessary to operate the equipment or system are on hand. | Ancillary components are missing. |
| 2 | A | Ensure that all components and accessories issued with the equipment or system are on hand. | Components or accessories are not readily available. |
| 3 | A | Ensure that all TMDE required to perform CVC and PMCS are on hand and calibrated. | TMDE is not available. |
| 4 | A | Inspect for corrosion, rust, physically damaged parts, deteriorated materials, and damage to protective coatings. | Rust on outer surface parts determined by the Infection Control Nurse to be a health hazard. |
| 5 | A | Ensure the operator and maintenance manuals or documentation are on hand. Identify the location of such material if it is not packed with the equipment. | Operator and maintenance manuals are not readily available. |
| 6 | A | Verify that the equipment or system has no broken parts or accessories, i.e., switches, knobs, casters, plastic coverings, hoses, casings, etc. | Equipment is not functional due to broken parts. |
| 7 | A | Ensure that fluid levels, lubricants, physical limits or settings for operation are correct. | Levels are below those established in the TM or manufacturer's literature. |
| 8 | A | During prolonged exercises or missions involving patient treatment, scheduled testing of electrically operated medical equipment designated for use in critical care areas will be performed. | Equipment fails the electrical safety test. |
| 9 | A | Verify operation of the equipment or system in accordance with published TMs and the manufacturer's literature. | Equipment does not function according to the TM or manufacturer's literature. |
| 10 | A | Perform CVC and PMCS as necessary indicating compliance with standards. Place appropriate labels on equipment. | Equipment cannot be calibrated to TM or manufacturer's specifications. |
| 11 | A | Inspect for unusual operation, noises, leakage, or other unexpected results. | Noticeable fluid leaks or unexpected noises are detected. |
| 12 | A | Shut down equipment, and clean and dry parts or components that were subjected to liquid contact. Use of compressed air and disassembly of components to remove liquid or reagent materials may be necessary. | Unit or components are not clean or dry. |
| 13 | A | Check the electrical power cord for cuts, fraying, or deterioration. | Electrical plug is missing a pin/blade or the cord insulation is cut through the outer coating. |
| 14 | A | Ensure that alarms and visual indicators are functioning properly. | Alarms and indicators are not functioning properly. |
| 15 | A | Verify proper battery condition. | Battery will not charge or is visibly defective (when applicable). |

6520-01-272-4531
 Dental Operating Unit, ADEC Model 3406 Porta-Cart

[M-Monthly, Q-Quarterly, S-Semiannually, and A-Annually]

| ITEM NO | INTERVAL | ITEM TO BE INSPECTED AND PROCEDURE | IS NOT MISSION CAPABLE IF: |
|---------|----------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1 | S | <p>Dental Unit</p> <p>a. Conduct an inventory to ensure that the items listed in the Equipment Parts or Accessories List are on hand.</p> <p>b. Inspect components for damage, discoloration, or excessively worn components.</p> | <p>Missing components or accessories prevent the operation of the dental unit.</p> <p>Unserviceable components prevent the use of the dental unit.</p> |
| 2 | S | <p>Operational Checks</p> <p>a. Review the general service information as provided in the manufacturer's literature.</p> <p>b. Check the air and water filters as directed in the manufacturer's literature.</p> <p>c. Check the air and water regulator as directed in the manufacturer's literature.</p> <p>d. Verify the operation of the "Century II Control System" as directed in the manufacturer's literature.</p> <p>e. Verify the operation of the three-way micro valves as directed in the manufacturer's literature.</p> <p>f. Verify the operation of the foot control valve as directed in the manufacturer's literature.</p> <p>g. Verify the operation of the signal relay valve as directed in the manufacturer's literature.</p> <p>h. Verify the operation of the chip blower valve as directed in the manufacturer's literature.</p> <p>i. Verify the operation of the three-way toggle valve as directed in the manufacturer's literature.</p> <p>j. Verify the operation of the needle valves as directed in the manufacturer's literature.</p> <p>k. Verify the operation of the syringe as directed in the manufacturer's literature.</p> | <p>The air pressure drops more than 15 psi or the water pressure drops more than 10 psi.</p> <p>The air regulator does not regulate between 60 psi to 80 psi or the water regulator does not regulate between 30 psi to 40 psi.</p> <p>There are air or water leaks that prevent the use of the dental unit.</p> <p>The three-way micro valves do not control the flow of coolant air or coolant water.</p> <p>The foot control valve does not operate the handpieces.</p> <p>The signal relay valve does not initiate the coolant air or coolant water.</p> <p>The chip blower valve does not provide chip-air flow to the handpieces.</p> <p>The three-way toggle valve does not pressurize or de-pressurize to water tank.</p> <p>The syringe leaks air or water or does not pass air or water.</p> |

(continued) Appendix B. Repairer PMCS

6520-01-272-4531
 Dental Operating Unit, ADEC Model 3406 Porta-Cart

[M-Monthly, Q-Quarterly, S-Semiannually, and A-Annually]

| ITEM NO | INTERVAL | ITEM TO BE INSPECTED AND PROCEDURE | IS NOT MISSION CAPABLE IF: |
|---------|----------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------|
| 3 | S | l. Verify the operation of the air vacuum system as directed in the manufacturer's literature. m. Verify the operation of the air saliva ejector as directed in the manufacturer's literature. Storage Case Inspect the storage case for cracks, dents, or broken latches. | The air vacuum system does not provide vacuum. The air saliva ejector does not provide vacuum. |

6520-01-398-4613
Compressor Dehydrator, Dental, Model PAC 6.7

[M-Monthly, Q-Quarterly, S-Semiannually, and A-Annually]

| ITEM NO | INTERVAL | ITEM TO BE INSPECTED AND PROCEDURE | IS NOT MISSION CAPABLE IF: |
|---------|----------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1 | A | <p>Compressor Dehydrator</p> <p>a. Conduct an inventory to ensure that the items listed in the Equipment Parts or Accessories List are on hand.</p> <p>b. Inspect and verify that the compressor-dehydrator operates as directed by the "Operational Checkout" procedures in the Operator Preventative Maintenance Checks and Services.</p> <p>c. Verify electrical safety.</p> | <p>Missing interconnecting air hoses, with appropriate connectors, which connect compressor to dental operating and treatment unit.</p> <p>The unit does not operate as directed by the operational checkout procedures.</p> <p>The compressor-dehydrator fails any of the safety tests.</p> |
| 2 | A | <p>Air Storage Tank</p> <p>a. Inspect air tank for leaks, damage, or excessive rust.</p> <p>b. Inspect hoses and ensure that the hose(s) can be properly connected.</p> <p>c. Ensure pressure relief/drain valve opens and closes properly.</p> | <p>Air tank leaks or damage or rust accumulation precludes operation.</p> <p>The hose(s) cannot be connected to the storage tank.</p> <p>The valve cannot be opened or it leaks when closed.</p> |
| 3 | A | <p>Pressure Gauge</p> <p>Check for dents, a cracked or broken dial cover, or gauge indications beyond the normal range.</p> | <p>The pressure gauge does not function.</p> |
| 4 | A | <p>Dryness Indicator</p> <p>a. Inspect for dents, a cracked or missing indicator cover, or the lack of any color indication.</p> <p>b. Ensure that the indicator is blue.</p> | <p>The damaged indicator is unserviceable.</p> <p>The dryness indicator is other than blue.</p> |
| 5 | A | <p>Case</p> <p>a. Inspect the case for signs of excessive wear.</p> <p>b. Check the air relief valve.</p> | |

(continued) Appendix B. Repairer PMCS

6525-01-312-6411

X-Ray Apparatus, Radiographic/Fluoroscopic, Model CS-8952

[M-Monthly, Q-Quarterly, S-Semiannually, and A-Annually]

| ITEM NO | INTERVAL | ITEM TO BE INSPECTED AND PROCEDURE | IS NOT MISSION CAPABLE IF: |
|---------|--------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1 | S | <p>X-Ray Apparatus</p> <p>a. Conduct an inventory to ensure that the items listed in the Equipment Parts and Accessories List are on hand.</p> <p>b. Unpack and install as directed by manufacturer’s literature.</p> <p>c. Ensure retrofit kit (consists of heavy steel brackets under each end of table) is installed for possible shipment.</p> <p>d. Inspect unit for damage, excessive rust to critical parts, bearing tracks and races, etc., or excessively worn components.</p> | <p>Missing components prevent the use of the x-ray unit.</p> <p>The unit cannot be installed.</p> <p>The unit is unable to deploy.</p> <p>The unserviceable components prevent the use of the unit.</p> |
| 2 | S A | <p>X-Ray Operational Test</p> <p>a. Ensure each component is operational as directed by the manufacturer’s literature.</p> <p>b. Ensure daily pre-operational systems checks were performed as directed by manufacturer’s literature.</p> <p>c. Verify the pre-calibration checks as directed by manufacturer’s literature.</p> <p>d. Verify calibration before attempting the calibration procedures.</p> <p>NOTE: Perform manufacturer's calibration procedures <u>ONLY</u> if x-ray apparatus does not meet manufacturer’s specifications.</p> <p>WARNING: FOLLOW X-RAY TUBE WARM UP PROCEDURE AS DIRECTED BY MANUFACTURER’S LITERATURE.</p> <p>e. Calibrate the unit as directed by the manufacturer’s literature.</p> <p>(1) Calibrate the generator as directed by manufacturer’s literature.</p> <p>(2) Calibrate the spot film device as directed by manufacturer’s literature.</p> <p>(3) Calibrate the under-table collimator as directed by manufacturer's literature.</p> | <p>Components not operational prevent the use of the x-ray unit.</p> <p>The unit is not prepared for calibration.</p> <p>The unit is in need of calibration.</p> |

(continued) Appendix B. Repairer PMCS

6525-01-312-6411

X-Ray Apparatus, Radiographic/Fluoroscopic, Model CS-8952

[M-Monthly, Q-Quarterly, S-Semiannually, and A-Annually]

| ITEM NO | INTERVAL | ITEM TO BE INSPECTED AND PROCEDURE | IS NOT MISSION CAPABLE IF: |
|---------|----------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------|
| | | <p>(4) Calibrate the over-table collimator as directed by manufacturer's literature.</p> <p>(5) Calibrate the automatic exposure control as directed by manufacturer's literature.</p> <p>(6) Verify the image intensifier as directed by manufacturer's literature.</p> <p>f. Update the Medical Equipment Verification / Certification sticker (DD Form 2163).</p> | <p>The unit has not been verified or calibrated within the last 12 months.</p> |

(continued) Appendix B. Repairer PMCS

6525-01-325-3740
Portable X-Ray System, Model 1200

[M-Monthly, Q-Quarterly, S-Semiannually, and A-Annually]

| ITEM NO | INTERVAL | ITEM TO BE INSPECTED AND PROCEDURE | IS NOT MISSION CAPABLE IF: |
|---------|---------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1 | S | <p>X-Ray System</p> <p>a. Conduct an inventory to ensure that the items listed on the Equipment Parts and Accessories List are on hand.</p> <p>b. Inspect unit for damage, discoloration, or excessively worn components.</p> <p>c. Verify assembly of unit as directed by the manufacturer's literature.</p> <p>d. Verify the electrical safety.</p> | <p>Missing components prevent the use of the X-Ray.</p> <p>Unserviceable components prevent the use of x-ray.</p> <p>The unit cannot be assembled.</p> <p>The x-ray system fails any of the electrical safety tests.</p> |
| 2 | M, Q M Q Q Q Q | <p>Periodic Maintenance</p> <p>Perform the "Periodic Maintenance Schedule and Procedure" as directed by manufacturer's literature.</p> <p>a. Clean the unit.</p> <p>b. Visually inspect unit; check electrical cables and connectors for bent, broken, or loose pins, cracked or broken insulators, weak, broken or loose pin connections, dirt, and corrosion; repair as required.</p> <p>c. Verify that unit meets all of the pre-operational check out procedures.</p> <p>d. Tighten any loose hardware.</p> <p>e. Touch up paint, any scratches, chips or exposed metal.</p> | <p>The maintenance cannot be completed.</p> |
| 3 | S | <p>Alignment, Adjustment, Calibration and Checkout Procedures</p> <p>a. Perform the "Alignment, Adjustment, Calibration and Checkout" procedures as directed by the manufacturer's literature:</p> <p>(1) Line Voltage</p> <p>(2) Line Set</p> <p>(3) Calibration Set-Up</p> <p>(4) mA/kVp Calibration</p> <p>(5) Verify mA/kVp with 220 VAC/50Hz</p> <p>(6) Timer Test Data</p> <p>(7) Exposure Indication</p> <p>(8) Line Current</p> | <p>The unit cannot be calibrated or verified as directed.</p> |

(continued) Appendix B. Repairer PMCS

6525-01-325-3740
 Portable X-Ray System, Model 1200

[M-Monthly, Q-Quarterly, S-Semiannually, and A-Annually]

| ITEM NO | INTERVAL | ITEM TO BE INSPECTED AND PROCEDURE | IS NOT MISSION CAPABLE IF: |
|---------|----------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------|
| | | (9) mAs Meter (10) Reproducibility (11) Half Value Layer (12) Leakage Test (13) Light Luminance (14) Beam Alignment (15) Final Step b. Update the Medical Equipment Verification/Certification label (DD Form 2163). | The unit has not been verified within the last 12 months. |

(continued) Appendix B. Repairer PMCS

6525-01-370-7552
Portable Dental X-Ray System, Model ALPHA MPDX

[M-Monthly, Q-Quarterly, S-Semiannually, and A-Annually]

| ITEM NO | INTERVAL | ITEM TO BE INSPECTED AND PROCEDURE | IS NOT MISSION CAPABLE IF: |
|---------|----------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1 | S | <p>X-Ray System</p> <p>a. Verify that the items listed on the Equipment Parts and Accessories List are on hand.</p> <p>b. Unpack and assemble the x-ray unit as directed by manufacturer's literature.</p> | <p>Missing components or accessories prevent the operation of the dental unit.</p> <p>The unit cannot be assembled.</p> |
| 2 | Q | <p>Preventive Maintenance Schedule and Procedures</p> <p>a. Inspection/check procedures</p> <p>(1) Visually inspect the unit as directed by the manufacturer's literature.</p> <p>(2) Verify that the unit meets all of the pre-operational requirements according to the Operator Preventive Maintenance Checks and Services.</p> <p>(3) Check all hardware connections for security. Tighten any loose connections.</p> <p>(4) Inspect the unit for chips, scratches or exposed metal. Use touch-up paint to repair any damage to paint or finish.</p> <p>(5) Perform corrective, adjustment or calibration procedures as required to resolve a malfunction, or perform periodic alignment adjustment and calibration functions in accordance with the schedule provided in manufacturer's literature.</p> | <p>The check out cannot be accomplished.</p> |
| | Q | <p>b. Perform the cleaning procedures as directed by the manufacturer's literature.</p> | |
| | S | <p>c. Perform "Adjustment, Calibration and Test" as directed by the manufacturer's literature.</p> <p>(1) Hi-Pot Test</p> <p>(2) Leakage Current</p> <p>(3) Line Voltage Meter</p> <p>(4) mA/kVp Calibration</p> <p>(a) Calibration Set-up</p> <p>(b) Line Voltage</p> <p>(c) mA Calibration</p> <p>(d) kVp Calibration</p> <p>(5) Timer Test Data</p> | <p>The adjustments and calibration cannot be accomplished.</p> <p>Leakage or breakdown occurs at 1500V within 60 seconds.</p> <p>Leakage is more than 100 microamps.</p> <p>X-ray will not calibrate to 7mA +/-10%. X-ray will not calibrate to 70kVp +/-10%.</p> |

6525-01-370-7552
 Portable Dental X-Ray System, Model ALPHA MPDX

[M-Monthly, Q-Quarterly, S-Semiannually, and A-Annually]

| ITEM NO | INTERVAL | ITEM TO BE INSPECTED AND PROCEDURE | IS NOT MISSION CAPABLE IF: |
|---------|----------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| | S | (a) Calibration (b) Verification (6) Exposure Indication (7) Line Current (8) Half Value Layer (9) Reproducibility (10) Leakage Test (11) Beam Limiting Device (12) Final Step (13) Update the Medical Equipment Verification/Certification label (DD Form 2163). (14) Verify electrical safety. d. Perform long term storage maintenance procedures as directed by the manufacturer's literature. | Will not calibrate within +/-10% and +/-4ms. Any indicators prevent safe operation. The current is not less than 7Amps The results are not greater than 0.51. The results are not less than 0.02. Any reading exceeds 50mR. Tolerance is not within 5.8 – 6.2cm. The unit has not been verified within the last 12 months. The x-ray system fails any of the electrical safety tests. The unit cannot complete the degassing process. |

(continued) Appendix B. Repairer PMCS

6525-01-384-9296
X-Ray Apparatus, Model LCROKS

[M-Monthly, Q-Quarterly, S-Semiannually, and A-Annually]

| ITEM NO | INTERVAL | ITEM TO BE INSPECTED AND PROCEDURE | IS NOT MISSION CAPABLE IF: |
|---------|-------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1 | M | <p>X-Ray Apparatus</p> <p>a. Conduct an inventory to ensure that the items listed on the Equipment Parts and Accessories List are on hand.</p> <p>b. Inspect unit for damage, discoloration, or excessively worn components</p> | <p>Missing components prevent the use of the x-ray.</p> <p>Unserviceable components prevent the use of x-ray.</p> |
| 2 | <p>A</p> <p>S</p> | <p>X-Ray Operational Test</p> <p>NOTE: Install the unit as direct by manufacturer's literature. Prepare x-ray tube for radiographic use in accordance with the manufacturer's break-in instructions.</p> <p>NOTE: An unseasoned tube will not calibrate and may develop hot spots.</p> <p>a. Calibrate the unit as directed by the manufacturer's literature.</p> <p>b. Perform the maintenance schedule checks as directed by the manufacturer's literature</p> <p>(1) Perform external visual checks as directed by the manufacturer's literature.</p> <p>(a) Check control panel stand, if so equipped, for nicks, scratches, or dents.</p> <p>(b) Check for proper seating of APR labels.</p> <p>(c) Inspect unit for all warning labels, serial tags, UL and CSA tags.</p> <p>(2) Perform mechanical checks as directed by the manufacturer's literature.</p> <p>(a) Check mechanical operation of control panel on/off and prep/expose switches.</p> <p>(b) Remove H.T. cables from transformer ports and check for proper level of oil. Check that H.T. cables are securely tightened.</p> <p>(c) Check connections on all cables on top of H.T. transformer.</p> <p>(d) Check connections on all cables in electronics cabinet.</p> <p>(e) Check connections on all cables in operator control panel.</p> | <p>The unit cannot be installed.</p> <p>Unit cannot be calibrated.</p> <p>The labels are missing, unreadable, or out-dated.</p> <p>X-ray does not operate or an electrical hazard exists.</p> <p>Oil level is low or H.T. cables are not securely tightened.</p> <p>The cables are not secure.</p> <p>The cables are not secure.</p> <p>The cables are not secure.</p> |

6525-01-384-9296
X-Ray Apparatus, Model LCROKS

[M-Monthly, Q-Quarterly, S-Semiannually, and A-Annually]

| ITEM NO | INTERVAL | ITEM TO BE INSPECTED AND PROCEDURE | IS NOT MISSION CAPABLE IF: |
|---------|----------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------|
| | | <p>(3) Perform operational checks as directed by the manufacturer's literature.</p> <ul style="list-style-type: none"> (a) Check for power-up sequence. (b) Check for operation of control panel switches; run fault diagnostics. (c) Check for operation of control panel LEDs; run fault diagnostics. (d) Check for operation of control panel display; run fault diagnostics. (e) Check for operation of control panel to generator communications; run fault diagnostics. (f) Check +5V power supply. (g) Check +15V power supply. (h) Check +24V power supply. (i) Depress "PREP" switch and check that control panel display reads "READY." (j) Depress "EXPOSURE" switch; listen for audible indicator to sound and check control panel for exposure indicator light. (k) Check that "BUT" logic works – "BUT" LED should light. (l) Check for actual mAs indication in display. (m) Check that another AEC exposure cannot be made. (n) Check that the reset button resets the "BUT" and another exposure can be made. (o) Check kV, mA, and time accuracy. (p) Check PT station(s) for density. <p>(4) Regrease high tension cables as directed by manufacturer's literature</p> <p>(5) Replace NVRAM every 72 months as directed by the manufacturer's literature.</p> | <p>X-ray does not operate or an electrical hazard exists.</p> |

(continued) Appendix B. Repairer PMCS

6525-01-384-9296
 X-Ray Apparatus, Model LCROKS

[M-Monthly, Q-Quarterly, S-Semiannually, and A-Annually]

| ITEM NO | INTERVAL | ITEM TO BE INSPECTED AND PROCEDURE | IS NOT MISSION CAPABLE IF: |
|---------|----------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------|
| | A | <p>(6) Perform "Final Appearance Checks" as directed by the manufacturer's literature.</p> <p>(a) Clean all exposed exterior surfaces of the Clinix VP4 Generator.</p> <p>(b) Check that all mounting hardware is secure and all covers are in place.</p> <p>c. Update Medical Equipment Verification/Certification label (DD Form 2361)</p> | <p>The mounting hardware is not secured.</p> <p>The unit has not been verified within the last 12 months.</p> |

6530-00-926-2151
Sterilizer, Surgical Dressing 16x36 in., Model M-138

[M-Monthly, Q-Quarterly, S-Semiannually, and A-Annually]

| ITEM NO | INTERVAL | ITEM TO BE INSPECTED AND PROCEDURE | IS NOT MISSION CAPABLE IF: |
|---------|----------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1 | S | <p>Sterilizer</p> <p>a. Verify the components and accessories according to the Operator Preventive Maintenance Checks and Services.</p> <p>b. Inspect the unit for obvious signs of damage such as cracks, dents, leaks, or broken components.</p> | <p>The shelves are missing.</p> |
| 2 | S | <p>Sterilizer Operational</p> <p>a. Ensure that the unit is set up and assembled properly as directed by the Operator Preventive Maintenance Checks and Services.</p> <p>b. Ensure unit is wired per data plate diagram to conform to incoming power.</p> <p>c. Inspect door for proper operation. Ensure hinges are properly lubricated. Inspect door gasket for damage or deterioration.</p> <p>d. Inspect the case for damage. Ensure hinges and latches are properly lubricated.</p> | <p>Unit cannot be wired according to diagram.</p> <p>Sterilizer door does not close and seal.</p> <p>Damage prevents operation of the unit.</p> |
| 3 | S | <p>Sterilizer Jacket</p> <p>Verify operation of the sterilizer jacket according to the Operator Preventive Maintenance Checks and Services.</p> <p>WARNING: LIFT THE RELIEF HANDLE OF THE SAFETY VALVE OR TURN OPERATING VALVE TO THE DRY POSITION TO RELEASE ANY PRESSURE IN THE JACKET BEFORE REMOVING THE PLUG FROM THE FILLING FUNNEL. FILL THE STERILIZER JACKET WITH THE PUREST WATER AVAILABLE AND INSPECT FOR WATER LEAKS. INSPECT THE WATER LEVEL INDICATOR GAUGE AND ENSURE WATER IS AT LEAST AT ¼ MARK.</p> | <p>Jacket leaks or cannot be filled with water.</p> <p>Water level indicator gauge is broken or excessive mineral deposits obscure the reading of the water level.</p> |
| 4 | S | <p>Operation Valve</p> <p>a. Conduct operating valve checks.</p> <p>b. Verify the increase in pressure and test the safety valve by depressing the safety lever.</p> | <p>Operating valve leaks or does not operate properly.</p> <p>Pressure does not increase or if the safety valve does not release pressure when depressed.</p> |

Appendix C. Maintenance Allocation Chart

6520-01-272-4531
Dental Operating Unit, ADEC Model 3406 Porta-Cart

[C – Operator or Crew; O – Unit Maintenance; DS – Direct Support Maintenance; GS – General Support Maintenance; D – Depot Maintenance]

| (1) GROUP NUMBER | (2) ASSEMBLY GROUP | (3) MAINTENANCE FUNCTION | (4) MAINTENANCE LEVEL | (5) TOOLS AND EQUIPMENT | (6) REMARKS |
|------------------------|---------------------------|--------------------------------|-----------------------------|----------------------------------|----------------|
| 00 | Dental Unit | Inspect | C | | |
| | | Service | O | 01 | |
| | | Test | O | 01 | |
| | | Repair | O | 01 | |
| | | Overhaul | O | 01 | |
| 01 | Air and Water Filters | Inspect | O | 01 | |
| | | Replace | O | 01 | |
| 02 | Air and Water Regulators | Inspect | O | 01 | |
| | | Service | O | 01 | |
| | | Test | O | 01 | |
| | | Repair | O | 01 | |
| | | Overhaul | O | 01 | |
| 03 | Century II Control System | Inspect | O | 01 | |
| | | Service | O | 01 | |
| | | Test | O | 01 | |
| | | Repair | O | 01 | |
| | | Overhaul | O | 01 | |
| 04 | Three-Way Micro Valves | Inspect | O | 01 | |
| | | Service | O | 01 | |
| | | Test | O | 01 | |
| | | Repair | O | 01 | |
| | | Overhaul | O | 01 | |
| 05 | Foot Control Valve | Inspect | O | 01 | |
| | | Service | O | 01 | |
| | | Test | O | 01 | |
| | | Repair | O | 01 | |
| | | Overhaul | O | 01 | |
| 06 | Signal Relay Valve | Inspect | O | 01 | |
| | | Service | O | 01 | |
| | | Test | O | 01 | |
| | | Repair | O | 01 | |
| | | Overhaul | O | 01 | |
| 07 | Chip Blower Valve | Inspect | O | 01 | |
| | | Service | O | 01 | |

(continued) Appendix C. Maintenance Allocation Chart

6520-01-272-4531
 Dental Operating Unit, ADEC Model 3406 Porta-Cart

[C – Operator or Crew; O – Unit Maintenance; DS – Direct Support Maintenance; GS – General Support Maintenance; D – Depot Maintenance]

| (1) GROUP NUMBER | (2) ASSEMBLY GROUP | (3) MAINTENANCE FUNCTION | (4) MAINTENANCE LEVEL | (5) TOOLS AND EQUIPMENT | (6) REMARKS |
|------------------------|--------------------------|--------------------------------|-----------------------------|----------------------------------|----------------|
| 08 | Three Way Toggle Valves | Test | O | 01 | |
| | | Repair | O | 01 | |
| | | Overhaul | O | 01 | |
| | | Inspect | O | 01 | |
| | | Service | O | 01 | |
| | | Test | O | 01 | |
| 09 | Needle Valves | Repair | O | 01 | |
| | | Overhaul | O | 01 | |
| | | Inspect | O | 01 | |
| | | Service | O | 01 | |
| | | Test | O | 01 | |
| 10 | Syringe | Repair | O | 01 | |
| | | Overhaul | O | 01 | |
| | | Inspect | O | 01 | |
| | | Service | O | 01 | |
| | | Test | O | 01 | |
| 11 | Air Vacuum System | Repair | O | 01 | |
| | | Overhaul | O | 01 | |
| | | Inspect | O | 01 | |
| | | Service | O | 01 | |
| | | Test | O | 01 | |
| 12 | Air Saliva Ejector | Repair | O | 01 | |
| | | Overhaul | O | 01 | |
| | | Inspect | O | 01 | |
| | | Service | O | 01 | |
| | | Test | O | 01 | |
| 13 | Storage Case | Repair | O | 01, 02 | |
| | | Overhaul | O | 01, 02 | |
| | | Inspect | O | 01 | |

(continued) Appendix C. Maintenance Allocation Chart

6520-01-398-4613
Compressor Dehydrator, Dental, Model PAC 6.7

[C – Operator or Crew; O – Unit Maintenance; DS – Direct Support Maintenance; GS – General Support Maintenance; D – Depot Maintenance]

| (1) GROUP NUMBER | (2) ASSEMBLY GROUP | (3) MAINTENANCE FUNCTION | (4) MAINTENANCE LEVEL | (5) TOOLS AND EQUIPMENT | (6) REMARKS |
|------------------------|------------------------------|--------------------------------|-----------------------------|----------------------------------|----------------|
| 00 | Compressor Dehydrator | Inspect | C | | |
| | | Test | C | | |
| | | Service | C | | |
| | | Adjust | O | 01, 02 | |
| | | Safety | O | 29 | |
| | | Repair | O | 01, 02 | |
| | | Overhaul | DS | 01, 02, 03 | |
| | | Rebuild | D | 01, 02, 03, 04 | |
| 01 | Compressor | Test | O | 01, 02 | |
| | | Service | O | 01, 02 | |
| | | Replace | O | 01, 02 | |
| | | Repair | O | 01, 02 | |
| | | Rebuild | D | 01, 02, 03, 04 | |
| 0101 | Intake Filter Element | Inspect | C | | |
| | | Replace | O | | |
| 0102 | Power Cord | Inspect | C | | |
| | | Replace | O | 01, 02 | |
| 0103 | Air Hose | Inspect | C | | |
| | | Replace | C | | |
| 02 | Fan | Test | O | 09 | |
| | | Replace | O | 01, 02 | |
| 03 | Cooling Coil | Inspect | O | 01, 02 | |
| | | Replace | O | 01, 02 | |
| 04 | Drying Chamber | Inspect | O | 01, 02 | |
| | | Replace | O | 01, 02 | |
| 05 | Storage Tank | Inspect | O | 01, 02 | |
| | | Replace | O | 01, 02 | |
| 0501 | Presssure Relief Drain Valve | Inspect | O | 01, 02 | |
| | | Replace | O | 01, 02 | |
| 06 | Dryness Indicator Disk | Inspect | C | | |
| | | Replace | O | 01, 02 | |
| 07 | Pressure Gauge | Test | O | 01, 02 | |
| | | Replace | O | 01, 02 | |
| | | Repair | O | 01, 02 | |

(continued) Appendix C. Maintenance Allocation Chart

6520-01-398-4613
Compressor Dehydrator, Dental, Model PAC 6.7

[C – Operator or Crew; O – Unit Maintenance; DS – Direct Support Maintenance; GS – General Support Maintenance; D – Depot Maintenance]

| (1) GROUP NUMBER | (2) ASSEMBLY GROUP | (3) MAINTENANCE FUNCTION | (4) MAINTENANCE LEVEL | (5) TOOLS AND EQUIPMENT | (6) REMARKS |
|------------------------|--------------------------|--------------------------------|-----------------------------|----------------------------------|----------------|
| 08 | Water Separator | Inspect Replace | O O | 01, 02 01, 02 | |
| 0801 | Muffler | Inspect Replace | C O | 01, 02 | |
| 09 | Case | Inspect Repair | C DS | 01, 02 01, 02, 03 | |
| 0901 | Latches | Inspect Replace | C O | 01, 02 | |

(continued) Appendix C. Maintenance Allocation Chart

6525-01-312-6411
X-Ray Apparatus, Radiographic/Fluoroscopic, Model CS-8952

[C – Operator or Crew; O – Unit Maintenance; DS – Direct Support Maintenance; GS – General Support Maintenance; D – Depot Maintenance]

| (1) GROUP NUMBER | (2) ASSEMBLY GROUP | (3) MAINTENANCE FUNCTION | (4) MAINTENANCE LEVEL | (5) TOOLS AND EQUIPMENT | (6) REMARKS |
|------------------------|--------------------------|--------------------------------|-----------------------------|----------------------------------|----------------|
| 00 | X-Ray Apparatus | Inspect | O | | |
| | | Test | O | 01, 02, 20, 21, 05, 14, 26 | |
| | | Service | O | 01, 02, 20, 21, 05, 14, 26 | |
| | | Repair | O | 01, 02, 20, 21 | |
| | | Replace | O | 01, 02 | |
| | | Overhaul | O | 01, 02 | |
| | | Rebuild | D | 01, 02, 04 | |
| 01 | Generator | Inspect | O | | |
| | | Test | O | 01, 02, 20, 21, 05, 14, 26 | |
| | | Service | O | 01, 02, 20, 21, 05, 14, 26 | |
| | | Repair | O | 01, 02, 20, 21 | |
| | | Replace | O | 01, 02, | |
| | | Overhaul | O | 01, 02 | |
| | | Rebuild | D | 01, 02, 04 | |
| 02 | Spot Film Device | Inspect | O | | |
| | | Test | O | 01, 02, 20, 21 | |
| | | Service | O | 01, 02 | |
| | | Repair | O | 01, 02, | |
| | | Replace | O | 01, 02 | |
| | | Overhaul | O | 01, 02, 20, 21 | |
| | | Rebuild | D | 01, 02, 04 | |
| 03 | Under-Table Collimator | Inspect | O | | |
| | | Test | O | 01, 02 | |
| | | Service | O | 01, 02 | |
| | | Repair | O | 01, 02 | |
| | | Replace | O | 01, 02 | |
| | | Overhaul | O | 01, 02, 04 | |
| | | Rebuild | D | 01, 02, 04 | |
| 04 | Over-Table Collimator | Inspect | O | | |
| | | Test | O | 01, 02, 20, 21 | |
| | | Service | O | 01, 02, 20, 21 | |
| | | Repair | O | 01, 02, 20 | |
| | | Replace | O | 01, 02 | |
| | | Overhaul | O | 01, 02, 20 | |

(continued) Appendix C. Maintenance Allocation Chart

6525-01-312-6411

X-Ray Apparatus, Radiographic/Fluoroscopic, Model CS-8952

[C – Operator or Crew; O – Unit Maintenance; DS – Direct Support Maintenance; GS – General Support Maintenance; D – Depot Maintenance]

| (1) GROUP NUMBER | (2) ASSEMBLY GROUP | (3) MAINTENANCE FUNCTION | (4) MAINTENANCE LEVEL | (5) TOOLS AND EQUIPMENT | (6) REMARKS |
|------------------------|----------------------------|--------------------------------|-----------------------------|----------------------------------|----------------|
| 05 | Automatic Exposure Control | Rebuild | D | 01, 02, 04 | |
| | | Inspect | O | | |
| | | Test | O | 01, 02, 20, 21 | |
| | | Service | O | 01, 02, 20, 21 | |
| | | Repair | O | 01, 02, 20 | |
| | | Replace | O | 01, 02 | |
| | | Overhaul | O | 01, 02 | |
| | | Rebuild | D | 01, 02, 04 | |
| 06 | Image Intensifier | Inspect | O | | |
| | | Test | O | 01, 02, 20, 21 | |
| | | Service | O | 01, 02, 20, 21 | |
| | | Repair | O | 01, 02, 20 | |
| | | Replace | O | 01, 02, 20 | |
| | | Overhaul | O | 01, 02, 20 | |
| | | Rebuild | D | 01, 02, 04, 20 | |
| | | 07 | Over-Table Tube | Inspect | O |
| Test | O | | | 01, 02, 20, 21 | |
| Service | O | | | 01, 02, 20, 21 | |
| Repair | D | | | 01, 02, 20 | |
| Replace | O | | | 01, 02, 20 | |
| Overhaul | D | | | 01, 02, 20 | |
| Rebuild | D | | | 01, 02, 04, 20 | |
| 08 | Under-Table Tube | | | Inspect | O |
| | | Test | O | 01, 02, 20, 21 | |
| | | Service | O | 01, 02, 20, 21 | |
| | | Repair | D | 01, 02, 20 | |
| | | Replace | O | 01, 02, 20 | |
| | | Overhaul | D | 01, 02, 20 | |
| | | Rebuild | D | 01, 02, 04, 20 | |
| | | 09 | Table | Inspect | O |
| Test | O | | | 01, 02, 20, 21 | |
| Service | O | | | 01, 02, 20, 21 | |
| Repair | O | | | 01, 02, 20 | |
| | | | | | |

(continued) Appendix C. Maintenance Allocation Chart

6525-01-312-6411
 X-Ray Apparatus, Radiographic/Fluoroscopic, Model CS-8952

[C – Operator or Crew; O – Unit Maintenance; DS – Direct Support Maintenance; GS – General Support Maintenance; D – Depot Maintenance]

| (1) GROUP NUMBER | (2) ASSEMBLY GROUP | (3) MAINTENANCE FUNCTION | (4) MAINTENANCE LEVEL | (5) TOOLS AND EQUIPMENT | (6) REMARKS |
|------------------------|--------------------------|--------------------------------|-----------------------------|----------------------------------|----------------|
| 010 | Tubestand | Replace | O | 01, 02, 20 | |
| | | Overhaul | O | 01, 02, 20 | |
| | | Rebuild | D | 01, 02, 04, 20 | |
| | | Inspect | O | | |
| | | Test | O | 01, 02, 20 | |
| | | Service | O | 01, 02, 20 | |
| | | Repair | O | 01, 02, 20 | |
| | | Replace | O | 01, 02, 20 | |
| | | Overhaul | O | 01, 02, 20 | |
| | | Rebuild | D | 01, 02, 04, 20 | |

(continued) Appendix C. Maintenance Allocation Chart

6525-01-325-3740
Portable X-Ray System, Model 1200

[C – Operator or Crew; O – Unit Maintenance; DS – Direct Support Maintenance; GS – General Support Maintenance; D – Depot Maintenance]

| (1) GROUP NUMBER | (2) ASSEMBLY GROUP | (3) MAINTENANCE FUNCTION | (4) MAINTENANCE LEVEL | (5) TOOLS AND EQUIPMENT | (6) REMARKS |
|------------------------|----------------------------|--------------------------------|-----------------------------|-------------------------------|----------------|
| 00 | X-Ray System | Inspect | O | | |
| | | Test | O | 01, 02, 20, 21, 05, 14, 26 | |
| | | Service | O | 01, 02, 20, 21, 05, 14, 26 | |
| | | Repair | O | 01, 02, 20, 21 | |
| | | Replace | O | 01, 02, 20 | |
| | | Overhaul | O | 01, 02, 20 | |
| | | Rebuild | D | 01, 02, 04, 20 | |
| 01 | Control Assembly | Inspect | O | | |
| | | Test | O | 01, 02, 20, 21 | |
| | | Service | O | 01, 02, 20, 21 | |
| | | Repair | O | 01, 02, 20 | |
| | | Replace | O | 01, 02, 20 | |
| | | Overhaul | O | 01, 02, 20, 21 | |
| | | Rebuild | D | 01, 02, 04, 20, 21 | |
| 02 | Panel Assembly, Control | Inspect | O | | |
| | | Test | O | 01, 02, 20 | |
| | | Service | O | 01, 02, 20 | |
| | | Repair | O | 01, 02, 20 | |
| | | Replace | O | 01, 02, 20 | |
| | | Overhaul | O | 01, 02, 20 | |
| | | Rebuild | D | 01, 02, 04 | |
| 03 | PCB Assembly, Mother Board | Inspect | O | | |
| | | Test | O | 01, 02, 20, 21 | |
| | | Service | O | 01, 02, 20, 21 | |
| | | Repair | O | 01, 02, 20, 21 | |
| | | Replace | O | 01, 02, | |
| | | Overhaul | O | 01, 02, 20, 21 | |
| | | Rebuild | D | 01, 02, 04, 20, 21 | |
| 04 | PCB Assembly, Timer Board | Inspect | O | | |
| | | Test | O | 01, 02, 20, 21 | |
| | | Service | O | 01, 02, 20, 21 | |
| | | Repair | O | 01, 02, 20, 21 | |
| | | Replace | O | 01, 02, 20, 21 | |
| | | Overhaul | O | 01, 02, 20, 21 | |
| | | Rebuild | D | 01, 02, 04, 20, 21 | |

(continued) Appendix C. Maintenance Allocation Chart

6525-01-325-3740
 Portable X-Ray System, Model 1200

[C – Operator or Crew; O – Unit Maintenance; DS – Direct Support Maintenance; GS – General Support Maintenance; D – Depot Maintenance]

| (1) GROUP NUMBER | (2) ASSEMBLY GROUP | (3) MAINTENANCE FUNCTION | (4) MAINTENANCE LEVEL | (5) TOOLS AND EQUIPMENT | (6) REMARKS |
|------------------------|-----------------------------------|--------------------------------|-----------------------------|-------------------------------|----------------|
| 05 | PCB Assembly, Line Set | Inspect | O | | |
| | | Test | O | 01, 02, 20, 21 | |
| | | Service | O | 01, 02, 20, 21 | |
| | | Repair | O | 01, 02, 20, 21 | |
| | | Replace | O | 01, 02, 20, 21 | |
| | | Overhaul | O | 01, 02, 20, 21 | |
| | | Rebuild | D | 01, 02, 04, 20, 21 | |
| 06 | PCB Assembly, MAS Interface Board | Inspect | O | | |
| | | Test | O | 01, 02, 20, 21 | |
| | | Service | O | 01, 02, 20, 21 | |
| | | Repair | O | 01, 02, 20, 21 | |
| | | Replace | O | 01, 02, 20, 21 | |
| | | Overhaul | O | 01, 02, 20, 21 | |
| | | Rebuild | D | 01, 02, 04, 20, 21 | |
| 07 | Plate Assembly, Base | Inspect | O | | |
| | | Test | O | 01, 02 | |
| | | Service | O | 01, 02 | |
| | | Repair | O | 01, 02 | |
| | | Replace | O | 01, 02 | |
| | | Overhaul | O | 01, 02 | |
| | | Rebuild | D | 01, 02, 04 | |
| 08 | Switch Assembly, Exposure | Inspect | O | | |
| | | Test | O | 01, 02, 20 | |
| | | Service | O | 01, 02, 20 | |
| | | Repair | O | 01, 02, 20 | |
| | | Replace | O | 01, 02, 20 | |
| | | Overhaul | O | 01, 02, 20 | |
| | | Rebuild | D | 01, 02, 04, 20 | |
| 09 | Case, Control Assembly | Inspect | O | | |
| | | Test | O | 01, 02 | |
| | | Service | O | 01, 02 | |
| | | Repair | O | 01, 02 | |
| | | Replace | O | 01, 02 | |
| | | Overhaul | O | 01, 02 | |

(continued) Appendix C. Maintenance Allocation Chart

6525-01-325-3740
Portable X-Ray System, Model 1200

[C – Operator or Crew; O – Unit Maintenance; DS – Direct Support Maintenance; GS – General Support Maintenance; D – Depot Maintenance]

| (1) GROUP NUMBER | (2) ASSEMBLY GROUP | (3) MAINTENANCE FUNCTION | (4) MAINTENANCE LEVEL | (5) TOOLS AND EQUIPMENT | (6) REMARKS |
|------------------------|--------------------------|--------------------------------|-----------------------------|-------------------------------|----------------|
| 10 | Cord, Line, Assembly | Rebuild | D | 01, 02, 04 | |
| | | Inspect | O | | |
| | | Test | O | 01, 02, 20 | |
| | | Service | O | 01, 02, 20 | |
| | | Repair | O | 01, 02, 20 | |
| | | Replace | O | 01, 02, 20 | |
| | | Overhaul | O | 01, 02, 20 | |
| | | Rebuild | D | 01, 02, 04, 20 | |
| 11 | Cable Assembly | Inspect | O | | |
| | | Test | O | 01, 02, 20 | |
| | | Service | O | 01, 02, 20 | |
| | | Repair | O | 01, 02, 20 | |
| | | Replace | O | 01, 02, 20 | |
| | | Overhaul | O | 01, 02, 20 | |
| | | Rebuild | D | 01, 02, 04, 20 | |
| 12 | Harness Assembly | Inspect | O | | |
| | | Test | O | 01, 02, 20 | |
| | | Service | O | 01, 02, 20 | |
| | | Repair | O | 01, 02, 20 | |
| | | Replace | O | 01, 02, 20 | |
| | | Overhaul | O | 01, 02, 20 | |
| | | Rebuild | D | 01, 02, 04, 20 | |
| 13 | Generator Assembly | Inspect | O | | |
| | | Test | O | 01, 02, 20, 21 | |
| | | Service | O | 01, 02, 20, 21 | |
| | | Repair | O | 01, 02, 20 | |
| | | Replace | O | 01, 02, 20 | |
| | | Overhaul | O | 01, 02, 20, 21 | |
| | | Rebuild | D | 01, 02, 04, 20, 21 | |
| 14 | Collimator Assembly | Inspect | O | | |
| | | Test | O | 01, 02, 20 | |
| | | Service | O | 01, 02, 20 | |
| | | Repair | O | 01, 02, 20 | |
| | | Replace | O | 01, 02, 20 | |
| | | Overhaul | O | 01, 02, 20 | |

(continued) Appendix C. Maintenance Allocation Chart

6525-01-325-3740
Portable X-Ray System, Model 1200

[C – Operator or Crew; O – Unit Maintenance; DS – Direct Support Maintenance; GS – General Support Maintenance; D – Depot Maintenance]

| (1) GROUP NUMBER | (2) ASSEMBLY GROUP | (3) MAINTENANCE FUNCTION | (4) MAINTENANCE LEVEL | (5) TOOLS AND EQUIPMENT | (6) REMARKS |
|------------------------|------------------------------|--------------------------------|-----------------------------|-------------------------------|----------------|
| 15 | Tube Head Assembly | Rebuild | D | 01, 02, 04, 20 | |
| | | Inspect | O | | |
| | | Test | O | 01, 02, 20 | |
| | | Service | D | 01, 02, 04, 20 | |
| | | Repair | D | 01, 02, 04, 20 | |
| | | Replace | O | 01, 02, 20 | |
| | | Overhaul | D | 01, 02, 04, 20 | |
| | | Rebuild | D | 01, 02, 04, 20 | |
| 16 | Yoke Assembly | Inspect | O | | |
| | | Test | O | 01, 02 | |
| | | Service | O | 01, 02 | |
| | | Repair | O | 01, 02 | |
| | | Replace | O | 01, 02 | |
| | | Overhaul | O | 01, 02 | |
| | | Rebuild | D | 01, 02, 04 | |
| 17 | Chassis, Generator, Assembly | Inspect | O | | |
| | | Test | O | 01, 02 | |
| | | Service | O | 01, 02 | |
| | | Repair | O | 01, 02 | |
| | | Replace | O | 01, 02 | |
| | | Overhaul | O | 01, 02 | |
| | | Rebuild | D | 01, 02, 04 | |
| 18 | Stand Assembly | Inspect | O | | |
| | | Test | O | 01, 02 | |
| | | Service | O | 01, 02 | |
| | | Repair | O | 01, 02 | |
| | | Replace | O | 01, 02 | |
| | | Overhaul | O | 01, 02 | |
| | | Rebuild | D | 01, 02, 04 | |
| 19 | Frame Assembly, Stand | Inspect | O | | |
| | | Test | O | 01, 02 | |
| | | Service | O | 01, 02 | |
| | | Repair | O | 01, 02 | |
| | | Replace | O | 01, 02 | |

(continued) Appendix C. Maintenance Allocation Chart

6525-01-325-3740
 Portable X-Ray System, Model 1200

[C – Operator or Crew; O – Unit Maintenance; DS – Direct Support Maintenance; GS – General Support Maintenance; D – Depot Maintenance]

| (1) GROUP NUMBER | (2) ASSEMBLY GROUP | (3) MAINTENANCE FUNCTION | (4) MAINTENANCE LEVEL | (5) TOOLS AND EQUIPMENT | (6) REMARKS |
|------------------------|--------------------------|--------------------------------|-----------------------------|-------------------------------|----------------|
| 20 | Cross Arm Assembly | Overhaul | O | 01, 02 | |
| | | Rebuild | D | 01, 02, 04 | |
| | | Inspect | O | | |
| | | Test | O | 01, 02 | |
| | | Service | O | 01, 02 | |
| | | Repair | O | 01, 02 | |
| | | Replace | O | 01, 02 | |
| | | Overhaul | O | 01, 02 | |
| | | Rebuild | D | 01, 02, 04 | |
| 21 | Gear Box Assembly | Inspect | O | | |
| | | Test | O | 01, 02 | |
| | | Service | O | 01, 02 | |
| | | Repair | O | 01, 02 | |
| | | Replace | O | 01, 02 | |
| | | Overhaul | O | 01, 02 | |
| | | Rebuild | D | 01, 02, 04 | |
| 22 | Pipe Assembly | Inspect | O | | |
| | | Test | O | | |
| | | Service | O | 01, 02 | |
| | | Repair | O | 01, 02 | |
| | | Replace | O | 01, 02 | |
| | | Overhaul | O | 01, 02 | |
| | | Rebuild | D | 01, 02, 04 | |
| 23 | Container, Reuseable | Inspect | O | | |
| | | Test | O | 01, 02 | |
| | | Service | O | 01, 02 | |
| | | Repair | O | 01, 02 | |
| | | Replace | O | 01, 02 | |
| | | Overhaul | O | 01, 02 | |
| | | Rebuild | D | 01, 02, 04 | |

(continued) Appendix C. Maintenance Allocation Chart
 25-01-370-7552
 Portable Dental X-Ray System, Model ALPHA MPDX

[C – Operator or Crew; O – Unit Maintenance; DS – Direct Support Maintenance; GS – General Support Maintenance; D – Depot Maintenance]

| (1) GROUP NUMBER | (2) ASSEMBLY GROUP | (3) MAINTENANCE FUNCTION | (4) MAINTENANCE LEVEL | (5) TOOLS AND EQUIPMENT | (6) REMARKS |
|------------------------|----------------------------|--------------------------------|-----------------------------|-------------------------------|----------------|
| 00 | X-Ray System | Inspect | C | | |
| | | Test | O | 01, 02, 20, 21, 05, 14, 26 | |
| | | Service | O | 01, 02, 20, 21, 05, 14, 26 | |
| | | Repair | O | 01, 02 | |
| | | Replace | O | 01, 02 | |
| | | Overhaul | D | 01, 02 | |
| | | Rebuild | D | 01, 02, 04 | |
| 01 | X-Ray Control Assembly | Inspect | C | | |
| | | Test | O | 01, 02, 20, 21, 05, 14, 26 | |
| | | Service | O | 01, 02, 20, 21, 05, 14, 26 | |
| | | Repair | O | 01, 02 | |
| | | Replace | O | 01, 02 | |
| | | Overhaul | O | 01, 02 | |
| | | Rebuild | D | 01, 02, 04 | |
| 02 | X-Ray Source Assembly | Inspect | O | | |
| | | Test | O | 01, 02, 20, 21, 05, 14, 26 | |
| | | Service | O | 01, 02, 20, 21, 05, 14, 26 | |
| | | Repair | O | 01, 02 | |
| | | Replace | O | 01, 02 | |
| | | Overhaul | O | 01, 02 | |
| | | Rebuild | D | 01, 02, 04 | |
| 0201 | X-Ray Tubehead Subassembly | Inspect | O | | |
| | | Test | O | 01, 02 | |
| | | Service | O | 01, 02 | |
| | | Repair | O | 01, 02 | |
| | | Replace | O | 01, 02 | |
| | | Overhaul | D | 01, 02, 04 | |
| | | Rebuild | D | 01, 02, 04 | |
| 0202 | Dental Cone Subassembly | Inspect | C | | |
| | | Test | O | 01, 02 | |
| | | Service | O | 01, 02 | |
| | | Repair | O | 01, 02 | |
| | | Replace | O | 01, 02 | |
| | | Overhaul | O | 01, 02 | |
| | | Rebuild | D | 01, 02, 04 | |
| 03 | Scissor Arm Assembly | Inspect | C | | |

(continued) Appendix C. Maintenance Allocation Chart

6525-01-370-7552
 Portable Dental X-Ray System, Model ALPHA MPDX

[C – Operator or Crew; O – Unit Maintenance; DS – Direct Support Maintenance; GS – General Support Maintenance; D – Depot Maintenance]

| (1) GROUP NUMBER | (2) ASSEMBLY GROUP | (3) MAINTENANCE FUNCTION | (4) MAINTENANCE LEVEL | (5) TOOLS AND EQUIPMENT | (6) REMARKS |
|------------------------|--------------------------|--------------------------------|-----------------------------|-------------------------------|----------------|
| 04 | Chair Unit | Test | C | | |
| | | Service | O | 01, 02 | |
| | | Repair | O | 01, 02 | |
| | | Replace | O | 01, 02 | |
| | | Overhaul | O | 01, 02 | |
| | | Rebuild | D | 01, 02, 04 | |
| | | Inspect | C | | |
| | | Test | C | | |
| | | Service | O | 01, 02 | |
| | | Repair | O | 01, 02 | |
| 0401 | Headrest Assembly | Replace | O | 01, 02 | |
| | | Overhaul | O | 01, 02 | |
| | | Rebuild | D | 01, 02, 04 | |
| | | Inspect | C | | |
| | | Test | C | | |
| | | Service | O | 01, 02 | |
| | | Repair | O | 01, 02 | |
| | | Replace | O | 01, 02 | |
| | | Overhaul | O | 01, 02 | |
| | | Rebuild | D | 01, 02, 04 | |
| 05 | Carrying Case | Inspect | C | | |
| | | Test | C | | |
| | | Service | O | 01, 02 | |
| | | Repair | O | 01, 02 | |
| | | Replace | O | 01, 02 | |
| | | Overhaul | O | 01, 02 | |
| | | Rebuild | D | 01, 02, 04 | |

(continued) Appendix C. Maintenance Allocation Chart

6525-01-384-9296
X-Ray Apparatus, Model LCROKS

[C – Operator or Crew; O – Unit Maintenance; DS – Direct Support Maintenance; GS – General Support Maintenance; D – Depot Maintenance]

| (1) GROUP NUMBER | (2) ASSEMBLY GROUP | (3) MAINTENANCE FUNCTION | (4) MAINTENANCE LEVEL | (5) TOOLS AND EQUIPMENT | (6) REMARKS |
|------------------------|--------------------------|--------------------------------|-----------------------------|-------------------------------|----------------|
| 00 | X-Ray Apparatus | Inspect | C, O | | |
| | | Test | O | 01, 02, 20, 21, 05, 14, 26 | |
| | | Service | O | 01, 02, 20, 21, 05, 14, 26 | |
| | | Calibrate | O | 01, 02, 20, 21, 05, 14, 26 | |
| | | Repair | O | 01, 02 | |
| | | Replace | O | 01, 02 | |
| | | Overhaul | O | 01, 02 | |
| | | Rebuild | D | 01, 02, 04 | |
| 01 | Operator Console | Inspect | C, O | | |
| | | Test | O | 01, 02, 20, 21 | |
| | | Service | O | 01, 02, 20 | |
| | | Calibrate | O | 01, 02, 20, 21 | |
| | | Repair | O | 01, 02 | |
| | | Replace | O | 01, 02 | |
| | | Overhaul | O | 01, 02 | |
| | | Rebuild | D | 01, 02, 04 | |
| 02 | Electronics Cabinet | Inspect | O | | |
| | | Test | O | 01, 02, 20, 21 | |
| | | Service | O | 01, 02, 20, 21 | |
| | | Calibrate | O | 01, 02, 20, 21 | |
| | | Repair | O | 01, 02 | |
| | | Replace | O | 01, 02 | |
| | | Overhaul | O | 01, 02 | |
| | | Rebuild | D | 01, 02, 04 | |
| 03 | Electronics Chassis | Inspect | O | | |
| | | Test | O | 01, 02, 20, 21 | |
| | | Service | O | 01, 02, 20 | |
| | | Calibrate | O | 01, 02, 20 | |
| | | Repair | O | 01, 02 | |
| | | Replace | O | 01, 02 | |
| | | Overhaul | O | 01, 02 | |
| | | Rebuild | D | 01, 02, 04 | |
| 04 | Inverter Chassis | Inspect | O | | |
| | | Test | O | 01, 02, 20, 21 | |
| | | Service | O | 01, 02, 20, 21 | |
| | | Calibrate | O | 01, 02, 20, 21 | |

(continued) Appendix C. Maintenance Allocation Chart

6525-01-384-9296
X-Ray Apparatus, Model LCROKS

[C – Operator or Crew; O – Unit Maintenance; DS – Direct Support Maintenance; GS – General Support Maintenance; D – Depot Maintenance]

| (1) GROUP NUMBER | (2) ASSEMBLY GROUP | (3) MAINTENANCE FUNCTION | (4) MAINTENANCE LEVEL | (5) TOOLS AND EQUIPMENT | (6) REMARKS |
|------------------------|--------------------------|--------------------------------|------------------------------|-------------------------------|----------------|
| 05 | High-Tension Transformer | Replace | O | 01, 02 | |
| | | Overhaul | O | 01, 02 | |
| | | Rebuild | D | 01, 02, 04 | |
| | | Inspect | O | | |
| | | Test | O | 01, 02, 20, 21 | |
| | | Service | O | 01, 02, 20, 21 | |
| | | Calibrate | O | 01, 02, 20, 21 | |
| | | Repair | O | 01, 02 | |
| | | Replace | O | 01, 02 | |
| | | Overhaul | O | 01, 02 | |
| | | Rebuild | D | 01, 02, 04 | |
| | | 06 0601 | Cables Power Cable | Inspect | C, O |
| Test | O | | | | |
| Service | O | | | | |
| Repair | O | | | 01, 02 | |
| Replace | O | | | 01, 02 | |
| Overhaul | O | | | 01, 02 | |
| Rebuild | D | | | 01, 02 | |
| 0602 | Interconnecting Cables | | | Inspect | O |
| | | Test | O | 01, 02, 20 | |
| | | Service | O | 01, 02, 20 | |
| | | Repair | O | 01, 02 | |
| | | Replace | O | 01, 02 | |
| | | Overhaul | O | 01, 02 | |
| | | Rebuild | D | 01, 02, 04 | |
| | | 0603 | High-Tension Generator Cable | Inspect | O |
| Test | O | | | 01, 02, 20 | |
| Service | O | | | 01, 02, 20 | |
| Repair | O | | | 01, 02 | |
| Replace | O | | | 01, 02 | |
| Overhaul | O | | | 01, 02 | |
| Rebuild | D | | | 01, 02, 04 | |
| 0604 | Rotor Drive Cable | | | Inspect | O |
| | | Test | O | 01, 02 | |
| | | Service | O | 01, 02 | |
| | | Repair | O | 01, 02 | |
| | | Replace | O | 01, 02 | |

(continued) Appendix C. Maintenance Allocation Chart

6525-01-384-9296
 X-Ray Apparatus, Model LCROKS

[C – Operator or Crew; O – Unit Maintenance; DS – Direct Support Maintenance; GS – General Support Maintenance; D – Depot Maintenance]

| (1) GROUP NUMBER | (2) ASSEMBLY GROUP | (3) MAINTENANCE FUNCTION | (4) MAINTENANCE LEVEL | (5) TOOLS AND EQUIPMENT | (6) REMARKS |
|------------------------|--------------------------|--------------------------------|-----------------------------|-------------------------------|----------------|
| | | Overhaul Rebuild | O | 01, 02 01, 02, 04 | |

(continued) Appendix C. Maintenance Allocation Chart

6530-00-926-2151

Sterilizer, Surgical Dressing 16x36 in., Model M-138

[C – Operator or Crew; O – Unit Maintenance; DS – Direct Support Maintenance; GS – General Support Maintenance; D – Depot Maintenance]

| (1) GROUP NUMBER | (2) ASSEMBLY GROUP | (3) MAINTENANCE FUNCTION | (4) MAINTENANCE LEVEL | (5) TOOLS AND EQUIPMENT | (6) REMARKS |
|------------------------|----------------------------------|-------------------------------------------|-----------------------------|----------------------------------|----------------|
| 00 | Sterilizer | Inspect Test Electrical Safety Test | C O O | 01, 02 01, 02, 20, 29 | |
| 01 | Heater Assembly | Test Replace | O O | 01, 02, 29 01, 02 | |
| 02 | Control Box Assembly | Inspect Repair | O O | 01, 02 01, 02, 29 | |
| 0201 | Relay, Armature | Test Replace | O O | 01, 02, 29 01, 02 | |
| 0202 | Pressure Control | Test Replace | O O | 01, 02, 29 01, 02 | |
| 0203 | Pilot Light | Test Replace | O O | 01, 02, 29 01, 02 | |
| 0204 | Lamp, Neon | Test Replace | O O | 01, 02, 29 01, 02 | |
| 0205 | Switch, Toggle | Test Replace | O O | 01, 02, 29 01, 02 | |
| 0206 | Switch, Low Water Cut-off | Test Replace | O O | 01, 02, 29 01, 02 | |
| 0207 | Block, Terminal | Replace | O | 01, 02 | |
| 03 | Operating Valve Assembly | Test Repair Replace | O O O | 01, 02 01, 02 01, 02 | |
| 04 | Door Assembly | Test Service Repair | O O O | 01, 02 01, 02 01, 02 | |
| 0401 | Packing, Preformed (Door Gasket) | Inspect Replace | C O | 01, 02 | |

(continued) Appendix C. Maintenance Allocation Chart

6530-00-926-2151
 Sterilizer, Surgical Dressing 16x36 in., Model M-138

[C – Operator or Crew; O – Unit Maintenance; DS – Direct Support Maintenance; GS – General Support Maintenance; D – Depot Maintenance]

| (1) GROUP NUMBER | (2) ASSEMBLY GROUP | (3) MAINTENANCE FUNCTION | (4) MAINTENANCE LEVEL | (5) TOOLS AND EQUIPMENT | (6) REMARKS |
|------------------------|--------------------------|--------------------------------|-----------------------------|----------------------------------|----------------|
| 05 | Vacuum Dryer Assembly | Test Replace | O O | 01, 02 01, 02 | |
| 06 | Gauges | Inspect Test Replace | C O O | 01, 02 01, 02 | |
| 07 | Timer | Test Replace | O O | 01, 02 01, 02 | |
| 08 | Case, Transport | Inspect Repair Overhaul | C O D | 01, 02 04 | |
| 09 | Shelves | Inspect Replace | C C | | |
| 10 | Chamber | Inspect Replace Overhaul | C O D | 01, 02 01, 02, 04 | |

(continued) Appendix D. Equipment Parts and Accessories List

6520-01-272-4531
Dental Operating Unit, ADEC Model 3406 Porta-Cart

| ITEM NO | Part or Accessory Description | Part or Accessory Number | Basic Issue | Operator | Repairer |
|---------|--------------------------------------------------------|--------------------------|-------------|----------|----------|
| 1 | Operator/Service Manual | 65-634 | 2 | 1 | 1 |
| 2 | Century II Automatic Control System for Two Handpieces | 38-0191-00 | 1 | 1 | 1 |
| 3 | Air Coolant Flow Control | 13-0361-00 | 1 | 1 | 1 |
| 4 | Water Coolant On-Off Toggle | 33-0048-00 | 1 | 1 | 1 |
| 5 | Water Coolant Flow Control for each Handpiece | 13-0361-00 | 1 | 1 | 1 |
| 6 | Drive Air Pressure Gauge | 026-009-00 | 1 | 1 | 1 |
| 7 | Disc Type Foot Control with Chip Blower Button | 38-0251-00 | 1 | 1 | 1 |
| 8 | Self Contained Two-Quart Water Tank | 36-0023-00 | 1 | 1 | 1 |
| 9 | Water Pressure On-Off Toggle | 33-0048-00 | 1 | 1 | 1 |
| 10 | Air Vacuum System | 10-0729-00 | 1 | 1 | 1 |
| 11 | AVS Handpiece Assembly | 10-0716-00 | 1 | 1 | 1 |
| 12 | Oral Evacuator Flow Control | 13-0361-00 | 1 | 1 | 1 |
| 13 | Air Saliva Ejector with Solids Separator | 12-0070-00 | 1 | 1 | 1 |
| 14 | Polypropylene Waste Bottle | 17-0270-00 | 1 | 1 | 1 |
| 15 | Soft Touch Button Syringe with Coiled Tubing | 23-0088-00 | 1 | 1 | 1 |
| 16 | Quick-Disconnect Water Outlet | 026-065-00 | 1 | 1 | 1 |
| 17 | Adjustable Height Frame | 36-0015-00 | 1 | 1 | 1 |
| 18 | Fiberglass Carrying Case | 36-0070-00 | 1 | 1 | 1 |
| 19 | Stainless Steel Tray, 15 1/8" x 10 5/8" | 043-003-00 | 1 | 1 | 1 |
| 20 | Two Coiled Handpiece Tubing with Midwest-Connectors | 98-0448-00 | 1 | 1 | 1 |
| 21 | 10 Foot Air Supply Tubing with Quick Disconnects | 45-0182-00 | 1 | 1 | 1 |
| 22 | Water Tank Filler Funnel | 009-003-00 | 1 | 1 | 1 |
| 23 | Stainless Steel Dry Oral Cup | 11-0450-00 | 1 | 1 | 1 |
| 24 | Stainless Steel Oral Evacuator Tips | 10-0010-00 | 1 | 1 | 1 |
| 25 | Porta-Cart Accessory Kit | 36-0089-00 | 1 | 1 | 1 |
| | a. Case | 36-0084-00 | 1 | 1 | 1 |
| | b. Caster | 16-0080-00 | 4 | 4 | 4 |
| | c. Service Kit, Air-Filter/Regulator | 90-0030-00 | 1 | 1 | 1 |
| | d. Service Kit, Century II | 90-0308-00 | 1 | 1 | 1 |
| | e. Stainless Steel Tip | 10-0010-00 | 2 | 2 | 2 |
| | f. Male 1/4" Quick-Disconnect | 026-035-00 | 1 | 1 | 1 |
| | g. Water Cup Filler Tube | 17-0240-00 | 1 | 1 | 1 |
| | h. 3/4" Open End Wrench | 009-004-00 | 1 | 1 | 1 |

(continued) Appendix D. Equipment Parts and Accessories List

6520-01-272-4531
 Dental Operating Unit, ADEC Model 3406 Porta-Cart

| ITEM NO | Part or Accessory Description | Part or Accessory Number | Basic Issue | Operator | Repairer |
|---------|---------------------------------------|--------------------------|-------------|----------|----------|
| | i. ½" Diameter Brush | 049-001-00 | 1 | 1 | 1 |
| | j. Offset Screwdriver | 009-001-00 | 1 | 1 | 1 |
| | k. Plastic Sleeve Tool | 98-0072-00 | 1 | 1 | 1 |
| | l. Syringe Tips | 23-0872-00 | 2 | 2 | 2 |
| | m. Syringe Service Kit | 90-0310-00 | 1 | 1 | 1 |
| | n. AVS Locking Button Kit | 10-0600-00 | 1 | 1 | 1 |
| 26 | Star Futura F303 High Speed Handpiece | 53874 | 1 | 1 | 1 |
| 27 | Star Titan 2 TA202M | 3055 | 1 | 1 | 1 |

(continued) Appendix D. Equipment Parts and Accessories List

6520-01-398-4613
Compressor Dehydrator, Dental, Model PAC 6.7

| ITEM NO | Part or Accessory Description | Part or Accessory Number | Basic Issue | Operator | Repairer |
|---------|--------------------------------------------------------------------------------------------------------------------------------------|--------------------------|-------------|----------|----------|
| 1 | Interconnecting Air Hoses; 10-foot section with appropriate connectors (connects Compressor to Dental Operating and Treatment Unit). | PAC6.7-035 | 2 | 1 | 1 |
| 2 | Technical Manual; Complete Operating and Maintenance Instructions. | PAC6.7M | 2 | 1 | 1 |

(continued) Appendix D. Equipment Parts and Accessories List

6525-01-312-6411

X-Ray Apparatus, Radiographic/Fluoroscopic, Model CS-8952

| ITEM NO | Part or Accessory Description | Part or Accessory Number | Basic Issue | Operator | Repairer |
|---------|----------------------------------------------|---------------------------|-------------|----------|----------|
| 1 | Tools | | | | |
| 2 | Manuals, Service (2 Volumes), and Operator | 9023.400 | 1 set | 1 | 1 |
| 3 | Tubestand | CTC | 1 | 1 | 1 |
| 4 | X-Ray Table | 5736.062.02 | 1 | 1 | 1 |
| 5 | Radiation Shield | | 1 | 1 | 1 |
| 6 | X-Ray Generator control | MXR-350 | 1 | 1 | 1 |
| 7 | Mobile Cassette Stand, with Shipping Crate | | 1 | 1 | 1 |
| 8 | Spot Film Device, with Shipping Crate | EXT-950 (DPSC) | 1 | 1 | 1 |
| 9 | Transformer, High Voltage | | 1 | 1 | 1 |
| 10 | Auxiliary Cabinet | | 1 | 1 | 1 |
| 11 | X-Ray Tube, Over-Table, with Shipping Crate | RAD-13 | 1 | 1 | 1 |
| 12 | X-Ray tube, Under-Table, with Shipping Crate | RAD-14 | 1 | 1 | 1 |
| 13 | Collimator, Over-Table, with Shipping Crate | 70-08040 LINEAR II (DPSC) | 1 | 1 | 1 |
| 14 | Collimator, Under-Table, with Shipping Crate | LINEAR FSE | 1 | 1 | 1 |
| 15 | Image Intensifier, with Shipping Crate | | 1 | 1 | 1 |
| 16 | Interconnecting Cables | | | | |
| 17 | Patient Handgrips | 5236.500.02 | 1 set | 1 set | 1 set |
| 18 | Urological Knee Crutches | 5436.504.02 | 1 set | 1 set | 1 set |
| 19 | Compression Band Device | 5536.500.01 | 1 | 1 | 1 |
| 20 | Shoulder Rest | 5536.504.01 | 1 set | 1 set | 1 set |
| 21 | Foot Rest | 5536.504.02 | 1 set | 1 set | 1 set |
| 22 | Head Clamp | 9491.201 | 1 | 1 | 1 |
| 23 | Over-Table HV Cables | | 1 set | 1 set | 1 set |
| 24 | Under-Table HV Cables | | 1 set | 1 set | 1 set |

(continued) Appendix D. Equipment Parts and Accessories List

6525-01-325-3740
 Portable X-Ray System, Model 1200

| ITEM NO | Part or Accessory Description | Part or Accessory Number | Basic Issue | Operator | Repairer |
|---------|--------------------------------|---------------------------|-------------|----------|----------|
| 1 | Operator/Service Manual | Model 1200 Service Manual | 2 | 1 | 1 |
| 2 | Control Assembly | 500507 | 1 | 1 | 1 |
| 3 | Exposure Switch Assembly | 500512 | 1 | 1 | 1 |
| 4 | X-Ray Generator Assembly | 500542 | 1 | 1 | 1 |
| 5 | Stand Assembly | 500592 | 1 | 1 | 1 |
| 6 | Reusable Storage Container | 201047 | 1 | 1 | 1 |
| 7 | Line Cord | 500081 | 1 | 1 | 1 |
| 8 | Interconnecting Cable Assembly | 500595 | 1 | 1 | 1 |

(continued) Appendix D. Equipment Parts and Accessories List

6525-01-370-7552
Portable Dental X-Ray System, Model ALPHA MPDX

| ITEM NO | Part or Accessory Description | Part or Accessory Number | Basic Issue | Operator | Repairer |
|---------|-------------------------------------|--------------------------|-------------|----------|----------|
| 1 | Carrying Case Unit | 500883 | 1 | 1 | 1 |
| | Parts Shipped in Cover | | | | |
| 2 | a. Chair Unit | 500892 | 1 | 1 | 1 |
| 3 | 1. Seat | 218110 | 1 | 1 | 1 |
| 4 | 2. Backrest | 218111 | 1 | 1 | 1 |
| 5 | 3. Headrest | 218108 | 1 | 1 | 1 |
| 6 | 4. Supporting assemblies | 800327 | 1 | 1 | 1 |
| 7 | b. Truss Arm | 500914 | 1 | 1 | 1 |
| 8 | c. Telescopic Leg | 218117 | 1 | 1 | 1 |
| | Parts Shipped in Case Lower Section | | | | |
| 9 | a. X-Ray Unit | 500876 | 1 | 1 | 1 |
| 10 | 1. X-Ray Control Assembly | 500877 | 1 | 1 | 1 |
| 11 | 2. X-Ray Source Assembly | 500891 | 1 | 1 | 1 |
| 12 | 3. Scissor Arm Assembly | 500882 | 1 | 1 | 1 |
| 13 | b. Steel Support Braces, Identical | 102860 | 2 | 2 | 2 |
| 14 | c. Extension Tube | 500885 | 1 | 1 | 1 |
| 15 | d. Leveling Brackets/Leveling Pads | 500884 | 4 | 4 | 4 |
| 16 | e. Dental Cone | 500897 | 1 | 1 | 1 |
| 17 | f. Line Cord | 500139 | 1 | 1 | 1 |
| 18 | g. Exposure Switch | 500902 | 1 | 1 | 1 |
| 19 | h. Operation Manual | 500893/OM | 1 | 1 | 1 |
| 20 | i. Maintenance Manual | 500893/MM | 1 | 1 | 1 |

(continued) Appendix D. Equipment Parts and Accessories List

6525-01-384-9296

X-Ray Apparatus, Model LCROKS

| ITEM NO | Part or Accessory Description | Part or Accessory Number | Basic Issue | Operator | Repairer |
|---------|-------------------------------------|--------------------------|-------------|----------|----------|
| 1 | Operators Literature | | 2 | 1 | 1 |
| 2 | Service Literature | | 2 | 1 | 1 |
| 3 | Generix- SynerGen- Console | 1173 | 1 | 1 | 1 |
| 4 | Control Unit CLINIX VP4 | 1184 | 1 | 1 | 1 |
| 5 | HT- Generator | 1243 | 1 | 1 | 1 |
| 6 | CLINIX VP4 | 1342 | 1 | 1 | 1 |
| 7 | Column Assembly | 1738 | 1 | 1 | 1 |
| 8 | Control Handle | 15675 | 1 | 1 | 1 |
| 9 | Collimator | 20 3072 | 1 | 1 | 1 |
| 10 | Cable Set | 1738-4-6 | 1 | 1 | 1 |
| 11 | Addition Kit for EP-Bucky Assembly | 1342-5 | 1 | 1 | 1 |
| | a. Rule | 088183 | 2 | 2 | 2 |
| | b. Roller | 086390 | 4 | 4 | 4 |
| | c. Cable | 087756 | 1 | 1 | 1 |
| | d. Tesamoll Tape, 30 mm x 6mm,White | 070304 | 1 | 1 | 1 |
| | e. Bracket | 040808 | 2 | 2 | 2 |
| | f. Cover for Bucky | 087521 | 1 | 1 | 1 |
| | g. Bracket, Left | 086579 | 1 | 1 | 1 |
| | h. Bracket, Right | 086580 | 1 | 1 | 1 |
| | i. Spring | 053548 | 2 | 2 | 2 |
| | j. Bracket | 088125 | 1 | 1 | 1 |
| | k. Knurled Screw | 055278 | 1 | 1 | 1 |
| 12 | Accessories | | | | |
| | a. Lateral Cassette Holder | 097789 | 1 | 1 | 1 |
| | b. Compression Immobilizing Device | 097787 | 1 | 1 | 1 |
| | c. Patient Hand Grip | 097786 | 2 | 2 | 2 |
| | d. Radiolucent Bands | 448606 | 2 | 2 | 2 |
| 13 | X- Ray Transport Device | 083282 | 1 | 1 | 1 |
| | a. Base Frame | 087836 | 1 | 1 | 1 |
| | b. Support | 087839 | 1 | 1 | 1 |
| | c. Support | 087840 | 1 | 1 | 1 |
| | d. Cover | 087841 | 1 | 1 | 1 |
| | e. Spring plug | 448703 | 4 | 4 | 4 |
| | f. Ring | 038900 | 4 | 4 | 4 |
| | g. Column Support | 087842 | 2 | 2 | 2 |
| | h. Base | 087865 | 1 | 1 | 1 |

(continued) Appendix D. Equipment Parts and Accessories List

6525-01-384-9296

X-Ray Apparatus, Model LCROKS

| ITEM NO | Part or Accessory Description | Part or Accessory Number | Basic Issue | Operator | Repairer |
|---------|-------------------------------|--------------------------|-------------|----------|----------|
| | i. Stirrup | 088055 | 2 | 2 | 2 |
| | j. Z- Locking Pliers | 083534 | 1 | 1 | 1 |
| | k. Bracket | 088067 | 1 | 1 | 1 |
| 14 | Bucky Assembly | 15669, 15676, 15682 | 1 | 1 | 1 |
| 15 | Cassette Size Sensing Tray | 15665 | 1 | 1 | 1 |

(continued) Appendix D. Equipment Parts and Accessories List

6530-00-926-2151

Sterilizer, Surgical Dressing 16x36 in., Model M-138

| ITEM NO | Part or Accessory Description | Part or Accessory Number | Basic Issue | Operator | Repairer |
|---------|-----------------------------------|--------------------------|-------------|----------|----------|
| 1 | Operators/Service Manual | TM 8-6530-004-24&P | 2 | 1 | 1 |
| 2 | Door Gasket | B300-264-23/2 | 2 | 1 | 1 |
| 3 | Gasket, Heating Element, Asbestos | C300-249-46/109 | 1 | 0 | 0 |
| 4 | Heating Element, Electrical | C300-249-48/109 | 3 | 0 | 0 |
| 5 | Sleeves, Heating Element | C300-250-50/109 | 1 | 1 | 1 |
| 6 | Washer, Flat, Teflon | 5310-00-926-9399 | 2 | 1 | 1 |
| 7 | Scraper Assembly | C300-906-139 | 1 | 1 | 1 |
| 8 | Shelf, Bottom | C300-906-119 | 1 | 1 | 1 |
| 9 | Shelf | C300-906-123 | 6 | 6 | 6 |

Appendix E. Tools and TMDE Code Listing for Maintenance Allocation Charts

| Reference Code | Item / Nomenclature | LIN | NSN | Model |
|----------------|--------------------------------------------------------------------------------|------------------|--------------------------------------|-----------------------------------------------|
| 01 | Tool Kit, Medical Equipment Maintenance Repairer | W45334 | 5180-00-611-7923 | Individual GSA Item (8001) |
| 02 | Tool Set, Medical Equipment Maintenance Unit Level | W45197 | 5180-01-483-1431 | Unit Level (alt: Org Maint) (8004) |
| | | | | |
| 04 | Tool Set, Medical Equipment and Maintenance General Support Level | T24386 | 6545-01-555-8683 | Shop Set Med Log Med Equip Maintenance (8007) |
| 05 | X Ray Calibration & Verification System | CO5856 | 6525-01-589-7774 | PN 1505061 |
| 06 | Calibrator Analyzer Hospital Equipment Computer Laptop 7010-01-502-5490 | C61523 | 6515-01-541-2864 | PN 2647336 |
| 07 | Analyzer Anesthetic Gas (AAG) | A00098 | 6630-01-530-7959 | Riken F1-21 |
| | | | | |
| 10 | Analyzer NonInvasive Blood Pressure (ANBP) | A27104 | 6515-01-449-1423 | PN 4515-0502 |
| 11 | Infusion Pump Analyzer (IPA) | J24245 | 6515-01-536-2122 | PN 2250063 |
| 12 | Defibrillator & Transcutaneous Pacemaker Analyzer (DATPA) | A83433 | 6515-01-575-2825 | PN 3366794 |
| 13 | Counter, Electronic Digital B64 | C19266 | 6515-01-406-7390 | |
| | | | | |
| 15 | Simulator, Medical Function | S56720 | 6515-01-548-3352 | MPS450 Army Medical Kit |
| 16 | Signal Generator B64 | S48323 | 6625-01-276-9421 | SG-1288/G |
| | | | | |
| 18 | Meter Foot Candle | M38443 | 6695-01-574-8935 | PN 07-621 |
| | | | | |
| 20 | Multimeter, AN/PSM-45A (CTA item) B64 Multimeter, AN/USM-486 (CTA item) B64 | M60449 M23954 | 6625-01-265-6000 6625-01-145-2430 | 27/FM 8050A |
| 21 | Oscilloscope, Digital Hand Held (ODHH) | P43667 | 6625-01-598-3865 | PN 199XRay |
| 22 | Wattmeter Ultrasonic | R95994 | 6695-01-593-8218 | USP-50 |
| 24 | Simulator, Pulse Oximetry (SPO) | S57953 | 6515-01-541-0432 | INDEX 2XLFE |
| 25 | Tachometer Stroboscopic | T07421 | 6625-01-550-3339 | |
| | | | | |
| 28 | Test Set, Electrosurgical Apparatus | T90883 | 6515-01-564-8554 | PN 2649769 |
| 29 | Tester, Current Leakage | T61791 | 6625-01-577-6744 | PN 3367232 |
| 31 | Tester, Ventilator Portable (TVP) | T05633 | 6515-01-449-1421 | PN 14500 |

INDEX OF DA SB 8-75-S8

| PARAGRAPH SUBJECT | PAGE |
|---------------------------------------------------------------------|----------------|
| Customers ----- | 2-2 |
| Equipment Parts and Accessories List----- | D-1 |
| Equipment Support Plans Availability ----- | 2-2 |
| Index ----- | IN-1 |
| Maintenance Allocation Chart (MAC) ----- | 1-3 and C-1 |
| Operator Level PMCS ----- | 1-1 |
| Operator PMCS - Generic Standards ----- | A-1 |
| Preventive Maintenance Checks and Services (PMCS) ----- | 1-1 |
| Repairer PMCS - Generic Standards----- | B-1 |
| Repairer Level PMCS ----- | 1-2 |
| Requirements----- | 2-1 |
| Responsibility----- | 2-1 |
| Supportability Strategy Overview ----- | 2-1 |
| Tools and TMDE Code Listing for Maintenance Allocation Charts ----- | E-1 |

By Order of the Secretary of the Army:

Official:



GERALD B. O'KEEFE

*Administrative Assistant to the
Secretary of the Army*

1322409

RAYMOND T. ODIERNO
*General, United States Army
Chief of Staff*

Distribution:

This publication is available through Electronic Media Only (EMO) by accessing the Army Knowledge Online (AKO) website and is intended for the use of command levels C, D, and E for the Active Army, the Army National Guard/Army National Guard of the United States, and the US Army Reserve.

PIN: 064785-000