

Operator Preventive Maintenance Checks and Services

6525-01-303-6235

X-Ray Process Machine, Model AFP14X3MIL

[B-Before Operation, D-During Operation, A-After Operation, Q-Quarterly, and S-Semiannually]

ITEM NO	INTERVAL	ITEM TO BE INSPECTED AND PROCEDURE	IS NOT MISSION CAPABLE IF:
1	<p>B</p> <p>B</p> <p>B, D</p> <p>B</p> <p>B</p> <p>B</p>	<p>X-Ray Processor</p> <p>a. Conduct an inventory to ensure that the items listed on the Equipment Parts and Accessories List are on hand.</p> <p>b. Inspect the processor for obvious signs of damage such as cracks, dents, leaks or broken components.</p> <p>c. Perform the "Periodic Maintenance Checks" as directed in the manufacturer's literature. (1) Weekly (2) Monthly</p> <p>d. Conduct the "Algae Control" procedure as directed in the manufacturer's literature.</p> <p>e. Conduct the "Whenever Chemistry is Changed" procedure as directed in the manufacturer's literature.</p> <p>f. Lubricate the processor as directed in the manufacturer's literature. (1) Weekly (2) Monthly (3) Quarterly (4) Semiannually (5) Every five years</p>	<p>Missing components or accessories prevent the operation of the processor.</p> <p>The damage prevents operation of the processor.</p> <p>The processor does not meet the "Periodic Maintenance Checks."</p> <p>The solutions are contaminated.</p> <p>The processor does not operate.</p>