

## Repairer Preventive Maintenance Checks and Services

6630-01-376-9823  
Analyzer, Clinical Chemistry, DT 60

[M-Monthly, Q-Quarterly, S-Semiannually, and A-Annually]

ITEM NO	INTERVAL	ITEM TO BE INSPECTED AND PROCEDURE	IS NOT MISSION CAPABLE IF:
1	S	<p><b>Analyzer, Clinical Chemistry</b></p> <p>a. Conduct an inventory to ensure that the items listed on the Equipment Parts and Accessories List are on hand.</p> <p>b. Inspect the unit for dust, dirt, damage, or excessively worn components.</p> <p>c. Verify the installation of the equipment according to the Operator Preventative Maintenance Checks and Services.</p> <p>d. Perform the procedures listed under Item 3, "Operating Instructions" in the Operator Preventative Maintenance Checks and Services.</p> <p>e. Verify the "Calibration" procedure according to the Operator Preventative Maintenance Checks and Services.</p> <p>f. Perform the "Instrument Care and Cleaning" procedures according to the Operator Preventative Maintenance Checks and Services.</p> <p>g. Verify electrical safety.</p> <p>h. Update the Medical Equipment Verification/Certification label (DD Form 2361).</p>	<p>Missing components or accessories prevent the operation of the DT60 system.</p> <p>Damage or deteriorated components prevent the operation of the unit.</p> <p>The analyzer fails any of the electrical safety tests.</p> <p>The unit has not been verified within the last six (6) months.</p>