

(continued) Appendix C. Maintenance Allocation Chart

6515-01-453-4003  
Defibrillator ECG Monitor/Recorder, LIFEPAK 10

[C – Operator or Crew; O – Unit Maintenance; DS – Direct Support Maintenance; GS – General Support Maintenance; D – Depot Maintenance

(1) GROUP NUMBER	(2) ASSEMBLY GROUP	(3) MAINTENANCE FUNCTION	(4) MAINTENANCE LEVEL	(5) TOOLS AND EQUIPMENT	(6) REMARKS
00	Defibrillator	Inspect	C		No equipment required for Crew level test. Defibrillator Programming Key required.
		Service	C		
		Test	C	01, 12, 20	
		Calibrate	O	01, 12, 20	
		Repair	O	01, 12, 20	
		Safety	O	12	
01	Batteries	Inspect	C		
		Service	C		
		Test	C		
		Replace	C		
02	Battery Connector Pins	Inspect	C		
		Service	C		
		Test	O	01, 20	
		Replace	O	01	
03	Recorder	Inspect	C		Cotton swab and isopropyl alcohol required.  Beware static sensitive components.
		Service	O		
		Test	O	12	
		Replace	O	01	
		Repair	O	01	
04	CRT Display	Inspect	C		No equipment required for Crew level test. Read all manufacture warnings before beginning.
		Service	O		
		Test	C	12	
		Replace	O	01, 20	
		Repair	O	01	
05	Defibrillator Paddles	Inspect	C		
		Service	C		
		Test	C	12	
		Replace	C		
06	ECG and Pacer Leads	Inspect	C		No equipment required for Crew level test.
		Service	C		
		Test	C	12	
		Replace	C		

(continued) Appendix C. Maintenance Allocation Chart

6515-01-453-4003  
Defibrillator ECG Monitor/Recorder, LIFEPAK 10

[C – Operator or Crew; O – Unit Maintenance; DS – Direct Support Maintenance; GS – General Support Maintenance; D – Depot Maintenance

(1) GROUP NUMBER	(2) ASSEMBLY GROUP	(3) MAINTENANCE FUNCTION	(4) MAINTENANCE LEVEL	(5) TOOLS AND EQUIPMENT	(6) REMARKS
07	Battery Support System	Inspect Test Service Repair Overhaul	C O O O O	01, 12, 20  01 12	
08	Printed Circuit Boards	Inspect Test Service Replace Repair	O O O O D	01, 12, 20 01, 12, 20 01, 12, 20 01, 12, 20 03, 12, 20	
09	Auxillary Power Module	Inspect Test Service Replace	C O O O	01, 12, 20  01	