

**ELECTROCARDIOGRAPH CP200 TOC**

**WELCH ALLYN**

**MODEL #: CP2A-1EM**

**Distributor's PHONE: (503) 530-7856 FAX: (503) 526-4919**

**Website: <http://www.welchallyn.com/products/en-us/x-11-ac-100-000000001026.htm>**

**AAC "L" ITEM NSN 6515-01-505-3025**

**UA's: F308, N308, N513, N713**

**USAMMA POC: Customer Relations Management, 301-619-4301/4316, DSN 343-4301/4316  
UNIT PRICE: \$2,920.64**

**System Description:** A multichannel electrocardiograph that is flexible enough to work in any setting. The CP200 is a powerful, easy-to-use device featuring an intuitive design and color display that make taking ECGs (electrocardiograms) easier.

**Features:**

- An interpretive 12-lead resting multichannel ECG that provides accurate measurements
- Full alphanumeric keyboard and color LCD screen to enter patient data quickly
- High-resolution thermal printer generates standard full-sized reports
- Help feature to access patient education tools as well as answer system operation questions
- Frequency range of 0.3 to 150 Hz
- Digital sampling rate of greater than 1,000 samples/second/channel
- Prints up to 100 continuous ECGs per battery charge




**ELECTROCARDIOGRAPH CP200  
MODEL # CP200 KIT W/ CASE  
NSN 6515-01-505-3025**



**CABLE, PATIENT DIGITAL ECG  
PART # 400293  
NSN 6515-01-560-9895**




**10-LEAD REPLACEABLE LEAD SET  
PART # 401129  
NSN 6515-01-561-0335**



**PAPER, THERMAL  
PART # 94018-0000  
NSN 6515-01-552-7128**



**ALLIGATOR CLIP  
PART # 58581-0000  
NSN 6515-01-505-3043**



**TAB ELECTRODES  
PART # 45008-0000  
NSN 6515-01-505-3037**

**ELECTROCARDIOGRAPH CP200  
PHYSICAL & ENVIRONMENTAL SPECIFICATIONS**

<b>Dimensions</b>	<b>Height</b>	<b>6.2 inches</b>
	<b>Width</b>	<b>16.2 inches</b>
	<b>Depth</b>	<b>15.6 inches</b>
<b>Cube</b>		<b>0.9 cu ft</b>
<b>Weight</b>	<b>Including battery</b>	<b>11.6 lb</b>
<b>Operating Temperature</b>		<b>32° to 104° F</b>
<b>Relative humidity</b>		<b>15 – 95% noncondensing (30 – 70% for printing)</b>
<b>Storage Temperature</b>		<b>-4°F to 120°F</b>
<b>Relative humidity</b>		<b>15 – 95% noncondensing</b>
<b>Power Requirements</b>		<b>100-240 VAC 50/60Hz</b>
<b>Amps</b>		<b>2.2 A maximum</b>
<b>Watts</b>		<b>13.2 W</b>
<b>Battery</b>	<b>12 hour recharging</b>	<b>Rechargeable lead acid gel, 6VDC, 5 AH (ampere-hour)</b>
<b>Paper</b>	<b>200 sheets/PG</b>	<b>8.25 x 11 inches, Z-fold thermal</b>

**ELECTROCARDIOGRAPH CP200  
(START-UP)**

NSN	NOTES	Nomenclature	PN	QTY	UI	Unit Price	Total Price	DIST/MFG	Shelf Life, months	Refrigeration Y/N
6515015053025	Included w/ NSN	Electrocardiograph CP200 kit w/case	CP2A-1EM	1	EA	\$2,720.64	\$2,720.64	Welch Allyn	indef	N
6515015609895		Cable, patient digital ECG, 10-lead, replaceable, AHA	400293	1	EA	Included w/ 3025 NSN	Included w/ 3025 NSN	Welch Allyn	indef	N
6515015610335		10-lead replaceable lead set, AHA, banana clip, 4 limb leads, 6 precordial (for use w/ 400293)	401129	1	SE	Included w/ 3025 NSN	Included w/ 3025 NSN	Welch Allyn	indef	N
	Use resupply	Tab electrodes, disposable	45008	1	100/PG	Included w/ 3025 NSN	Included w/ 3025 NSN	Welch Allyn	indef	N
	Use resupply	Paper, thermal, 200 sheets	94018	1	RM	Included w/ 3025 NSN	Included w/ 3025 NSN	Welch Allyn	60	N
6540015009058		AC power cord	761076-0	1	EA	Included w/ 3025 NSN	Included w/ 3025 NSN	Welch Allyn	indef	N
6515015609958		Operator manual, training guide & physician's guide, CD	401151	1	EA	Included w/ 3025 NSN	Included w/ 3025 NSN	Welch Allyn	indef	N
7610015436723		Service manual, CD	704340	1	EA	Included w/ 3025 NSN	Included w/ 3025 NSN	Welch Allyn	indef	N
6515015611869		Case, storm, hard-sided	IM2720	1	EA	Included w/ 3025 NSN	Included w/ 3025 NSN	Welch Allyn	indef	N
NO NSN REQUIRED		Case, soft-side/Dust Cover	701586	1	EA	Included w/ 3025 NSN	Included w/ 3025 NSN	Welch Allyn	indef	N

**NOTE: [CLICK FOR SERVICE MANUAL /](http://www.welchallyn.com/documents/Cardiopulmonary/Electrocardiographs/CP200/ServiceManual_2007-8-3_CP100&CP200.pdf)**

**[http://www.welchallyn.com/documents/Cardiopulmonary/Electrocardiographs/CP200/ServiceManual\\_2007-8-3\\_CP100&CP200.pdf](http://www.welchallyn.com/documents/Cardiopulmonary/Electrocardiographs/CP200/ServiceManual_2007-8-3_CP100&CP200.pdf)**