

CONCENTRATOR OXYGEN SYSTEM PORTABLE W/CASE TOC

AIRSEP NEWLIFE ELITE

MODEL # AS119-1D

Manufacturer's PHONE: (716) 691-0202 FAX: (716) 691 0174

Website: <http://www.airsep.com/medical/newlifeelite.html>

AAC "J" ITEM NSN 6515-01-552-9967

UA's: 267B, 901A, 905A, 912A

USAMMA POC: Customer Relations Management 301-619-4301/4316, DSN 343-4301/4316

UNIT PRICE: \$2,372.00

System Description: Produces from 1 to 3 liters per minute of therapeutic oxygen, purity levels of up to 92%, 4 LPM at 89 %, and 5 LPM at 82 %. This system has built-in alarms for power failure and high/low pressure. Unit includes filters to remove dust, bacteria and other particulates. Unit provides for the attachment of an in-line humidifier. The gas mixture is delivered to the patient for use by producing oxygen rich gas mixture by drawing in room air and extracting the oxygen. This is a large, easy-to-read flow meter with back-lines for accurate flow settings. The built-in handles make it easy to move or lift. The large casters glide easily over virtually any flooring surface.



**CONCENTRATOR OXYGEN SYSTEM
PORTABLE W/CASE
MODEL # AS119-1D
NSN 6515-01-552-9967**

CONCENTRATOR OXYGEN SYSTEM PORTABLE W/CASE PHYSICAL AND ENVIRONMENTAL SPECIFICATIONS

Dimensions	Height	27.5 inches
	Width	16.7 inches
	Depth	14.5 inches
Weight		54 lbs
Cube		3.85 cu ft
Operating Temperature Range		40° to 105 ° F
Relative Humidity (storage)		4° to 140° F
Relative Humidity (in use)		0 to 95%
Power Requirements		120 VAC 50/60 Hz
Amps		4 amps
Watts		480 watts
Performance:		92% at 1-3 LPM, 89% at 1-4 LPM, 82% at 1-5 LPM
Alarms: Power failure, High and low pressure, Temperature, Battery test		
Low purity (Oxygen monitor with optional Eco Check.)		

CONCENTRATOR OXYGEN SYSTEM PORTABLE W/CASE (START-UP)

NSN	Nomenclature	PN	QTY	UI	Unit Price	Total Price	DIST/MFG	Shelf Life, months	Refrig Y/N	Notes
6515015529967	Oxygen System Is The Forward Surg. Team Port W/Case	AS119-1D	1	EA	\$2,372.00	\$2,372.00	AirSep	Indef.	N	
6515014344629	Concentrator Oxygen 1-5 LPM 115/230V 50/60Hz Built-In Alarms	AS005-4	1	EA	Included W/9967 NSN	Included W/9967 NSN	AirSep	N/A	N	
	Humidifier Bottle	HU003-1	1	EA	Included W/9967 NSN	Included W/9967 NSN	AirSep	N/A	N	Mounted on the Oxygen Concentrator
6515015527414	Case Oxygen Concentrator Padded For 6515014344629	CA253-1	1	EA	Included W/9967 NSN	Included W/9967 NSN	AirSep	N/A	N	
6515014724663	Filter Bacteria End Item ID 6515-01-327-6798	F1001-1	1	EA	Included W/9967 NSN	Included W/9967 NSN	AirSep	N/A	N	Installed inside of the concentrator; replace after 50,000 hours of use.
7610015379355	Manual, Technical/Service Oxygen Concentrator 6515014344629	MN105-1	2	EA	Included W/9967 NSN	Included W/9967 NSN	AirSep	N/A	N	Must request two manuals, otherwise just one will ship