

Maintenance Allocation Chart

6630-01-316-5085

Centrifugal Hematology Analyzer System with QBC II Reader, Model 4477 and QBC Centrifuge, Model 4207

[C – Operator or Crew; O – Unit Maintenance; DS – Direct Support Maintenance; GS – General Support Maintenance; D – Depot Maintenance]

(1) GROUP NUMBER	(2) ASSEMBLY GROUP	(3) MAINTENANCE FUNCTION	(4) MAINTENANCE LEVEL	(5) TOOLS AND EQUIPMENT	(6) REMARKS
00	QBC II Reader Model 4477	Install	C		Requires calibration check tubes – venous and capillary modes.
		Inspect	C		
		Service	C	01, 20, 29	
		Test	O	01, 20, 29	
		Calibrate	O	01	
		Repair	O	01, 02, 20, 29	
		Refurbish	D	01, 02, 03, 20, 29	
01	QBC Centrifuge Model 4207	Install	C		Requires a stopwatch, 1-second accuracy.
		Inspect	C		
		Service	C		
		Test	O	01, 25, 29	
		Calibrate	O	01, 25, 21	
		Repair	O	01, 02, 25, 21, 29	
		Refurbish	O	01, 02, 03, 25, 21, 29	
02	Work Station	Install	C		
		Inspect	C		
		Replace	C		
03	Voltage Step-down Transformer	Install	C		
		Inspect	C		
		Test	O	01, 20	
		Replace	O		