

ORTHO BLOOD BANK GEL TOC

Distributed by: GOVERNMENT SCIENTIFIC SOURCE
Manufactured by: ORTHO CLINICAL DIAGNOSTIC INC
Distributor's MODEL# FSCOKMTSPACKAGEEA
Manufacturer's MODEL# DSTDSI120806100EA
Distributor's PHONE: (703) 734-1805 FAX: (703) 734-3321
Manufacturer's PHONE: (585) 453-3400 FAX: (585) 453-3660
Distributor's Website: <https://mozart.govsci.com/>
Manufacturer's Website: www.orthoclinical.com

ACC "J" ITEM NSN 6640-01-549-8235

UA's: N304, N504, N704

USAMMA POC: Customer Relations Management 301-619-4301/4316, DSN 343-4301/4316
UNIT PRICE: \$12,005.57

System Description: Ortho provides the tools you need to manage your pre-transfusion testing with dramatically improved workflow efficiency and standardization. Our product offering includes the ID-Micro Typing System™ (ID-MTS) Gel Test™, with flexible automation to meet your changing needs: ABO & Rh typing, antibody screening & identification, and antigen typing.

Because the ID-Micro Typing System Gel Test requires fewer procedural steps, it's easier to perform and more cost effective than other serological methods. Plus, reduced handling means fewer operator errors and more objective interpretation of results.

Testing is performed in a prefilled card containing gel mixed with the appropriate reagent. Simply pipette, incubate, spin and read. Results are ready to report, without the need for cell washing, red cell resuspension or tube shaking. The result is optimized utilization of personnel, simplified training and increased cross-training opportunities for improved productivity and greater workflow efficiency.



**ORTHO BLOOD BANK GEL
MODEL# FSCOKMTSPACKAGEEA
NSN 6640-01-549-8235**

ORTHO BLOOD BANK GEL PHYSICAL AND ENVIRONMENTAL SPECIFICATIONS

Dimensions of Incubator	Height	4.25 inches
	Width	10.75 inches
	Depth	12.5 inches
	Cube	0.33 cu ft
Weight	Free Standing	11 lbs
Operating Temperature Range		50° F to 86° F
Relative Humidity (in use)		Up to 90%
Power Requirements		110 VAC 50/60 Hz
Amps		1 amp
Watts		110 watts

Dimensions of Centrifuge	Height	5.75 inches
	Width	13.50 inches
	Depth	19.25 inches
Cube		0.86 cu ft
Weight	Free Standing	14 lbs
Operating Temperature Range		41° F to 104° F
Relative Humidity (in use)		Up to 90%
Power Requirements		110-120VAC 50/60 Hz
Amps		3 amps
Watts		330 watts

Caution: The main supply voltage should not exceed $\pm 10\%$ of the nominal voltage

ORTHO BLOOD BANK GEL (START-UP)

NSN	Notes	Nomenclature	PN	QTY	UI	Unit Price	Total Price	DIST/MFG	Shelf Life, months	Refrig Y/N
6640015498235	ECAT	Ortho Blood Bank Gel System	FSCOKMTSPACK AGEEA	1	EA	\$12,005.57	\$12,005.57	Government Scientific Source	indef.	N
6640014953487		Blood Bank Incubator gel System 110v AC 50/60Hz 12.5x10.75x4.25in	MTS9680	1	EA	Included W/8235 NSN	Included W/8235 NSN	Ortho Clinical Diagnostic Inc	N/A	N
7610015436723		Operators Manual	50000616	1	EA	Included W/8235 NSN	Included W/8235 NSN	Ortho Clinical Diagnostic Inc	N/A	N
6640014953309		Centrifuge laboratory 110/120VAC 60hz single phase timer preprogramed	MT515060	1	EA	Included W/8235 NSN	Included W/8235 NSN	Ortho Clinical Diagnostic Inc	N/A	N
7610015437003		Operators Manual	50000615	1	EA	Included W/8235 NSN	Included W/8235 NSN	Ortho Clinical Diagnostic Inc	N/A	N
6515014937272		Table work Ortho MTS plastic 14"x14"x3" use with Ortho blood bank system	MTS-9650	1	EA	Included W/8235 NSN	Included W/8235 NSN	Ortho Clinical Diagnostic Inc	N/A	N
	ECAT Only	Hardened case to fit above items	DSTD5I120806100 EA	1	EA	Included W/8235 NSN	Included W/8235 NSN	Government Scientific Source	N/A	N