

**MEDICAL EQUIPMENT
SET PRIMARY GYN
CARE**



**Support and Consumables
Handbook**

6545-01-527-8013

UA 301A

LIN M29701

Semper Responsus

REVISED

January 2013

U.S. Army Medical Materiel Agency

The information contained in this manual was prepared by the

U.S. Army Medical Materiel Agency
Office of PM, Medical Devices
Acquisition Information Division
693 Neiman Street
Fort Detrick, MD 21702-5001

This manual is a living document and may contain inaccuracies due to changes in mission requirements or equipment updates. To report corrections or additions please contact:

Customer Relationship Management
301-619-1288/4301, DSN 343-1288/4301, or at:

USAMMACRM@amedd.army.mil

**For additional equipment item information, visit the
USAMMA Medical Equipment Literature section on the web at
<<http://www.usamma.amedd.army.mil/med equip lit cdroms.cfm>>**

**Review the Master Table of Contents for all of the
published USAMMA Medical Equipment CDs at
<<http://www.usamma.amedd.army.mil/assets/docs/Master%20Table%20of%20Contents.pdf>>**

**You can also order a copy of each published CD by completing and submitting the form at
<<http://www.usamma.amedd.army.mil/CDROM Requests form.cfm>>**

**NOTE: CLICK ON TITLE
FOR DIRECT LINK**

TABLE OF CONTENTS

SECTION	DEVICE	NSN	PAGE
1	<u>LIGHT FLOOR GENERAL EXAM GREEN SERIES</u>	(L) 6530-01-589-4229	1
2	<u>TABLE EXAMINING AND TREATMENT W STIRRUPS</u>	(L) 6530-01-531-2896	5

Please provide all comments at the website link listed below.

https://ice.disa.mil/index.cfm?fa=card&site_id=750&service_provider_id=102377

Select: EQUIPMENT PUBLICATIONS BUTTON AT TOP OF FORM

LIGHT FLOOR GENERAL EXAM GREEN SERIES

[TOC](#)

MANUFACTURER: WELCH ALLYN

Manufacturer Model# GS 300

Manufacturer's PHONE: (315) 685-4560 FAX: (315) 683-3361

Manufacturer's Website: <http://www.welchallyn.com>

Item Web Link:

<http://www.welchallyn.com/products/en-us/x-19-br-96-1255625713648.htm>

AAC 'L' ITEM NSN 6530-01-589-4229

UAs: 301A, 324A, 623A, 625A, N305, N312, N315, N316, N413, N432, N436, N503, N513, N703, N713, P301, P307, P308, P309, P310, P311, P314, P433

USAMMA POC: Customer Relations Management 301-619-1288/4301 - DSN 343-1288/4301

UNIT PRICE: \$304.15

System Description: This general exam light is a perfect combination of intense light output and cool operation, all delivered in a compact design and with rugged performance.

- Features:**
- Touchless On/Off
 - Multiple mounting options for floor space maximization
 - Low cost of ownership due to elimination of bulb replacement and lower energy consumption
 - Easy Assembly



LIGHT FLOOR GENERAL EXAM GREEN SERIES
PART # 44400
NSN 6530015894229

**LIGHT FLOOR GENERAL EXAM GREEN SERIES
PHYSICAL & ENVIRONMENTAL SPECIFICATIONS**

Item Dimensions	Height (Inches)	57.59
	Width (Inches)	20.87
	Length (Inches)	20.87
	Cube (Cubic Feet)	14.516
	Weight	18.40 lbs
Ship Dimensions	Height (Inches)	0
	Width (Inches)	0
	Length (Inches)	0
	Cube (Cubic Feet)	0
	Weight	
Storage Temp		-4 to 120 °F
Operating Temp		50 to 95 °F
Storage Humidity		
Operating Humidity		
Amps		0.35 (max)
Watts		40.25
Power Requirements	100/240 VAC, 50/60 Hz	
Notes		

**LIGHT FLOOR GENERAL EXAM GREEN SERIES
ACCESSORIES & CONSUMABLES - STARTUP**

NSN	Notes	Nomenclature	MFR PN Welch Allyn	DIST PN	QTY	UI	Unit Price	Total Price	Shelf Life Mon.
6530015894229		Light Floor General Exam Green Series	44400		1	EA	\$304.15	\$304.15	N/A
NSN Not Required		GS300 W/Control Box	406790		1	EA	Incl w/NSN 4229	Incl w/NSN 4229	N/A
NSN Not Required		Literature, Operator	713317		1	EA	Incl w/NSN 4229	Incl w/NSN 4229	N/A
NSN Not Required		Power Cord 8 Foot Long	714030		1	EA	Incl w/NSN 4229	Incl w/NSN 4229	N/A
NSN Not Required		Stand Mobile GS Series	405515		1	EA	Incl w/NSN 4229	Incl w/NSN 4229	N/A

**LIGHT FLOOR GENERAL EXAM GREEN SERIES
ACCESSORIES & CONSUMABLES - RESUPPLY**

NSN	Notes	Nomenclature	MFR PN Welch Allyn	DIST PN	QTY	UI	Unit Price	Total Price	Shelf Life Mon.
NSN Not Required		Power Cord 8 Foot Long	714030		1	EA	\$15.00	\$15.00	N/A

TABLE EXAMINING AND TREATMENT W STIRRUPS

[TOC](#)

MANUFACTURER: NORTH COAST OUTFITTERS

Manufacturer Model# SR901-OBGYN

Manufacturer's PHONE: (631) 727-5580 FAX: (631) 727-5584

Manufacturer's Website: <http://www.CharliesHorse.com>

Item Web Link:

<http://www.CharliesHorse.com/html/comm/sr901rt.html>

AAC 'L' ITEM NSN 6530-01-531-2896


UAs: 301A, N316

USAMMA POC: Customer Relations Management 301-619-1288/4301 - DSN 343-1288/4301


UNIT PRICE: \$6,518.90

System Description: The SR901 is truly unique. This Multi-Purpose Utility Table can be utilized to provide a functional operating table or OB GYN table, which is lightweight and compact to far forward remote areas allowing doctors to stabilize casualties and operate when necessary prior to transport. This table assists nurses in customizing the table to the doctor's specific needs for treatment in seconds. End item includes Accessory Kit B (stirrups). Unit has a capacity of 500 lbs.

- Features:**
- Mounted items can be repositioned horizontally along the side of table
 - Table has Trendelenburg capability, even while patient is on the table
 - Patient can be transitioned to the right or left
 - Each end of table can be raised a total of 6 inches in 1-inch intervals
 - 1 inch of vertical micro adjustments available for leveling
 - Patient can be lowered 4 inches
 - IV poles can be mounted at any corner
 - Accessory Kits come in individual carry bags



**TABLE EXAMINING WITH STIRRUPS (SHOWN WITH ACCESSORY KITS A, B & C)
PART # SR901-OBGYN
NSN 6530015312896**



**TABLE, EXAMINING & TREATMENT WITH STIRRUPS (IN OPEN CARRYING BAG)
PART # SR901-OBGYN
NSN 6530015312896**



CARRYBAG LEVEL 1 - 1.5 INCH SUPPORT TUBES AND 1 INCH LEG SUPPORTS



CARRYBAG LEVEL 2 - LEGS AND STRETCHER CLAMP ASSEMBLIES



CARRYBAG LEVEL 3 - LOWER SHELF, MESS SHELIVING (2) AND FLUID CAPTURE

**TABLE EXAMINING AND TREATMENT W STIRRUPS
PHYSICAL & ENVIRONMENTAL SPECIFICATIONS**

Item Dimensions	Height (Inches)	84
	Width (Inches)	33
	Length (Inches)	29
	Cube (Cubic Feet)	46.5208
	Weight	90 lbs
Ship Dimensions	Height (Inches)	48
	Width (Inches)	30
	Length (Inches)	10
	Cube (Cubic Feet)	8.3333
	Weight	110 lbs
Storage Temp		
Operating Temp		
Storage Humidity		
Operating Humidity		
Amps		
Watts		
Power Requirements	N/A	
Notes		

**TABLE EXAMINING AND TREATMENT W STIRRUPS
ACCESSORIES & CONSUMABLES - STARTUP**

NSN	Notes	Nomenclature	MFR PN North Coast Outfitters	DIST PN	QTY	UI	Unit Price	Total Price	Shelf Life Mon.
6530015312896		Table Examining and Treatment w Stirrups	SR901- OBGYN		1	EA	\$6,518.90	\$6,518.90	N/A
6530015294859	See Table Below for Kit Contents	Accessory, OR Table, Kit B	Accessory Kit B		1	KT	Incl w/NSN 2896	Incl w/NSN 2896	N/A

**TABLE EXAMINING AND TREATMENT W STIRRUPS
ACCESSORIES & CONSUMABLES - RESUPPLY**

NSN	Notes	Nomenclature	MFR PN North Coast Outfitters	DIST PN	QTY	UI	Unit Price	Total Price	Shelf Life Mon.
6530015294859	See Table Below for Kit Contents	Accessory, OR Table, Kit B	Accessory Kit B		1	KT	\$2,275.00	\$2,275.00	N/A

**TABLE EXAMINING AND TREATMENT W STIRRUPS
ADDITIONAL NOTE:**

The SR901 has been packaged into the carry bag in the order it is assembled. This table is returned into the carry bag in the reverse order that it is assembled. If you take the time to review the instruction manual prior to assembly your efforts will be rewarded. Once you have reviewed the manual, set-up should be less than 20 minutes. Subsequent assemblies should take no more than 10 minutes.

For price quotes on individual items refer to user manual, call 631-727-5580, or e-mail
info@CharliesHorse.com.

TABLE EXAMINING AND TREATMENT W STIRRUPS
ITEM TABLES

Contents of Examining and Treatment Table, Accessory Kit A
NSN 6530-01-529-4842

Part	Part#	Qty/UI
Arm board assembly	NOC5067A	1/PG
Mayo tray assembly	NCO5083	1/PG
Accessory clamp assembly	NCO5090	2/PG
IV pole w/mount	6000228	2/PG

Contents of Arm Board Assembly Part# NOC5067A
(Contained in Kit A)

Part	Part#	Qty/UI
Tension knob	6000227	1/EA
Arm board shaft	NCO1199A	2/EA

Contents of Mayo Tray Assembly Part# NOC5083
(Contained in Kit A)

Part	Part#	Qty/UI
Tension knob	6000227	1/EA
Mayo tray shaft	NCO1199B	1/EA
Mayo tray shaft	6000226	1/EA

Contents of Accessory Clamp Assembly Part# NCO5090 (Contained in Kit A)

Part	Part#	Qty/UI
1.97 knurled knob 5/16-18 2.000"	60001000/156	2/EA

Contents of Examining and Treatment Table, Accessory Kit B NSN
6530-01-529-4859

Part	Part#	Qty/UI
Stirrup crossmember assembly	NCO5084	1/KT

Contents of Stirrup Crossmember Assembly Part# NCO5084 (Contained in
Kit B)

Part	Part#	Qty/UI
1.97 knurled knob 5/16-18 1.875"	6000100/153	4/PG
1.97 knurled knob 5/16-18 1.500"	6000100/154	2/PG
1.97 knurled knob 3/8-16 1.250"	6000101/155	2/PG
Stirrup	NCO1147	2/KT
Stirrup key	NCO0027	2/PG
Stirrup IV block	NCO0028	2/PG
Stirrup rail assembly	NCO5091	2/PG

Contents of Examining and Treatment Table, Accessory Kit C
NSN 6530-01-529-4857

Part	Part#	Qty/UI
Examination light assembly	NCO5086	2/PG

**Contents of Examination Light Assembly part# NCO5086
(Contained in Kit C)**

Part	Part#	Qty/UI
Examination light shaft	NCO1199C	2/PG

