

# **MMS PHARMACY - 84 BED EQUIPMENT SET**



## **Support and Consumables Handbook**

**6545-01-526-5130**

**UA N506**

**LIN M73254**

*Semper Responsus*

**ARCHIVED**

**August 2012**

**U.S. Army Medical Materiel Agency**



The information contained in this manual was prepared by the

U.S. Army Medical Materiel Agency  
Office of PM, Medical Devices  
Acquisition Information Division  
693 Neiman Street  
Fort Detrick, MD 21702-5001

This manual is a living document and may contain inaccuracies due to changes in mission requirements or equipment updates. To report corrections or additions please contact:

Customer Relationship Management  
301-619-1288/4301, DSN 343-1288/4301, or at:

[USAMMACRM@amedd.army.mil](mailto:USAMMACRM@amedd.army.mil)

**For additional equipment item information, visit the  
USAMMA Medical Equipment Literature section on the web at  
<<http://www.usamma.army.mil/med equip lit cdroms.cfm>>**

**Review the Master Table of Contents for all of the  
published USAMMA Medical Equipment CDs at  
<<http://www.usamma.army.mil/assets/docs/Master%20Table%20of%20Contents.pdf>>**

**You can also order a copy of each published CD by completing and submitting the form at  
<<http://www.usamma.army.mil/CDROM Requests form.cfm>>**

**NOTE: CLICK ON TITLE  
FOR DIRECT LINK**

## TABLE OF CONTENTS

SECTION	DEVICE	NSN	PAGE
<a href="#">1</a>	<a href="#">COUNTER TABLET AUTOMATIC</a>	<a href="#">(J) 6515-01-508-9046</a>	<a href="#">1</a>
<a href="#">2</a>	<a href="#">HOOD LAMINAR FLOW LAB</a>	<a href="#">(J) 6640-01-493-8582</a>	<a href="#">5</a>
<a href="#">3</a>	<a href="#">REFRIGERATOR-FREEZER</a>	<a href="#">(L) 4110-01-425-8009</a>	<a href="#">9</a>
<a href="#">4</a>	<a href="#">SINK UNIT SURGICAL SCRUB</a>	<a href="#">(L) 6530-01-572-6775</a>	<a href="#">13</a>

**Please provide all comments at the website link listed below.**

**[https://ice.disa.mil/index.cfm?fa=card&site\\_id=750&service\\_provider\\_id=102377](https://ice.disa.mil/index.cfm?fa=card&site_id=750&service_provider_id=102377)**

**Select: EQUIPMENT PUBLICATIONS BUTTON AT TOP OF FORM**



## COUNTER TABLET AUTOMATIC

[TOC](#)

**DISTRIBUTOR: AMERISOURCE BERGEN**

**MANUFACTURER: KIRBY LESTER**

**Distributor Model# KL15DF**

**Manufacturer Model# KL15DF**

**Distributor's PHONE: (847) 495-6685 FAX: (847) 808-3322**

**Manufacturer's PHONE: (203) 323-1499 FAX:**

**Distributor's Website: <http://www.AmerisourceBergen.com>**

**Manufacturer's Website: <http://www.kirbylester.com>**

**Item Web Link:**

**<http://www.kirbylester.com/pharmacy-products-KL15df.htm>**

**AAC 'J' ITEM NSN 6515-01-508-9046**

**UAs: N306, P506, P706**

**USAMMA POC: Customer Relations Management 301-619-1288/4301 - DSN 343-1288/4301**


**UNIT PRICE: \$5,242.38**

**System Description:** This automatic counter transfers tablets and capsules directly into the final container. It counts the average prescription and feeds it into a vial in less than 4 seconds. The counting technology used guarantees the most accurate counting available in any automatic system – 99.97%, or less than 3 errors in 10,000 tablets.

**Features:** Accuracy: 99.97% (less than 3 errors in 10,000 tablets)  
Maximum count: 500 tablets or capsules or 85 drams  
Maximum throughput: 1,500 parts per minute  
Maximum size counted: not applicable. Counts the largest known tablets and capsules



**COUNTER TABLET AUTOMATIC  
PART # KL15df  
NSN 6515015089046**



**COUNTER TABLET AUTOMATIC DEMONSTRATED USE  
PART # KL15df  
NSN 6515015089046**

**COUNTER TABLET AUTOMATIC  
PHYSICAL & ENVIRONMENTAL SPECIFICATIONS**

<b>Item Dimensions</b>	<b>Height (Inches)</b>	<b>19.5</b>
	<b>Width (Inches)</b>	<b>9</b>
	<b>Length (Inches)</b>	<b>17</b>
	<b>Cube (Cubic Feet)</b>	<b>1.7266</b>
	<b>Weight</b>	<b>21 lbs</b>
<b>Ship Dimensions</b>	<b>Height (Inches)</b>	
	<b>Width (Inches)</b>	
	<b>Length (Inches)</b>	
	<b>Cube (Cubic Feet)</b>	
	<b>Weight</b>	
<b>Storage Temp</b>		
<b>Operating Temp</b>		
<b>Storage Humidity</b>		
<b>Operating Humidity</b>		
<b>Amps</b>		<b>6</b>
<b>Watts</b>		<b>690</b>
<b>Power Requirements</b>	<b>115VAC/60Hz for US &amp; Canada, 220VAC/50Hz for overseas</b>	
<b>Notes</b>		

**COUNTER TABLET AUTOMATIC  
ACCESSORIES & CONSUMABLES - STARTUP**

NSN	Notes	Nomenclature	MFR PN Kirby Lester	DIST PN Amerisour ce Bergen	QTY	UI	Unit Price	Total Price	Shelf Life Mon.
6515015089046		Counter Tablet Automatic	KL15df	KL15df	1	EA	\$5,242.38	\$5,242.38	N/A
	Awaiting NSN	Channel Bottom	8CAI-0315		1	EA	Incl w/NSN 9046	Incl w/NSN 9046	N/A
	Awaiting NSN	Channel Top	D8AI- 0305-A		1	EA	Incl w/NSN 9046	Incl w/NSN 9046	N/A
	Awaiting NSN	Funnel	D8AM- 4001-A		1	EA	Incl w/NSN 9046	Incl w/NSN 9046	N/A
7610015437003		Literature Operation			1	EA	Incl w/NSN 9046	Incl w/NSN 9046	N/A

**COUNTER TABLET AUTOMATIC  
ACCESSORIES & CONSUMABLES - RESUPPLY**

NSN	Notes	Nomenclature	MFR PN Kirby Lester	DIST PN Amerisour ce Bergen	QTY	UI	Unit Price	Total Price	Shelf Life Mon.
	Awaiting NSN	Channel Bottom	8CAI-0315		1	EA	\$45.00	\$45.00	N/A
	Awaiting NSN	Channel Top	D8AI- 0305-A		1	EA	\$25.00	\$25.00	N/A
	Awaiting NSN	Funnel	D8AM- 4001-A		1	EA	\$90.00	\$90.00	N/A

## HOOD LAMINAR FLOW LAB

TOC

**MANUFACTURER: GERMFREE LABORATORIES INC.**

**Manufacturer Model# BZ4SS-RX**

**Manufacturer's PHONE: (386) 677-7742 FAX: (386) 677-2651**

**Manufacturer's Website: <http://www.germfree.com>**

**Item Web Link:**

**<http://97.74.181.168/pharmacy/laminar/BZ/>**

**AAC 'J' ITEM NSN 6640-01-493-8582**

**UAs: N306, P506, P706**

**USAMMA POC: Customer Relations Management 301-619-1288/4301 - DSN 343-1288/4301**

**UNIT PRICE: \$8,148.27**


**System Description:** The Horizontal Laminar Flow (BZ) Series Hood is designed for the handling of materials in a sterile working environment. The work area is continuously bathed with positive pressure ISO 5/ Class 100 horizontal laminar flow air that has passed through a High Efficiency Particulate Air (HEPA) filter. This filter removes organisms and particulates 0.3 micron in size with an efficiency of 99.99%. The laminar flow principle consists of moving individual streams of unidirectional, ultra-clean air along parallel lines with minimal turbulence. The HEPA filter is positioned in the rear of the work area and is protected by a removable, perforated metal diffuser. The BZ is constructed of quality all-welded stainless steel with a pharmaceutical grade finish. The fluorescent light housing is externally mounted to provide optimum lighting into the work area. The acrylic top and side panels allow for additional ambient light and visibility into the work area. Our SpillGuard feature is utilized in the rear and sides of the full-width work area. The CurveFlex forearm support feature provides operator comfort.

**Features:**

- Outer cabinet and work surface are welded stainless steel with a pharmaceutical grade finish
- High capacity motor/blower system with speed control to extend HEPA filter life
- Top mounted pre-filters are easily changed
- HEPA filter is full width and height (36"), ensuring unidirectional airflow and is easily changed from both sides
- Top mounted fluorescent lighting is isolated from the clean environment
- IV hanging bar with 12 hooks in work area
- Separate lighted power ON/OFF indicator switches for blower and lighting
- Electrical is 115Volt, 60 Hz (220/50-60 Hz also available)
- Ten foot power cord with molded grounded plug



**HOOD LAMINAR FLOW LAB  
PART # BZ4SS-RX  
NSN 6640014938582**



**LAMINAR HOOD AIR FLOW DIAGRAM**

**HOOD LAMINAR FLOW LAB  
PHYSICAL & ENVIRONMENTAL SPECIFICATIONS**

<b>Item Dimensions</b>	<b>Height (Inches)</b>	<b>49.5</b>
	<b>Width (Inches)</b>	<b>49</b>
	<b>Length (Inches)</b>	<b>34</b>
	<b>Cube (Cubic Feet)</b>	<b>47.724</b>
	<b>Weight</b>	<b>425 lbs</b>
<b>Ship Dimensions</b>	<b>Height (Inches)</b>	<b>60</b>
	<b>Width (Inches)</b>	<b>58</b>
	<b>Length (Inches)</b>	<b>41</b>
	<b>Cube (Cubic Feet)</b>	<b>82.5694</b>
	<b>Weight</b>	<b>515 lbs</b>
<b>Storage Temp</b>		
<b>Operating Temp</b>		
<b>Storage Humidity</b>		
<b>Operating Humidity</b>		
<b>Amps</b>		<b>12</b>
<b>Watts</b>		<b>1380</b>
<b>Power Requirements</b>	<b>115VAC 60HZ</b>	
<b>Notes</b>		

**HOOD LAMINAR FLOW LAB  
ACCESSORIES & CONSUMABLES - STARTUP**

NSN	Notes	Nomenclature	MFR PN GermFree Laboratorie s Inc.	DIST PN	QTY	UI	Unit Price	Total Price	Shelf Life Mon.
6640014938582		Hood Laminar Flow Lab	BZ4SS-RX		1	EA	\$8,148.27	\$8,148.27	N/A
	Awaiting NSN	Stainless Steel I.V. Bar	BZ-A9		1	EA	Incl w/NSN 8582	Incl w/NSN 8582	N/A
6640015757759		BZ-4 HEPA Filter (36"x48"x6")	BZ-HF4		1	EA	Incl w/NSN 8582	Incl w/NSN 8582	N/A
6640015772078		Stand, Manually Adjustable	BZ-A8BX		1	EA	Incl w/NSN 8582	Incl w/NSN 8582	N/A
6640015772669		BZ-4 Prefilter (48.25"x12.25"x1")	BZ-PF4		4	EA	Incl w/NSN 8582	Incl w/NSN 8582	N/A
7610045437003		Operation Manual	BZ-Manual		1	EA	Incl w/NSN 8582	Incl w/NSN 8582	N/A
NSN Not Required		Casters, Stand	BZ-A11		4	EA	Incl w/NSN 8582	Incl w/NSN 8582	N/A

**HOOD LAMINAR FLOW LAB  
ACCESSORIES & CONSUMABLES - RESUPPLY**

NSN	Notes	Nomenclature	MFR PN GermFree Laboratorie s Inc.	DIST PN	QTY	UI	Unit Price	Total Price	Shelf Life Mon.
6640015757759		BZ-4 HEPA Filter (36"x48"x6")	BZ-HF4		1	EA	\$455.00	\$455.00	N/A
6640015772669		BZ-4 Prefilter (48.25"x12.25"x1")	BZ-PF4		1	EA	\$125.00	\$125.00	N/A

**HOOD LAMINAR FLOW LAB  
WEB LINKS**

**[Click Below for Instruction Manual](http://97.74.181.168/downloads/manuals/BZ._Seriespdf.pdf)**  
**<[http://97.74.181.168/downloads/manuals/BZ.\\_Seriespdf.pdf](http://97.74.181.168/downloads/manuals/BZ._Seriespdf.pdf)>**

## REFRIGERATOR-FREEZER

TOC

**MANUFACTURER: NORTHLAND-MARVEL**

**Manufacturer Model# 6CRF0610**

**Manufacturer's PHONE: (616) 754-5601 FAX: (616) 754-0970**

**Manufacturer's Website: [www.marvelscientific.com](http://www.marvelscientific.com)**

**Item Web Link:**

**<http://www.marvelscientific.com/prod-combo.php>**

**AAC 'L' ITEM NSN 4110-01-425-8009**

**UAs: 263C, 273A, 623A, 905C, 912B, 913B, 914B, N306, N413, N503, N513, N703, N713, N725, P308, P309, P310, P506, P706**

**USAMMA POC: Customer Relations Management 301-619-1288/4301 - DSN 343-1288/4301**

**UNIT PRICE: \$2,219.18**

**System Description:** Six cubic foot under counter refrigerator/freezer fits easily under lab benches and casework.

**Features:** Forced air ventilation eliminates need for clearance space in built-in applications. Manual defrost, 5.3 cubic foot refrigerator capacity and 0.7 cubic foot freezer capacity. Dual voltage capable, carrying handles and olive drab paint scheme.



**REFRIGERATOR-FREEZER MECHANICAL FIELD WARD  
PART # 6CRF0610  
NSN 4110014258009**

**REFRIGERATOR-FREEZER  
PHYSICAL & ENVIRONMENTAL SPECIFICATIONS**

<b>Item Dimensions</b>	<b>Height (Inches)</b>	<b>35.5</b>
	<b>Width (Inches)</b>	<b>24.75</b>
	<b>Length (Inches)</b>	<b>25.5</b>
	<b>Cube (Cubic Feet)</b>	<b>12.9658</b>
	<b>Weight</b>	<b>120 lbs</b>
<b>Ship Dimensions</b>	<b>Height (Inches)</b>	<b>35.5</b>
	<b>Width (Inches)</b>	<b>24.75</b>
	<b>Length (Inches)</b>	<b>25.5</b>
	<b>Cube (Cubic Feet)</b>	<b>12.9658</b>
	<b>Weight</b>	<b>120 lbs</b>
<b>Storage Temp</b>		
<b>Operating Temp</b>		
<b>Storage Humidity</b>		
<b>Operating Humidity</b>		
<b>Amps</b>		<b>2.66</b>
<b>Watts</b>		<b>320 Normal - 1200 Surge</b>
<b>Power Requirements</b>	<b>120/220 VAC, 50/60 Hz</b>	
<b>Notes</b>	<b>Uses R-134A Refrigerant</b>	

**REFRIGERATOR-FREEZER  
ACCESSORIES & CONSUMABLES - STARTUP**

NSN	Notes	Nomenclature	MFR PN Northland- Marvel	DIST PN	QTY	UI	Unit Price	Total Price	Shelf Life Mon.
4110014258009		Refrigerator-Freezer	6CRF0610		1	EA	\$2,219.18	\$2,219.18	N/A
7610015437003		Literature, Operator	41011992		1	EA	Incl w/NSN 8009	Incl w/NSN 8009	N/A

**REFRIGERATOR-FREEZER  
ACCESSORIES & CONSUMABLES**

**NO RESUPPLY ITEMS**

## SINK UNIT SURGICAL SCRUB

[TOC](#)

**MANUFACTURER: ASEPTICO INC.**

**Manufacturer Model# ADU-42CF**

**Manufacturer's PHONE: (425) 487-3157 FAX: (360) 668-8722**

**Manufacturer's Website: <http://www.aseptico.com>**

**Item Web Link:**

**<http://www.aseptico.com/ADU-42CF.html>**

**AAC 'L' ITEM NSN 6530-01-572-6775**

**UAs: 211B, 213B, 218B, 267B, 274A, 504B, 905C, 912B, N312, N316, N413, N436, N503, N504, N513, N523, N713, N725, P301, P307, P308, P309, P310, P311, P314, P506**

**USAMMA POC: Customer Relations Management 301-619-1288/4301 - DSN 343-1288/4301  
UNIT PRICE: \$4,786.88**

**System Description:** This portable scrub sink provides a practical water source for mobile medical and dental operations. Both warm and cold water is available for washing and instrument cleansing as part of the process of disinfecting medical devices. Hands free operation makes surgical scrubs simple. The rigid stainless steel basin and instrument basket are durable and easy to disinfect. And, the accessory tray allows sink-side consolidation of all instruments. Light weight and compact size make the ADU-42CF easy to set-up and transport.

**Features:**

- Sink flow rates: High - 0.95 gpm (3.6 Lpm), Low - 0.64 gpm. (2.4 Lpm)
- Max. pump discharge head: 8.3 ft (2.5 meters)
- Pump self-priming limit: 48 inches (1.2 meters)
- Tank capacity: 1.45 gal (5.5 L)
- Tank Fill Time: High Flow - 1.1 minutes, Low Flow - 1.7 minutes
- Heating capacity: 1400 Watts
- Tank Heating Time: High Flow - 6.0 minutes, Low Flow - 6.5 minutes
- Peak water temperature: 106 °F, ±7 °F (41 °C ±4 °C)
- Cold water foot pump flow rate: 0.66 gpm (2.5 Lpm)
- Sink Basin Size: 21.3 x 15.4 x 7.4 inches (54.1 x 39.1 x 18.8 cm)
- Instrument Basket Size: 20.75 x 12.75 x 6 inches (52.7 x 32.4 x 15.2 cm)

**NOTE: Avoid freezing. To protect internal components, purge the sink's tank and lines prior to unattended periods in subfreezing conditions.**



**SINK UNIT SURGICAL SCRUB  
PART # ADU-42CF  
NSN 6530015726775**

**Items not supplied, Jerry cans or alternate water supply.**



**CASE TRANSPORT MIL-SPEC OD GREEN  
PART # 410191-08**



**CASE TRANSPORT MIL-SPEC  
PART # 410191-08**

**Case Shown Open With Contents**



**POUCH STORAGE LEG ASSY ADU-42CF CMPL  
PART # 410192-08**



**POUCH STORAGE ADU-40CF HOSES CMPL  
PART # 410131-08**



**POUCH STORAGE ADU-40CF TRAY CORDS  
PART # 410130-08**



**BASKET ASSY INSTRUMENT STAINLESS  
PART # 330550**



**SOAP DISPENSER ASSEMBLY  
PART # AA-42SD**



**PAPER TOWEL ASSY  
PART # 330532**

**SINK UNIT SURGICAL SCRUB  
PHYSICAL & ENVIRONMENTAL SPECIFICATIONS**

<b>Item Dimensions</b>	<b>Height (Inches)</b>	<b>0</b>
	<b>Width (Inches)</b>	<b>0</b>
	<b>Length (Inches)</b>	<b>0</b>
	<b>Cube (Cubic Feet)</b>	<b>0</b>
	<b>Weight</b>	<b>50 lbs</b>
<b>Ship Dimensions</b>	<b>Height (Inches)</b>	<b>12.5</b>
	<b>Width (Inches)</b>	<b>22.4</b>
	<b>Length (Inches)</b>	<b>25.3</b>
	<b>Cube (Cubic Feet)</b>	<b>4.0995</b>
	<b>Weight</b>	<b>104 lbs</b>
<b>Storage Temp</b>		<b>34 to 140 °F</b>
<b>Operating Temp</b>		<b>34 to 113 °F</b>
<b>Storage Humidity</b>		<b>5 to 95% Non-condensing</b>
<b>Operating Humidity</b>		<b>5 to 95% Non-condensing</b>
<b>Amps</b>		<b>12</b>
<b>Watts</b>		<b>1380</b>
<b>Power Requirements</b>	<b>115/230 VAC, 50/60Hz</b>	
<b>Notes</b>		

**SINK UNIT SURGICAL SCRUB  
ACCESSORIES & CONSUMABLES - STARTUP**

NSN	Notes	Nomenclature	MFR PN Aseptico Inc.	DIST PN	QTY	UI	Unit Price	Total Price	Shelf Life Mon.
6530015726775		Sink Unit Surgical Scrub	ADU-42CF		1	EA	\$4,786.88	\$4,786.88	N/A
	Awaiting NSN	Paper Towel Assy	330532		1	EA	Incl w/NSN 6775	Incl w/NSN 6775	N/A
	Awaiting NSN	Pouch Storage ADU-40CF Hoses Cmpl	410131-08		1	EA	Incl w/NSN 6775	Incl w/NSN 6775	N/A
	Awaiting NSN	Pouch Storage ADU-40CF Tray/ Cords	410130-08		1	EA	Incl w/NSN 6775	Incl w/NSN 6775	N/A
	Awaiting NSN	Pouch Storage Leg Assy ADU-42CF Cmpl	410192-08		1	EA	Incl w/NSN 6775	Incl w/NSN 6775	N/A
	Awaiting NSN	Pump & Tank Assembly	330522		1	EA	Incl w/NSN 6775	Incl w/NSN 6775	N/A
6530015883968		Regulator Assembly Water Pressure	330285		1	EA	Incl w/NSN 6775	Incl w/NSN 6775	N/A
6530015883979		Footswitch Switch Assembly Surgical Sink	AE-7PM		1	EA	Incl w/NSN 6775	Incl w/NSN 6775	N/A
6530015884124		Intake Hose Assembly	330252		1	EA	Incl w/NSN 6775	Incl w/NSN 6775	N/A
6530015884132		Waste Hose Assembly	330284		1	EA	Incl w/NSN 6775	Incl w/NSN 6775	N/A
6530015884170		Foot Pump Assembly	330458		1	EA	Incl w/NSN 6775	Incl w/NSN 6775	N/A
6530015885078		Soap Dispenser Assembly	AA-42SD		1	EA	Incl w/NSN 6775	Incl w/NSN 6775	N/A
6530015885138		Basket Assembly Surgical Sink CRS	330550		1	EA	Incl w/NSN 6775	Incl w/NSN 6775	N/A
6530015885163		Case Transport Surgical Sink	410191-08		1	EA	Incl w/NSN 6775	Incl w/NSN 6775	N/A
6530015888456		Leg Assembly	330521		1	EA	Incl w/NSN 6775	Incl w/NSN 6775	N/A
6530015888537		Tray Instrument	330246		1	EA	Incl w/NSN 6775	Incl w/NSN 6775	N/A
6530015888935		Spout Assembly	120354		1	EA	Incl w/NSN 6775	Incl w/NSN 6775	N/A
7610015437003		Literature Operation	420752		1	EA	Incl w/NSN 6775	Incl w/NSN 6775	N/A

**SINK UNIT SURGICAL SCRUB  
ACCESSORIES & CONSUMABLES - RESUPPLY**

NSN	Notes	Nomenclature	MFR PN Aseptico Inc.	DIST PN	QTY	UI	Unit Price	Total Price	Shelf Life Mon.
	Awaiting NSN	Paper Towel Assy	330532		1	EA	\$62.80	\$62.80	N/A
	Awaiting NSN	Pouch Storage ADU-40CF Hoses Cmpl	410131-08		1	EA	\$36.52	\$36.52	N/A
	Awaiting NSN	Pouch Storage ADU-40CF Tray/ Cords	410130-08		1	EA	\$35.39	\$35.39	N/A
	Awaiting NSN	Pouch Storage Leg Assy ADU-42CF Cmpl	410192-08		1	EA	\$35.14	\$35.14	N/A
	Awaiting NSN	Pump & Tank Assembly	330522		1	EA	\$1,133.94	\$1,133.94	N/A

## ACCESSORIES & CONSUMABLES - RESUPPLY

NSN	Notes	Nomenclature	MFR PN Aseptico Inc.	DIST PN	QTY	UI	Unit Price	Total Price	Shelf Life Mon.
6530015883968		Regulator Assembly Water Pressure	330285		1	EA	\$159.74	\$159.74	N/A
6530015883979		Footswitch Switch Assembly Surgical Sink	AE-7PM		1	EA	\$92.26	\$92.26	N/A
6530015884124		Intake Hose Assembly	330252		1	EA	\$80.05	\$80.05	N/A
6530015884132		Waste Hose Assembly	330284		1	EA	\$20.47	\$20.47	N/A
6530015884170		Foot Pump Assembly	330458		1	EA	\$146.64	\$146.64	N/A
6530015885078		Soap Dispenser Assembly	AA-42SD		1	EA	\$114.25	\$114.25	N/A
6530015885138		Basket Assembly Surgical Sink CRS	330550		1	EA	\$340.03	\$340.03	N/A
6530015885163		Case Transport Surgical Sink	410191-08		1	EA	\$828.96	\$828.96	N/A
6530015888456		Leg Assembly	330521		1	EA	\$187.50	\$187.50	N/A
6530015888537		Tray Instrument	330246		1	EA	\$73.71	\$73.71	N/A
6530015888935		Spout Assembly	120354		1	EA	\$198.21	\$198.21	N/A
7610015437003		Literature Operation	420752		1	EA			N/A



