

SCREEN X-RAY PROTECTIVE MOBILE

MANUFACTURER: FLUKE BIOMEDICAL

Manufacturer Model# 56-603

Manufacturer's PHONE: (440) 498-2565 FAX: (440) 498-0000

Manufacturer's Website: <<http://www.flukebiomedical.com/biomedical/usen/home/>>

Item Web Link:

<<http://www.flukebiomedical.com/biomedical/usen/radiation-safety/radiation-protection/56-663-series.htm?PID=55000>>

AAC 'L' ITEM NSN 6525-01-166-9033

UAs: N305

USAMMA POC: Customer Relationship Management 301-619-1288/4301 - DSN 343-1288/4301

UNIT PRICE: \$833.09

System Description: The regular and the jumbo mobile barriers can safely shield one or two persons. The updated stainless steel hardware provides the “look” to your discriminating healthcare designers and their health care staff. These models are available in 0.5 mm lead equivalence. This shield is ideal for Operating Room Special Procedure Suites and Cardiac Catheterization Labs. The stainless steel hardware mobile barriers are the perfect fit in any imaging room due to their 30 inch width and easy maneuverability.



**SCREEN X-RAY PROTECTIVE MOBILE 6 FT. X 2 FT.
PART # 56-603
NSN 6525011669033**

**SCREEN X-RAY PROTECTIVE MOBILE
PHYSICAL & ENVIRONMENTAL SPECIFICATIONS**

Item Dimensions	Height (Inches)	72
	Width (Inches)	30
	Length (Inches)	2
	Cube (Cubic Feet)	2.5
	Weight	110 lbs
Ship Dimensions	Height (Inches)	
	Width (Inches)	
	Length (Inches)	
	Cube (Cubic Feet)	
	Weight	
Storage Temp		
Operating Temp		
Storage Humidity		
Operating Humidity		
Amps		
Watts		
Power Requirements	N/A	
Notes	Window Height: 24 inches Window Width: 30 inches Lead equivalent window panel: 24" by 30", 0.5 mm Lead equivalent opaque panel: 48" by 30", Min 0.03 in (0.8mm)	

**SCREEN X-RAY PROTECTIVE MOBILE
ACCESSORIES & CONSUMABLES - STARTUP**

NSN	Notes	Nomenclature	MFR PN Fluke BioMedica 1	DIST PN	QTY	UI	Unit Price	Total Price	Shelf Life Mon.
6525011669033		Screen X-Ray Protective Mobile	56-603		1	EA	\$833.09	\$833.09	N/A

SCREEN X-RAY PROTECTIVE MOBILE
WEB LINKS

[Click Below for Literature, Operator Manual](#)

http://assets.fluke.com/manuals/56_614__umeng0000.pdf