

SINK UNIT SCRUB FIELD

MANUFACTURER: FISHER HAMILTON LLC

Manufacturer Model# 950S936

Manufacturer's PHONE: (800) 826-1166 FAX: (412) 490-8304

Manufacturer's Website: <<http://www.fishersci.com>>

AAC 'J' ITEM NSN 6530-01-308-7740

UAs: N305, N306, N432, N703, N704, P302, P433, P706

USAMMA POC: Customer Relationship Management 301-619-1288/4301 - DSN 343-1288/4301

UNIT PRICE: \$4,753.28


System Description: This sink supplies hot and cold water and is mounted on four detachable casters. It is a self contained sink with a tankless electric water heater and pump.

Features:

- Removable faucet with a gooseneck spout
- Paddle shaped handles to control the water flow
- Foot activated, externally mounted soap dispenser
- Two recessed carrying handles on each side
- Two doors
- A 10-foot water supply hose with auto shut-off and quick disconnect couplings
- Water temperature adjustable from 100°F to 145°F
- Corrosion-resistant steel (CRS) top and sink



**SINK UNIT SCRUB FIELD HOSPITAL CRS 30 INCHES
PART # 950S936
NSN 6530013087740**



- Self-contained sink base cabinet is 30" wide X 22-5/8" deep X 36-1/2" high
- Includes 24" X 14" X 12" stainless steel X-Ray sink, self-contained pump and heater, foot operated soap dispenser.
- 4" casters

**SINK UNIT SCRUB FIELD
PART # 950S936**

**SINK UNIT SCRUB FIELD
PHYSICAL & ENVIRONMENTAL SPECIFICATIONS**

Item Dimensions	Height (Inches)	36.5
	Width (Inches)	30
	Length (Inches)	22.625
	Cube (Cubic Feet)	14.337
	Weight	335 lbs
Ship Dimensions	Height (Inches)	
	Width (Inches)	
	Length (Inches)	
	Cube (Cubic Feet)	
	Weight	
Storage Temp		
Operating Temp		
Storage Humidity		
Operating Humidity		
Amps		38
Watts		9000
Power Requirements	120/230 VAC, 50/60Hz	
Notes		

**SINK UNIT SCRUB FIELD
ACCESSORIES & CONSUMABLES - STARTUP**

NSN	Notes	Nomenclature	MFR PN Fisher Hamilton LLC	DIST PN	QTY	UI	Unit Price	Total Price	Shelf Life Mon.
6530013087740		Sink Unit Scrub Field	950S936		1	EA	\$4,753.28	\$4,753.28	N/A

**SINK UNIT SCRUB FIELD
ADDITIONAL NOTE:**

To order Operator and Maintenance Literature for this sink, please refer to the USAMMA Medical Equipment Literature Web page (http://www.usamma.army.mil/med equip_lit_cdroms.cfm); click on "Order the CD-ROM here," fill in the requested information, and select "Operator & Maintenance Literature for Medical Equipment - Discs 3 & 4"