

STERILIZER SURGICAL INSTRUMENT STEAM

DISTRIBUTOR: GSS INC.

MANUFACTURER: TUTTNAUER USA CO.

Distributor Model# TUT2540MPACKAGEEA

Manufacturer Model# 2540M

Distributor's PHONE: (703) 734-1805 FAX: (703) 734-1803

Manufacturer's PHONE: (631) 737-4850 FAX: (631) 737-0720

Distributor's Website: <<http://www.govsci.com>>

Manufacturer's Website: <<http://www.tuttnauer.com>>

Item Web Link:

<<http://www.tuttnauer.com/products/laboratory-autoclaves/benchtop-laboratory-sterilizers/ml-semi-automatic-benchtop-sterilizers>>

AAC 'L' ITEM NSN 6530-01-531-4807

ECAT Item Number: TUT2540MPACKAGEEA

UAs: 267C, 274A, 512A, 905C, 912B

USAMMA POC: Customer Relationship Management 301-619-1288/4301 - DSN 343-1288/4301

UNIT PRICE: \$3,611.54

System Description: Tuttnauer manual autoclaves set the standard for steam sterilization. The Tuttnauer manual autoclave is designed to compliment any healthcare and laboratory facility. Our basic manual autoclave will satisfy all of your sterilization needs. Time tested reliability and durability are built into every Tuttnauer product. The Tuttnauer M Series is also available in a quick cycle model. This sterilizer offers all the proven, dependable features of our manual autoclave, with the added benefit of reducing the total sterilization cycle time by 50%. Speed and reliability are accomplished without sacrificing load size. This is the perfect choice for quick instrument turn-around or a high-volume office.

Features:

- Sterilizer is a table top steam autoclave
- Heat insulation and quiet operation for pleasant environment
- Durable 316L stainless steel chamber and door with electro-polish finish
- Double locking safety device prevents door from opening while chamber is pressurized
- Steam pressure activated lock prevents door from opening if chamber is pressurized
- Door failure protection prevents door from opening at high pressure
- Automatic safety shutoff prevents overheating
- Insufficient water detector with automatic shutoff
- Fast cycle time: 14 minutes from a hot start; 22 minutes from a cold start
- Drain valve located on the front, allowing for quick and easy draining of water reservoir
- Sterilization cycle and separate drying cycle feature automatic shut off




**STERILIZER SURGICAL INSTRUMENT STEAM
PART # 2540M
NSN 6530015314807**



**HARDENED CASE FOR THE TUTTNAUER STERILIZER
PART # HDD041805200
NSN 6530015392908**



STERILIZER SHOWN IN HARDENED CASE



**INSTRUMENT POUCH RACK
NSN 6530015408014**



DISTILLER UNIT WATER
PART # 9000
NSN 6530015329092
Found in UA Support

**STERILIZER SURGICAL INSTRUMENT STEAM
PHYSICAL & ENVIRONMENTAL SPECIFICATIONS**

Item Dimensions	Height (Inches)	14
	Width (Inches)	20
	Length (Inches)	20
	Cube (Cubic Feet)	3.2407
	Weight	33 lbs
Ship Dimensions	Height (Inches)	21
	Width (Inches)	25
	Length (Inches)	28
	Cube (Cubic Feet)	8.5069
	Weight	80 lbs
Storage Temp		59 to 95 °F
Operating Temp		59 to 95 °F
Storage Humidity		0 to 85%
Operating Humidity		0 to 85%
Amps		12/6
Watts		1400
Power Requirements	120/230 VAC 50/60 Hz 12 amps at 120VAC	
Notes	Chamber Volume: 23 liter Number of Trays: 4 Tray Dimensions: 6.6 inches wide X 0.8 inches height X 16.3 inches deep	

**STERILIZER SURGICAL INSTRUMENT STEAM
ACCESSORIES & CONSUMABLES - STARTUP**

NSN	Notes	Nomenclature	MFR PN Tuttnauer USA Co.	DIST PN GSS Inc.	QTY	UI	Unit Price	Total Price	Shelf Life Mon.
6530015314807		Sterilizer Surgical Instrument Steam		2540M	1	EA	\$3,611.54	\$3,611.54	N/A
6530015392908		Case, Hardened		HDD04180 5200	1	EA	Incl w/NSN 4807	Incl w/NSN 4807	N/A
6530015408014		Pouch Rack		AR910	2	EA	Incl w/NSN 4807	Incl w/NSN 4807	N/A
6530015408023		Tray		CT620010	4	EA	Incl w/NSN 4807	Incl w/NSN 4807	N/A
7610015408195		Manual, Operation/Service		EMHC01.	1	EA	Incl w/NSN 4807	Incl w/NSN 4807	N/A

STERILIZER SURGICAL INSTRUMENT STEAM
WEB LINKS

[Click Below for Product Brochure](#)

<http://www.tuttnauer.com/sites/default/files/Tabletop%20Line%20-%20Automatic%20and%20Semi-Automatic%20Autoclaves%20-%20English_0.pdf>

Click Below for Equipment Support Plans. Note: Only available to personnel who reside in a Military (.mil) domain. Please select the plans you wish to receive and enter your email address. Support Plans include 10/20 Standards, maintenance intervals, Maintenance Allowance Chart, and more.

<http://www.usamma.army.mil/assets/apps/email_doc/email_doc_index.cfm>