

## ANALYZER SET CHEMISTRY URINE

**MANUFACTURER: SIEMENS**

**Manufacturer Model# Clinitek Advantus**

**Manufacturer's PHONE: (914) 631-8000 FAX: (914) 524-2132**

**Manufacturer's Website: <<http://www.usa.siemens.com>>**

**Item Web Link:**

**[http://www.medical.siemens.com/webapp/wcs/stores/servlet/ProductDisplay~q\\_catalogId~e - 101~a\\_catTree~e\\_100004~a\\_langId~e\\_-101~a\\_productId~e\\_175341~a\\_storeId~e\\_10001.htm](http://www.medical.siemens.com/webapp/wcs/stores/servlet/ProductDisplay~q_catalogId~e-101~a_catTree~e_100004~a_langId~e_-101~a_productId~e_175341~a_storeId~e_10001.htm)**

**AAC 'J' ITEM NSN 6630-01-562-1648**

**UAs: N503, N703**

**USAMMA POC: Customer Relationship Management 301-619-1288/4301 - DSN 343-1288/4301**

**UNIT PRICE: \$11,020.65**

**System Description:** The Clinitek Advantus is the next generation of the Clinitek 500 Analyzer. This self-pacing compact analyzer is designed to increase productivity and enhance the efficiency in a semiautomated urinalysis laboratory. Innovative system offers automatic reagent strip detection, multiple printed report formats, rapid recall of patient results and enhanced data management, including consolidation of microscopy results. Increase accuracy and confidence delivered through automatic calibration with each sample and barcode entry of patient sample information.

**Features:**

- New, easy-to-use color touch screen (see picture of Input Screen below)
- Ability to consolidate microscopy results onto one patient record
- Increased control with automatic QC prompting and QC lockout capability to ensure compliance
- Testing flexibility with simple, one-touch switching between different strip configurations
- A total of three sieve functions – confirmatory and microscopic to flag abnormal results for quick review
- Convenient data retrieval for up to 500 patients by patient id, and 200 control results
- Streamlined operation with immediate start-up
- Improved efficiency with nonpacing processing
- Improve accuracy with automatic calibration for every strip
- Continued productivity boost with throughput of 7 seconds per sample
- Minimize errors with simple, barcode identification of sample, color and clarity values
- Flexibility in reporting options, on-screen display, printout and/or uploading to laboratory information system




**ANALYZER SET, CHEMISTRY URINE  
PART # Clinitek Advantus  
NSN 6630015621648**



**USER INTERFACE, TESTING AND PRINTING AREAS**



**CASE, RUGGEDIZED  
PART # DSI-30308-1500**




**PAPER ROLLS  
PART # 5773**



**LUBRIPLATE LUBRICANT  
PART # 50336008**



**READER BAR CODE HANDHELD  
PART # 1421**



INPUT SCREEN

**ANALYZER SET CHEMISTRY URINE  
PHYSICAL & ENVIRONMENTAL SPECIFICATIONS**

<b>Item Dimensions</b>	<b>Height (Inches)</b>	<b>12.75</b>
	<b>Width (Inches)</b>	<b>15.75</b>
	<b>Length (Inches)</b>	<b>13.75</b>
	<b>Cube (Cubic Feet)</b>	<b>1.5979</b>
	<b>Weight</b>	<b>16.30 lbs</b>
<b>Ship Dimensions</b>	<b>Height (Inches)</b>	
	<b>Width (Inches)</b>	
	<b>Length (Inches)</b>	
	<b>Cube (Cubic Feet)</b>	
	<b>Weight</b>	
<b>Storage Temp</b>		
<b>Operating Temp</b>		<b>64 to 86 °F</b>
<b>Storage Humidity</b>		
<b>Operating Humidity</b>		<b>20% - 80%</b>
<b>Amps</b>		<b>72 volt amps</b>
<b>Watts</b>		<b>80</b>
<b>Power Requirements</b>	<b>100/240 VAC, 50/60 HZ</b>	
<b>Notes</b>	<b>Maintainer Test Kit:</b> <b>NSN 6630-01-611-8506, Repair Kit Chemistry Analyzer</b> <b>Located in the following UA's:</b> <b>N523, N725, T500, T503</b>	

**ANALYZER SET CHEMISTRY URINE  
ACCESSORIES & CONSUMABLES - STARTUP**

NSN	Notes	Nomenclature	MFR PN Siemens	DIST PN	QTY	UI	Unit Price	Total Price	Shelf Life Mon.
6630015621648		Analyzer Set Chemistry Urine	Clinitek Advantus		1	EA	\$11,020.65	\$11,020.65	N/A
	Awaiting NSN	Card, Color/Clarity	111028		1	EA	Incl w/NSN 1648	Incl w/NSN 1648	N/A
	Awaiting NSN	Case, Ruggedized, Custom	DSI- 30308- 1500		1	EA	Incl w/NSN 1648	Incl w/NSN 1648	N/A
	Awaiting NSN	Cord Power AC (European)	00171415A		1	EA	Incl w/NSN 1648	Incl w/NSN 1648	N/A
	Awaiting NSN	Cord Power AC (US)	40330028		1	EA	Incl w/NSN 1648	Incl w/NSN 1648	N/A
	Awaiting NSN	Fixed Platform/Holdown/Waste Bin	133374		1	EA	Incl w/NSN 1648	Incl w/NSN 1648	N/A
	Awaiting NSN	Lubriplate Lubricant	50336008		1	EA	Incl w/NSN 1648	Incl w/NSN 1648	N/A
	Awaiting NSN	Manual Operator CD	1133594		1	EA	Incl w/NSN 1648	Incl w/NSN 1648	N/A
	Awaiting NSN	Moving Table	50577530		1	EA	Incl w/NSN 1648	Incl w/NSN 1648	N/A
	Awaiting NSN	Paper Roll	5773		1	5/BX	Incl w/NSN 1648	Incl w/NSN 1648	N/A
	Awaiting NSN	Push Bar	50540224		1	EA	Incl w/NSN 1648	Incl w/NSN 1648	N/A
	Awaiting NSN	Reader Bar Code handheld	1421		1	EA	Incl w/NSN 1648	Incl w/NSN 1648	N/A
6550010934010		Chex-Sticks +25,-25	1364		1	25/BT	Incl w/NSN 1648	Incl w/NSN 1648	18
6550011225540		Test Strips & Color Chart Bilirubin Bld Glucose	2161		1	100/BT	Incl w/NSN 1648	Incl w/NSN 1648	24
6630015619938		Analyzer Chemistry Urine	1420		1	EA	Incl w/NSN 1648	Incl w/NSN 1648	N/A
6630015623969		Waste Bin Liners	6472		1	5/PG	Incl w/NSN 1648	Incl w/NSN 1648	N/A

**ANALYZER SET CHEMISTRY URINE  
WEB LINKS**

**Click Below for Equipment Support Plans. Note: Only available to personnel who reside in a Military (.mil) domain. Please select the plans you wish to receive and enter your email address. Support Plans include 10/20 Standards, maintenance intervals, Maintenance Allowance Chart, and more.**

**<[http://www.usamma.army.mil/assets/apps/email\\_doc/email\\_doc\\_index.cfm](http://www.usamma.army.mil/assets/apps/email_doc/email_doc_index.cfm)>**