

FIXATION SET EXTERNAL ORTHOPEDIC HOFFMAN II TOC

STRYKER/HOWMEDICA

MODEL #: 4920-9-962

Manufacturer's PHONE: (201) 831-5726 FAX: (201) 831-6553

Website: <http://www.trauma.stryker.com/physicians/usa/hoffmann2.php>


AAC "L" ITEM NSN 6515-01-551-1198

UA's: 1366, 1370, 1374, 267B, F302, N302, N583, N783

USAMMA POC: Customer Relations Management, 301-619-4301/4316, DSN 343-4301/4316

UNIT PRICE: \$22,397.00

System Description: Device providing rigid immobilization through external skeletal fixation utilizing rods attached to pins placed in or through the bone. Used for both open fractures and septic nonunion. Pictured below are examples of this device in use.




FIXATION SET EXTERNAL ORTHOPEDIC HOFFMAN II PHYSICAL & ENVIRONMENTAL SPECIFICATIONS

Dimensions	Height	21.6 inches
	Width	10.2 inches
	Depth	4.3 inches
Cube		0.54 cu ft
Weight		20 lb

X Level 3/4/5 External Fixation System


- ◆ Components and instruments interchangeable with the Forward Deployed Field Pack Fixator
- ◆ Color-coded components for improved intra-operative and post-operative inventory management
- ◆ Integrated snap-fit mechanism allows easy assembly and temporary fixation during reduction maneuvers
- ◆ Clamps made of aluminum alloys and stainless steel, significantly reducing overall frame weight, without compromising stability

X Upper Insert



Component Listing			
Description	Qty	Catalog Number	
A Rod-to-Rod Coupling 8/8mm	16	4920-1-010	
B Removable Thumbwheel	2	4920-9-020	
C Pin-to-Rod Coupling 8mm/4+5mm	6	4920-1-020	
D 4mm x 300mm, 40mm Central Thread	2	5030-4-300	
E 5mm x 180mm, 50mm THD SD/ST Apex® Half Pins	26	5018-6-180	
F 5mm x 150mm, 40mm THD SD/ST Apex® Half Pins	12	5018-5-150	
G 4mm x 150mm, 40mm THD SD/ST Apex® Half Pins	12	5023-5-150	

X Lower Insert




Component Listing			
Description	Qty	Catalog Number	
H Multi-Pin Clamp Assembly	16	4920-2-020	
I 30° Angled Post, Stainless Steel	16	4920-2-140	
J Guide Sleeve with Trocar for 5mm/150mm pins	2	5054-5-150	
K Guide Sleeve with Trocar for 4mm/150mm pins	2	5054-4-150	

FIXATION SET EXTERNAL ORTHOPEDIC HOFFMAN II
MODEL # 4920-9-962
NSN 6515-01-551-1198

X Level 3/4/5 External Fixation System - continued

X Base



Component Listing		
Description	Qty	Catalog Number
L Curved Aluminum Rod Semi-Circular Large 214mm Width	2	5029-7-032
M T-Wrench	1	5054-3-051
N Stabilization/Reduction Wrench	2	4920-9-010
O Carbon Connecting Rod 8mm x 400mm	6	5029-8-840
P Carbon Connecting Rod 8mm x 300mm	6	5029-8-830
Q Drill Brace	1	5054-2-002
R 4mm Chuck	1	5054-2-004
S 5mm Chuck	1	5054-2-051

FIXATION SET EXTERNAL ORTHOPEDIC HOFFMAN II (continued)
 MODEL # 4920-9-962
 NSN 6515-01-551-1198

FIXATION SET EXTERNAL ORTHOPEDIC HOFFMAN II (START-UP)

NSN	NOTES	Nomenclature	PN	QTY	UI	Unit Price	Total Price	DIST/MFG	Shelf Life, months	Refrig Y/N
6515015511198		Fixation set external	4920-9-962	1	EA	\$22,397.00	\$22,397.00	Howmedica	indef	N
UPPER INSERT										
6530015224432		Case, reusable, steam sterilizable custom	4920-9-9	1	EA	Included w/ 0077 NSN	Included w/ 0077 NSN	Howmedica	indef	N
6515014631498		Coupling fixation ortho 8/8mm	4920-1-010	16	EA	Included w/ 0077 NSN	Included w/ 0077 NSN	Howmedica	indef	N
6515015224925		Thumbwheel, removable	4920-9-020	2	EA	Included w/ 0077 NSN	Included w/ 0077 NSN	Howmedica	indef	N
6515014631473		Coupling, pin to rod 8mm/4-5mm	4920-1-020	6	EA	Included w/ 0077 NSN	Included w/ 0077 NSN	Howmedica	indef	N
6515014631494		Pin ortho 4mm dia x 300mm Lg 40mm central thread	5030-4-300	2	EA	Included w/ 0077 NSN	Included w/ 0077 NSN	Howmedica	indef	N
6515015131261		Pin bone 5mm dia	5018-6-180	1	26/PG	Included w/ 0077 NSN	Included w/ 0077 NSN	Howmedica	indef	N
6515015039984		Pin bone CRS 5mm dia 150mm Lg half pin fix	5018-5-150	12	EA	Included w/ 0077 NSN	Included w/ 0077 NSN	Howmedica	indef	N
6515015061850		Pin bone 150mm Lg 4mm thread	5023-5-150	12	EA	Included w/ 0077 NSN	Included w/ 0077 NSN	Howmedica	indef	N
LOWER INSERT										
6515014631475		Pin clamp assembly w/ no post ortho	4920-2-020	16	EA	Included w/ 0077 NSN	Included w/ 0077 NSN	Howmedica	indef	N
6515014631497		Post angle fixation ortho 30 deg SS	4920-2-140	16	EA	Included w/ 0077 NSN	Included w/ 0077 NSN	Howmedica	indef	N
6515013532367		Fixator dbl sleeve Trocar ortho	5054-5-150	2	EA	Included w/ 0077 NSN	Included w/ 0077 NSN	Howmedica	indef	N
6515014469703		Fixator dbl sleeve Trocar 150mm	5054-4-150	2	EA	Included w/ 0077 NSN	Included w/ 0077 NSN	Howmedica	indef	N
BASE										
6515014631467		Rod semi circular alum ext	5029-7-032	2	EA	Included w/ 0077 NSN	Included w/ 0077 NSN	Howmedica	indef	N
6515013532554		Wrench fixation, T-wrench s/ 35mm ext SS	5054-3-051	1	EA	Included w/ 0077 NSN	Included w/ 0077 NSN	Howmedica	indef	N
6515015232438		Wrench, stabilization/reduction	4920-9-010	2	EA	Included w/ 0077 NSN	Included w/ 0077 NSN	Howmedica	indef	N
6515015040434		Rod carbon connect ext 400mm Lg 8mm W	5029-8-840	6	EA	Included w/ 0077 NSN	Included w/ 0077 NSN	Howmedica	indef	N
6515015040433		Rod carbon connect ext 300mm Lg 8mm W	5029-8-830	6	EA	Included w/ 0077 NSN	Included w/ 0077 NSN	Howmedica	indef	N
6515012085966		Brace drill bit w/o chuck, reusable	5054-2-002	1	EA	Included w/ 0077 NSN	Included w/ 0077 NSN	Howmedica	indef	N
6515012969606		Chuck drill plain 4mm arbor hole thread	5054-2-004	1	EA	Included w/ 0077 NSN	Included w/ 0077 NSN	Howmedica	indef	N
6515012085965		Chuck AGM 51 5mm 75mm ext	5054-2-051	1	EA	Included w/ 0077 NSN	Included w/ 0077 NSN	Howmedica	indef	N