

ELECTROCARDIOGRAPH CP200 KIT W CASE

MANUFACTURER: WELCH ALLYN

Manufacturer Model# CP2A-1EM

Manufacturer's PHONE: (315) 685-4560 FAX: (315) 683-3361

Manufacturer's Website: <<http://www.welchallyn.com>>

Item Web Link:

<<http://www.welchallyn.com/products/en-us/x-11-ac-100-000000001026.htm>>

AAC 'L' ITEM NSN 6515-01-563-7185

UAs: 009A, N513, N713, P308, P309

USAMMA POC: Customer Relationship Management 301-619-1288/4301 - DSN 343-1288/4301

UNIT PRICE: \$4,410.24

System Description: A multichannel electrocardiograph that is flexible enough to work in any setting. The CP200 is a powerful, easy-to-use device featuring an intuitive design and color display that make taking ECGs (electrocardiograms) easier.

- Features:**
- An interpretive 12-lead resting multichannel ECG that provides accurate measurements
 - Full alphanumeric keyboard and color LCD screen to enter patient data quickly
 - High-resolution thermal printer generates standard full-sized reports
 - Help feature to access patient education tools as well as answer system operation questions
 - Frequency range of 0.3 to 150 Hz
 - Digital sampling rate of greater than 1,000 samples/second/channel
 - Prints up to 100 continuous ECGs per battery charge




**ELECTROCARDIOGRAPH CP200
PART # CP2A-1EM
NSN 6515015637185**



**CABLE, PATIENT DIGITAL ECG
PART # 401126
NSN 6515015609895**



**10-LEAD REPLACEABLE LEAD SET
PART # 401129
NSN 6515015610335**



**PAPER, THERMAL
PART # 94018-0000
NSN 6515015527128**



ALLIGATOR CLIP
PART # 58581-0000
NSN 6515015053043



TAB ELECTRODES
PART # 45008-0000
NSN 6515015053037

**ELECTROCARDIOGRAPH CP200 KIT W CASE
PHYSICAL & ENVIRONMENTAL SPECIFICATIONS**

Item Dimensions	Height (Inches)	6.2
	Width (Inches)	16.2
	Length (Inches)	15.6
	Cube (Cubic Feet)	0.9068
	Weight	11.6 lbs
Ship Dimensions	Height (Inches)	0
	Width (Inches)	0
	Length (Inches)	0
	Cube (Cubic Feet)	0
	Weight	
Storage Temp		-4 to 120 °F
Operating Temp		32 to 104 °F
Storage Humidity		15 – 95% Noncondensing (30 – 70% for printing)
Operating Humidity		15 – 95% Noncondensing (30 – 70% for printing)
Amps		2.2 Max
Watts		13.2
Power Requirements	100-240 VAC, 50/60Hz	
Notes	Battery: Rechargeable lead acid gel, 6 VDC, 5 Ah (ampere-hour), 12 hour recharging Paper: 200 sheets/PG, 8.25 x 11 inches, Z-fold thermal	

**ELECTROCARDIOGRAPH CP200 KIT W CASE
ACCESSORIES & CONSUMABLES - STARTUP**

NSN	Notes	Nomenclature	MFR PN Welch Allyn	DIST PN	QTY	UI	Unit Price	Total Price	Shelf Life Mon.
6515015637185		Electrocardiograph CP200 Kit w Case	CP2A- 1EM		1	EA	\$4,410.24	\$4,410.24	N/A
6515015053025		Electrocardiograph CP200	CP2A		1	EA	Incl w/NSN 7185	Incl w/NSN 7185	N/A
6515015609895	Mfr uses PN 400293 for Resupply	Cable, Patient Digital ECG, 10-lead, Replaceable, AHA	401126	400293	1	EA	Incl w/NSN 7185	Incl w/NSN 7185	N/A
6515015609958		Manual Operator, Training Guide & Physician's Guide, CD	401151		1	EA	Incl w/NSN 7185	Incl w/NSN 7185	N/A
6515015610335	Mfr uses PN 401129 for Resupply	10-Lead Replaceable Lead Set, AHA, Banana Clip, 4 Limb Leads, 6 Precordial (for Use w/400293)	401672	401129	1	SE	Incl w/NSN 7185	Incl w/NSN 7185	N/A
6515015611869		Case, Storm, Hard-Sided	IM2720		1	EA	Incl w/NSN 7185	Incl w/NSN 7185	N/A
6540015009058		AC Power cord	761076-0		1	EA	Incl w/NSN 7185	Incl w/NSN 7185	N/A
NSN Not Required		Case, Soft-Sided/Dust Cover	701586		1	EA	Incl w/NSN 7185	Incl w/NSN 7185	N/A
NSN Not Required	See Additional Note 2 Below	Electrode ECG Disposable 100s	45008- 0001		1	100/PG	Incl w/NSN 7185	Incl w/NSN 7185	N/A
NSN Not Required		European Power Cord	761076-2		1	EA	Incl w/NSN 7185	Incl w/NSN 7185	N/A
NSN Not Required		Literature, Maintenance, CD	704340		1	EA	Incl w/NSN 7185	Incl w/NSN 7185	N/A
NSN Not Required	See Additional Note 1 Below	Paper, Thermal, 200 sheets	704367		1	RM	Incl w/NSN 7185	Incl w/NSN 7185	60

**ELECTROCARDIOGRAPH CP200 KIT W CASE
ADDITIONAL NOTE:**

1. Paper must be stored at 77 °F or less and 45 to 65% Relative Humidity.
2. Start-up quantity or unit of issue is different from Resupply quantity/unit of issue.

**ELECTROCARDIOGRAPH CP200 KIT W CASE
WEB LINKS**

[Click Below for Service Manual](#)

<http://www.welchallyn.com/documents/Cardiopulmonary/Electrocardiographs/CP200/ServiceManual_2007-8-3_CP100&CP200.pdf>

[Click Below for Equipment Support Plans. Note: Only available to personnel who reside in a Military \(.mil\) domain. Please select the plans you wish to receive and enter your email address. Support Plans include 10/20 Standards, maintenance intervals, Maintenance Allowance Chart, and more.](#)

<http://www.usamma.army.mil/assets/apps/email_doc/email_doc_index.cfm>