

WATER TEST KIT BACTERIO (NEEDS START-UP)

MANUFACTURER: HACH COMPANY

Manufacturer Model# 2569700

Manufacturer's PHONE: (970) 669-3050 FAX: (970) 962-6739

Manufacturer's Website: <<http://www.hach.com>>

Item Web Link:

<<http://www.hach.com/mel-mf-total-coliform-laboratory/product?id=7640250897>>

AAC 'L' ITEM NSN 6665-01-467-3718

UAs: 218B

USAMMA POC: Customer Relationship Management 301-619-1288/4301 - DSN 343-1288/4301

UNIT PRICE: \$3,615.26

System Description: From a remote field location to a small scale microbiology lab within your existing lab, this Microbiological Environmental Lab (MEL) offers all the equipment and accessories to perform microbiological analysis.

Our MEL/MF (Membrane Filtration) Total Coliform Lab comes as a complete membrane filtration system. This lab is ideal for supporting disinfection adjustments or modifications, short-term and seasonal monitoring, emergency field operations, or as a microbiology lab within an existing facility.

A charged battery for the portable incubator will last 12 hours at ambient temperature. Because microbiological tests normally require at least 24 hours incubation, the portable incubator is best used for transport of samples to a place where the incubator can be plugged into an electrical outlet using the battery eliminator. The portable incubator can also be plugged into the cigarette lighter of an automobile.

- Features:**
- Step-by-step, illustrated procedures manual
 - Portable Incubator
 - Can be powered from automobile cigarette lighter receptacle.)
 - 200 Whirl-Pak Bags with dechlorinating agent for sampling
 - 200 funnel assemblies individually packed for sterility; components:
 - Membrane filter
 - Petri dish with pad
 - Funnel
 - Portable, illuminated magnifier with 2.5X and 5X magnification
 - Forceps and alcohol burner
 - Armored alcohol thermometer, 14 to 230 °F (-10 to 110 °C)
 - Five germicidal cloths for disinfecting test surfaces
 - Colony counter
 - MF Petri Dish Rack
 - Media must be ordered separately (ensures longest shelf-life)



**WATER TESTING KIT BACTERIOLOGICAL
PART # 25697-00
NSN 6665014673718**



**WATER TESTING KIT BACTERIOLOGICAL
PART # 2569700
NSN 6665014673718**



**INCUBATOR PORTABLE 12VDC 30° TO 50° C
PART # 2569900
NSN 6640015523532**



**BATTERY PACK
PART # 2580300
NSN 6130015530374**

**WATER TEST KIT BACTERIO (NEEDS START-UP)
PHYSICAL & ENVIRONMENTAL SPECIFICATIONS**

Item Dimensions	Height (Inches)	12
	Width (Inches)	12
	Length (Inches)	10
	Cube (Cubic Feet)	0.8333
	Weight	4 lbs
Ship Dimensions	Height (Inches)	0
	Width (Inches)	0
	Length (Inches)	0
	Cube (Cubic Feet)	0
	Weight	
Storage Temp		-40 to 140 °F
Operating Temp		32 to 104 °F
Storage Humidity		
Operating Humidity		
Amps		
Watts		
Power Requirements	12 VDC - Vehicle cigarette lighter outlet; Portable Battery (12 hr) in optional Battery Pack; or 115 VAC with Battery Eliminator (AC adapter)	
Notes	Warm up time: 2 ± 1 hour Capacity: 42 - 50 mm petri dishes or 40 - MPN tubes (19 mm OD) or 6 - P/A disposable bottles	

**WATER TEST KIT BACTERIO (NEEDS START-UP)
ACCESSORIES & CONSUMABLES - STARTUP**

NSN	Notes	Nomenclature	MFR PN Hach Company	DIST PN	QTY	UI	Unit Price	Total Price	Shelf Life Mon.
6665014673718		Water Test Kit Bacterio (Needs Start-up)	2569700		1	EA	\$3,615.26	\$3,615.26	N/A
	Awaiting NSN	Bag Sterile Whirl-Pak w Dechlorinating Agent 177mL 100s	2075333		2	100/PG	Incl w/NSN 3718	Incl w/NSN 3718	N/A
	Awaiting NSN	Burner Alcohol	2087742		1	EA	Incl w/NSN 3718	Incl w/NSN 3718	N/A
	Awaiting NSN	Counter Hand Tally	1469600		1	EA	Incl w/NSN 3718	Incl w/NSN 3718	N/A
6640015523532		Incubator Portable 12VDC 30 to 50 C	2569900		1	EA	Incl w/NSN 3718	Incl w/NSN 3718	N/A
7610015437003	Download Below	Literature Operator	25699-88		1	EA	Incl w/NSN 3718	Incl w/NSN 3718	N/A
NSN Not Required		Rack Petri Dish General Purpose	2580502		1	EA	Incl w/NSN 3718	Incl w/NSN 3718	N/A

**WATER TEST KIT BACTERIO (NEEDS START-UP)
WEB LINKS**

[Click Below for Operator Manual](#)

<http://www.hach.com/hc/view.file.invoker/2569700/View=FILECAT_MAN_MICROBIOLOGY-MEL_PORTABLE_LABS/NewLinkLabel=Portable+Laboratory+Manuals/SESSIONID|BzVNekkwT1R-FME1qZzJKbWQxWlhOMFNGQIJRZz09QkVzeE1qUQ==|>>